Order information:

Catalog No.	Contents
4705	R1 4 x 50 ml R2 4 x 10 ml
	R4 1 x 5 ml

Intended use:

Enzymatic in vitro test for the quantitative determination of urea in human serum, plasma and urine.



Summary:

The determination of urea is the most widely used test for the evaluation of kidney function. The test is frequently used in conjunction with the determination of creatinine for the differential diagnosis of prerenal hyperuremia (cardiac decompensation, water depletion increased protein catabolism), renal hyperuremia (glomerulonephritis, chronic nephritis, polycystic kidney, nephrosclerosis, tubular necrosis) and postrenal hyperuremia (obstructions of the urinary tract). Urea is the final degradation product of protein and amino acid metabolism. In protein catabolism the proteins are broken down to amino acids and deaminated. The ammonia formed in this process is synthesized to urea in the liver. This is the most important catabolic pathway for eliminating excess nitrogen in the human body.

In 1914 Marshall introduced an assay based on the enzyme urease for determining urea in blood. The ammonia released from urea by urease was measured titrimetrically. Numerous other techniques have since been employed to measure the ammonia produced. These include Bertholot's indophenol assay and the reaction of ammonia with Nessler's reagent. Subsequent modifications have been published by Fawcett and Scott and by Chaney and Marbach. In 1995, Talke and Schubert published a totally enzymatic procedure for the determination of urea using the coupled urease/glutamate dehydrogenase (GLDH) enzyme system. The Analyticon UREA-UV assay is based on the completely enzymatic method. It has been optimized for automatic analyzers that permit kinetic measurements.

Test principle:

Urea is hydrolysed in presence of urease to produce ammonia and CO₂. The ammonia produced combines with 2-oxoglutarate and NADH in presence of GLDH to yield glutamate and NAD.

The decrease in absorbance due to the decrease of NADH concentration in unit time is proportional to the urea concentration.

Reagent concentration:

R1:	
Tris buffer pH 8.0	50 mmol/l
2-oxoglutarate	15 mmol/l
Urease	≥ 1000 U/l
GLDH	≥ 6000 U/l
R2:	
NADH	0.18 mmol/l
R4:	
Urea	50 mg/dl (8.325 mmol/l)

Preparation and stability:

R1: ready for use

R2: ready for use

R4: ready for use

The reagents are stable up to the expiry date on the label when stored at 2-8°C.

Onboard stability:	R1	28 days
	R2	28 days

Alternatively the reagent can be used as single reagent (manual procedure). Mix 5 volumes enzyme reagent/R1 with one volume starting buffer /R2. The working solution is stable up to:

14 day	at + 2°C to +8°C
3 days	at +20°C to +25°C

Specimen:

Serum

Collect serum using standard sampling tubes. Heparin, citrated plasma.

Stability:	7 days	at +20°C to +25°C
	7 days	at +2°C to +8°C
	1 year	at -20°C

Urine

Urine diluted 1:50 with distilled water.

Collect urine without using preservatives.

Stability:	2 days	at +20°C to +25°C
	7 days	at +2°C to +8°C
	1 month	at -20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

Xn - Harmful

All Reagents (R1, R2 and R4) contain harmful components (sodium azide). R22: Harmful if swallowed.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Ammonia produced on the cuvette during a GLDH or lactate UV determination interferes with the urea-UV assay. The urea-UV must therefore not be installed on the analyzers together with reagents for the GLDH or lactate UV test. In urine endogenous ammonium ions interfere with the urea-UV assay. Elevated concentrations can occur under acidic conditions (e.g. acidosis.)

Great care must be taken to prevent ammonia contamination of the specimens and calibrators to be analyzed for urea/urea nitrogen.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure

Wavelength:	340nm, Hg 334 nm
Temperature:	+25 °C or +37 °C
Cuvette:	1 cm light path
Zero adjustment:	distilled water

	Standard/Calibrator	Sample
Working reagent	1000 µl	1000 µl
Standard R4	10 µl	----
Sample	----	10 µl

Mix, and read A₁ after exactly 30 sec., read A₂ after additional 60 sec.

Calculation:

Determine the absorbance change as

$$\Delta A_{\text{sample}} = [A_1 (\text{sample}) - A_2 (\text{sample})]$$

$$\Delta A_{\text{calibrator}} = [A_1 (\text{calibrator}) - A_2 (\text{calibrator})]$$

and use this for calculation.

$$\frac{\Delta A_{\text{sample}}}{\Delta A_{\text{calibrator}}} \times \text{calibrator conc.} = \text{Urea conc. (mg/dl)}$$

Measuring/reportable range

Kinetic measuring: linear up to 400 mg/dl (66.7 mmol/l)

Conversion in SI-units, relation between urea and urea-nitrogen:

$$\text{mg/dl} \times 0.166 = \text{mmol/l (urea)}$$

$$\text{mg/dl urea} \times 0.467 = \text{mg/dl urea-nitrogen}$$

Expected values:

Serum

15 – 45 mg/dl (1.7 – 8.3 mmol/l) urea

4.7 – 2300 mg/dl (31.1-15.3 mmol/l) urea-nitrogen

Reference ranges for children are given in the brochure "Reference ranges for adults and children; preanalytical considerations" by Heil W. Koberstein R, Zawta B. (published by Boehringer Mannheim GmbH 1997).

24-hour urine

20 – 36 g/l (333 – 600 mmol/l) urea

9.4 – 16.9 g/l (6.3-11.3 mmol/l) urea-nitrogen

The expected values are influenced by daily take-up of proteins on relation to the body weight.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the urea results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 0.5 mg/dl

The lower detection limit represents the lowest measurable urea activity that can be distinguished from zero.

Imprecision:

Urea

Reproducibility was determined using human samples and controls between day.
The following results were obtained:

Serum

	Between day		
Sample	Mean mg/dl	SD mg/dl	CV %
Human serum	31	1.1	3.4
Contronorm	50	0.9	1.8
Controptath	145	1.6	1.1

Urine

	Between day		
Sample	Mean mg/dl	SD mg/dl	CV %
Human serum	692	13.1	1.9
Contronorm	1753	25.6	1.5
Controptath	1813	23.9	1.3

Urea-N

Reproducibility was determined using human samples and controls between day.
The following results were obtained:

Serum

	Between day		
Sample	Mean mg/dl	SD mg/dl	CV %
Human serum	17	0.56	3.3
Contronorm	20	0.50	2.4
Controptath	66	0.82	1.2

Urine

	Between day		
Sample	Mean mg/dl	SD mg/dl	CV %
Human serum	568	17.94	3.2
Contronorm	500	3.68	0.7
Controptath	769	10.58	1.4

Quality Control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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