Hemolyzer-Hypocleaner

Hemolyzer-Hypo-Reiniger

Order information:

Catalog-No.	Content
HE3401	100 ml

Intended use:

Concentrated hypochlorite cleaning reagent for automated analysis on Hemolyzer instruments.

Hemolyzer-Hypocleaner is a stabilized and micro-filtered concentrated hypochlorite solution for intensive oxidative automated cleaning, rinsing and washing of capillaries, tubing and chambers in Analyticon hematology analyzers, removing blood component precipitates and lipoprotein deposits ("Emergency cleaning").

Emergency cleaning procedure is advisable, if the blank measurement result is too high or unacceptable results appear during the measurement or calibration. Please refer to the instrument's User Manual for further information.

Reagent Concentration:

Detergents	< 1.0 %
Sodium hypochlorite	< 12.0 %
Sodium hydroxide	< 3.5 %
Stabilizers	< 0.5 %
in ion-free water	

Preparation and Stability:

Reagent is ready for use

Leave the reagent at room temperature for at least 12 hours before use.

Stability:		
Unopened container:	Up to the expiration date	at 15 – 30 °C
Opened container:	120 days	at 15 – 30 °C

Do not freeze!

Cleaning procedure:

During daily maintenance procedure put a proper sample tube containing Hemolyzer-Hypocleaner into the sample rotor and aspirate the cleanser through the probe needle. Perform cleaning cycle as specified in the instrument's User Manual.

If the blank is too high or unacceptable results appear pipette Hemolyzer-Hypocleaner directly into the appropriate chamber according to the instrument's User Manual.

Wait for 2 minutes to perform an efficient cleaning and measure the blank.

Specimen:

Hemolyzer-Hypocleaner is intended for use with blood specimens collected by vein puncture in EDTA anticoagulant. Specimens for hematological analysis may be stored for up to 5 hours at 15-30 °C or up to 24hours after collection when refrigerated (2-8 °C).

Notes:

For in vitro diagnostic use only.
For trained laboratory professional only.
Exercise the normal precautions required for handling all laboratory reagents.

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