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#### BIBLIOGRAFIA

- Thompson JM. The control of heparin therapy by the activated partial thromboplastin time: sensitivity of various thromboplastins to heparin, Triplett DA. Ed. Standardisation of coagulation assays: an overview. Skokie IL: College of American Pathologists 1982; 195-206
- Triplett DA. Heparin: Clinical use and monitoring. Ed. Triplett DA. Laboratory evaluation of coagulation. Chicago IL: American society of clinical pathologists, 1982; 217-313
- Poller L. Activated partial thromboplastin time (aPTT). Laboratory techniques in thrombosis a manual, 2a edizione, the Netherlands: Kluwer academic Publishers, 1999; 37-43
- NCCLS. Collection, transport and processing of blood specimens for coagulation testing and performance of coagulation assays, 2a edizione H21-A3 1991.
- Thompson JM, Poller L. The activated partial thromboplastin time. Thompson JM ed. Blood coagulation and haemostasis: a practical guide. Edinburgh: Churchill Livingstone, 1985; 301-9
- Hathaway WE, Bonnar J., Physiology of coagulation in the fetus and newborn infant., Perinatal coagulation. New York: Grune & Stratton, 1978; 69

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
#### REFERENCIAS

- Thompson JM. The control of heparin therapy by the activated partial thromboplastin time: sensitivity of various thromboplastins to heparin, Triplett DA. Ed. Standardisation of coagulation assays: an overview. Skokie IL: College of American Pathologists 1982; 195-206
- Triplett DA. Heparin: Clinical use and monitoring. Triplett DA ed. Laboratory evaluation of coagulation. Chicago IL: American society of clinical pathologists, 1982; 217-313
- Poller L. Activated partial thromboplastin time (aPTT). Laboratory techniques in thrombosis a manual, 2nd edition, the Netherlands: Kluwer academic Publishers, 1999; 37-43
- NCCLS. Collection, transport and processing of blood specimens for coagulation testing and performance of coagulation assays, 2nd edition H21-A3 1991.
- Thompson JM, Poller L. The activated partial thromboplastin time, Thompson JM ed. Blood coagulation and haemostasis: a practical guide. Edinburgh: Churchill Livingstone, 1985; 301-9
- Hathaway WE, Bonnar J., Physiology of coagulation in the fetus and newborn infant., Perinatal coagulation. New York: Grune & Stratton, 1978; 69

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