Calcium Chloride (0.025M) (English) CE

REF	CG0024	CONT	10x4mL
------------	---------------	-------------	---------------

INTENDED USE

Calcium Chloride (0.025M) is a reagent for use only in the *in vitro* determination of prothrombin time (PT), activated partial thromboplastin time (APTT) and re-calcification (plasma clotting) time, as well as in various coagulation factor assays ^(1,2).

REAGENT

Calcium Chloride (0.025M) is a liquid preparation of 0.025M calcium chloride and 0.095% sodium azide as a preservative. The lot number and expiry date of the reagent are shown on the vial's label. The reagent is provided ready to use.

PRECAUTIONS

- Do not ingest. Avoid contact with skin, eyes or clothing.
- Contains sodium azide. Sodium Azide may react with lead joints in copper drain lines to form explosive compounds. Drains should be well flushed with water when discarding the reagent.

STORAGE AND STABILITY

- The **Calcium Chloride (0.025M)** is stable until expiry date when stored in the tightly capped original container at 2 to 8°C.
- The reagent should not be used if cloudy or discoloured.
- This reagent contains 0.095% sodium azide as a preservative.

Stability of Calcium Chloride (0.025M)		
Unopened	2 – 8°C	Refer to vial label
Open vial	2 – 8°C	7 Days
On-board Coagulyzer® 100	18 - 26°C	5 Days
On-board Coagulyzer® 1,2,4 Instruments	37°C	7 Days*

*Opened vial must be resealed after measurements are completed.

PROCEDURE

Refer to your instrument manual, laboratory procedures guide and to the specific test insert for detailed instructions.

REFERENCES

- Poller L. Activated partial thromboplastin time (aPTT). Laboratory techniques in thrombosis a manual, 2nd edition, the Netherlands: Kluwer academic Publishers, 1999; 37-43
- Thompson JM, Poller L. The activated partial thromboplastin time, Thompson JM ed. Blood coagulation and haemostasis: a practical guide. Edinburgh: Churchill Livingstone, 1985; 301-9

WARRANTY

This product is warranted to perform in accordance with its labeling and literature. **Analyticon** disclaims any implied warranty of merchantability or fitness for any other purpose. Purchaser must calibrate and determine the suitability of **Analyticon**'s products for their specific applications. In no event will **Analyticon** be liable for any consequential damages arising out of aforesaid express warranty.

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REF	CG0024	CONT	10x4mL
------------	---------------	-------------	---------------

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