Combi-Screen® VET

Zur Veterinär-In-Vitro Diagnostik

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Combi-Screen® VET

For Veterinary In-Vitro Diagnostic Use

Urine Test Strips for the Rapid Determination of Ascorbic Acid, Bilirubin, Blood, Glucose, Ketones, Leucocytes, Nitrite, pH-value, Protein, Specific Gravity and Urobilinogen. Refer to the carton and label for specific parameter combination on the product you are using.

Intended Use

For use as a preliminary screening test for diabetes, liver diseases, haemolytic diseases, urogenital and kidney disorders and metabolic abnormalities for animals.

Procedure and Notes

- Use only well mixed, non-centrifuged urine, which should not be older than 4 hours. First morning urine is recommended. Protect the samples from light.
- If the samples cannot be tested immediately, they should be stored at 2 – 4 °C and brought to room temperature (15 – 25 °C) before testing.
- Collect specimens in clean, well rinsed containers, free of detergents. Do not add any preservatives.

Parameter	Horse	Cattle	Pig	Sheep	Goat	Dog	Cat	Bunny	Guinea pig
Ascorbic Acid)))))))))
Bilirubin	neg.	neg. – weak pos.	neg.	neg.	neg.	neg. – weak pos.	neg.	neg.	neg.
Blood	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.
Glucose	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.
Ketones	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.
Leucocytes	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.
Nitrite	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.
pH	7.6 – 9.0	7.0 – 8.4	5.5 – 8.0	7.5 – 8.5	7.5 – 8.5	5.5 – 7.0	5.0 – 7.0	8.2	8.0 – 9.0
Protein	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.
Specific Gravity	1.020 – 1.040	1.020 – 1.040	1.020 – 1.040	1.020 – 1.040	1.020 – 1.040	1.030 – 1.045	1.035 – 1.060	1.003 – 1.036	1.000 – 1.040
Urobilinogen	neg. – weak pos.	neg. – weak pos.	neg. – weak pos.	neg. – weak pos.	neg. – weak pos.	neg. – weak pos.	neg. – weak pos.	neg. – weak pos.	neg. – weak pos.
Color	Clay – ochre	Bright yellow – brown yellow	Bright yellow – dark yellow	Bright yellow – brown yellow	Bright yellow – brown yellow	Paly yellow – brown yellow	Yellow – dark yellow	Bright yellow – auburn	Yellow
Turbidity	Turbid	Clear	Clear	Clear	Clear	Clear	Clear	Clear – turbid	Clear – slightly turbid
24h-Volume (ml/kg)	8 – 30	16 – 50	20 – 80	10 – 40	10 – 40	24 – 50	18 – 25	20 – 350	–
Odor	Aromatic	Aromatic	Aromatic, unpleasant	Aromatic – neutral	Aromatic – neutral	Garlic, meat broth	Caustic	–	–

Special Information on veterinary urine diagnostics

Besides the information on Combi-Screen® VET test strips, there is some additional information especially related to veterinary diagnostics. The nutrition, habits of liquid intake, physical activity of the animal and conditions of collection of the sample must be considered during unalysis of animals.

⁾ **Ascorbic Acid**: Most of the animals are in the position to produce Ascorbic Acid with its body in limited quantities independently (Exceptions: Guinea pig, Primates and some species of birds). By the influence of heat, coldness, disease, stress and deficits in feeding it could come to an acute shortage of vitamin C for the animal. Vitamin C is indispensable for the creation of connective tissue, bones and cartilage as well as for normal function of tissue and several metabolism reactions. In connection with vitamin B and vitamin E it is stimulating the defense reaction of the body after the influence of stress and infections. Most of the animal feeds are containing only a few or no vitamin C. The vitamin could be added with supplemental animal feed, fresh fruits and vegetables.

Bilirubin: The test normally is negative for all animals. Dogs and cattle may have bilirubin in their urine without being ill.

Blood: The test normally is negative for all animals.

Glucose: The test normally is negative for all animals.

Ketones: The test normally is negative for all animals. Ruminants sometimes may show a weak positive test result.

Leucocytes: The test normally is negative for all animals. The validity of the determination of leucocytes in dogs and cats is very limited. A microscopic inspection is recommended.

Nitrite: The test normally is negative for all animals. Especially for carnivores, the test may show a limited sensitivity, because the necessary nitrate concentration is not present in their urine in case of a non-herbal diet.

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General information for urinalysis:

Clinical Utility, Test Principles, Expected Values, Limitations

Ascorbic Acid: Intended to measure the level of ascorbic acid (vitamin C) in urine. The detection is based on the decoloration of Tillmans reagent. In the presence of ascorbic acid a color change takes place from grey blue to orange. As ascorbic acid can disturb various test fields, the test must be repeated if the ascorbic acid reaction is positive, however, at the earliest 10 hours after the last vitamin C intake (medication, fruit, vegetables). Values of 5–10 mg/dl or 0.6–1.1 mmol/l are indicated.

Bilirubin: Intended to measure the levels of bilirubin conjugates in urine. Measurements of urinary bilirubin and its conjugates are used in the diagnosis and treatment of certain liver and bile diseases. A red azo compound is obtained in the presence of acid by coupling of bilirubin with a diazonium salt. Normally, no bilirubin is detectable in urine. Concentrations of 0.5 mg/dl and more lead to a color of red-orange peach and indicate the early stage of a liver disease. The reaction is unaffected by pH of urine. False low or negative results may be simulated by large amounts of vitamin C or Nitrite or by longer exposure of the sample to direct light. Increased concentrations of urobilinogen can reinforce the sensitivity of the test field. Different urine contents (e.g. urine indicans) can lead to atypical coloration. For metabolites of drugs see urobilinogen. The color fields correspond to the following values: 0 (negative), 1(+), 2(++), 4(+++) mg/dl or 0 (negative), 17(+), 35(++), 70(+++) µmol/l. Values of 0.5–1 mg/dl (8.5–17 µmol/l) Bilirubin are indicated.

Blood: Intended to detect occult blood in urine. Occult blood indicates serious urological or kidney diseases. Microhaematuria does not affect the colour of urine and is only detect

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