# ACP/ACP-P

## ACID/PROSTATIC PHOSPHATASE

### Order information:

Catalog No.	Contents					
927	R1	1 x	65 ml	R2	20 x	for 3 ml
	R3	1 x	30 ml	R4	1 x	5 ml
H0101 Hit I (ILab*)	R1	6 x	13 ml	R2	6 x	for 13 ml
	R3	3 x	13 ml	R4	1 x	5 ml
H0103 Hit 917 (AU*)	R1	6 x	13 ml	R2	6 x	for 13 ml
	R3	3 x	13 ml	R4	1 x	5 ml

(\* ) Kit contains only reagent barcodes for Hitachi systems.

### System information:

Hitachi 911: ACN 001 (ACP) / ACN 002 (NPP)

Hitachi 917: ACN 021 (ACP) / ACN 022 (NPP)

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

### Intended use:

In vitro test for the quantitative determination of acid phosphatase and prostate-specific acid phosphatase in human serum.

### Summary:

Serum acid phosphatase consists of 5 isoenzymes that originate mainly from erythrocytes, platelets, spleen and liver reticuloendothelial cells, the kidneys, bone, and prostate epithelial cells. Prostatic acid phosphatase isoenzyme 2 is formed mainly, but not exclusively, in the prostate.

In general, total acid phosphatase and prostatic acid phosphatase levels increase in the presence of progressive and metastasizing prostate carcinoma, the increase being dependent upon the disease stage in 80% of patients with metastasizing prostate cancer. The percentage increase at each stage depends on the classification (pathological or clinical). Increased acid phosphatase levels occur in Gaucher's disease, Niemann-Pick disease, 1-2 days after prostate surgery, biopsy, manipulation or catheterization, and in the presence of benign prostate hypertrophy, prostatitis, and prostate infarction.


The assay used here is a modification of the method described by Hillmann. Addition of 1.5-pentandiol increases the activity of prostatic acid phosphatase.

### Test principle:

1-naphthylphosphate is hydrolyzed by acid phosphatase to phosphate and 1-naphthol, which is converted with Fast-Red TR to an azo dye.

The increase of absorbance is proportional to the acid phosphatase (ACP) activity in the sample.

The prostatic fraction of the ACP can be blocked by tartrate. The activity of this fraction has to be calculated by the difference of the activity determined in reagent without and reagent with tartrate.

### Reagent Concentration:

<b>R1:</b>	
Citrate buffer pH 5.0	150 mmol/l
1,5-Pentane diol	200 mmol/l
Detergents	0.3%

<b>R2:</b>	
1-Naphthyl Phosphate	10 mmol/l
Fast Red TR	1.1 mmol/l

<b>R3:</b>	
Citrate buffer, pH5.0	150 mmol/l
Sodium tartrate	90 mmol/l
1,5-Pentane diol	200 mmol/l
Detergents	0.3%

<b>R4:</b>	
Citrate-Acetate buffer, pH 5.0	1.0 mol/l

### Preparation and stability:

For total acid phosphatase dissolve R2 with the corresponding volume of R1. For non-prostatic acid phosphatase dissolve R2 with the corresponding volume of R3. Gently swirl working solution until is completely dissolved. DO NOT SHAKE!

Unopened kit components: Up to the expiration dates at +2°C to +8°C

On board stability: 5 days Working reagent R1/R2  
5 days Working reagent R3/R2

Store the reagent refrigerated and protected from light.

### Specimen:

Collect serum with the standard sampling tubes. Hemolysis will interfere.

Perform determinations on the samples immediately. Samples which cannot be examined immediately should be stabilized as follows:

Add one drop (40µl) of R4 per 1 ml of fresh taken sample and mix.

Stability: 8 days at +2°C to +8°C

Centrifuge samples containing precipitates before performing the assay.

### Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

### Limitation interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 5 (approximate conjugated bilirubin concentration: 5 mg/dl)

Hemolysis: No significant interference up to an index H of 60 (approximate hemoglobin concentration: 60 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 400 (approximate triglycerides concentration: 800 mg/dl) There is poor correlation between turbidity and triglycerides concentration.

The addition of stabilizer to the sample interferes with the determination of other parameters.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

Applications for automated systems are available on request.

#### Materials provided

• Working solutions as described above

#### Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure (short time kinetic):			
Wavelength:	Hg 405 nm		
Temperature:	+37°C		
Cuvette:	1 cm light path		
Zero adjustment:	air or distilled water		
	Macro	Semi	Micro
Working reagent	2500 µl	1000 µl	500 µl
Sample	250 µl	100 µl	50 µl
Mix and incubate 3 min, then read A <sub>1</sub> . Incubate exactly 3 min. and read A <sub>2</sub> . Calculate ΔA/min. = (A <sub>2</sub> - A <sub>1</sub> ) /3.			
<b>Calculation for sample start:</b>			
For total ACP activity:			
U/l = (ΔA /min) x 643			
For prostatic phosphatase activity:			
U/l = (ΔA /min <sub>TOTAL</sub> - dA/min <sub>NONPROSTATIC</sub> ) x 643			

### Measuring range:

Short time kinetic: 1-150 U/l (0.017-2.50 µkat/l)

Roche/Hitachi: 0.1-180 U/l (0.002-3.00 µkat/l)

Determine samples having higher activities via the rerun function. On instruments without rerun function, manually dilute the samples with 0.9% NaCl or distilled/deionized water (e.g. 1 + 2). Multiply the result by the appropriate dilution factor (e.g. factor 3).

### Analytical sensitivity (lower detection limit):

0.1 U/l or 0.002 µkat/l

The lower detection limit represents the lowest measurable activity that can be distinguished from zero using a Hitachi.

### Imprecision:

Reproducibility was determined using human samples within run. The following results were obtained:

Sample	Within run		
	Mean U/l	SD U/l	CV %
Control serum I	21.8	0.07	0.34
Control serum II	30.5	0.22	0.46
Control serum III	36.7	0.17	0.71

Sample	Between day		
	Mean U/l	SD U/l	CV %
Control serum I	22.0	0.31	1.41
Control serum II	25.3	0.48	1.90
Control serum III	40.5	0.45	1.33

# ACP/ACP-P

## ACID/PROSTATIC PHOSPHATASE

### Reference value:

#### Total ACP:

Men <6.6 U/l (0.110 µkat/l)  
Women <6.5 U/l (0.108 µkat/l)

#### Prostatic Phosphatase:

Men <3.5 U/l (0.058 µkat/l)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic examinations and other purposes, albumin results should always be assessed in conjunction with the patient's medical history and clinical findings.

### Method comparison:

A comparison of the Analyticon ACP (y) with a commercial obtainable assay without Petandiole activation (x) gave the following result:

$$y = 0.867 x + 1.09; \quad r = 0.986$$

### Quality Control:

#### Human Control Serum:

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### Calibration:

S1: 0.9% NaCl		
S2: Bio Cal E	10 x 3 ml	#1430

### Calibration frequency:

A two-point-calibration is recommended in case of:

- change of lot
- quality control requirements

### Disposal:

Please note the legal regulations.

### Literature:

1. Bablok W. et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
2. Glick M.R., Ryder K.W., Jackson S.A.. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474.
3. Guder W.G., Narayanan S., Wisser H., Zawta B. List of Analytes Preanalytical Variables. Broschüre in: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag, 1996.
4. Heller J.E. Prostatic acid phosphatase: Its current clinical status.
5. Hillmann G. Z. klein. Chem. u. klin. Biochem. 1971;9:273.
6. J. Urol. 1987;137:1091-1103.
7. Junge W., Thormeyer I., Schlottmann A. et al. Determination of Reference Values for Acid Phosphatase using a New Photometric Assay. Pecs, Hungary: 3<sup>rd</sup> Alpe-Adria Congress on Clinical Chemistry and Laboratory Medicine September 7-9, 1994.
8. Passing H., Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720.
9. Tietz N.W. Clinical Guide to Laboratory Tests, 3<sup>rd</sup> Philadelphia: W.B. Saunders Company, 1995:516-519.

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