### Order information:

Catalog No.	Contents
H0501 Hit I (ILab*)	R1 12 x 50 ml   R2 6 x 22 ml
H0503 Hit 917 (AU*)	R1 6 x 60 ml   R2 6 x 15 ml
AU0503 AU	R1 6 x 60 ml   R2 6 x 15 ml

(\* ) Kit contains only reagent barcodes for Hitachi systems.

### System information:

Hitachi 911/917: ACN 413  
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

### Intended use:


In vitro test for the quantitative determination of albumin in human serum and plasma.

### Summary:

Albumin is a carbohydrate-free protein, which constitutes 55-65% of total plasma protein. It maintains oncotic plasma pressure, is also involved in the transport and storage of a wide variety of ligands and is a source of endogenous amino acids. Albumin binds and solubilizes various compounds, e.g. bilirubin, calcium and long-chain fatty acids. Furthermore albumin is capable of binding toxic heavy metals ions as well as numerous pharmaceuticals, which is the reason why lower albumin concentrations in blood have a significant effect on pharmacokinetics. Hyperalbuminemia is of little diagnostic significance except in the case of dehydration. Hypoalbuminemia occurs during many illnesses and is caused by several factors: Compromised synthesis due either to liver disease or as a consequence of reduced protein uptake, elevated catabolism due to tissue damage (severe burns) or inflammation, malabsorption of amino acids (Crohn's disease), proteinuria as a consequence of nephrotic syndrome; protein loss via the stool (neoplastic disease). In severe cases of hypoalbuminemia, the maximum albumin concentration of plasma is 2.5g/dl. Due to the low osmotic pressure of the plasma water permeates through blood capillaries into tissue (edema). The determination of albumin allows monitoring of a controlled patient dietary supplementation and serves also as an excellent test of liver function.

### Test principle:

Colorimetric assay, endpoint method  
• Sample and addition of R1 (buffer)  
• Addition of R2 (substrate) and start of the reaction:  
At a pH value of 4.1 albumin displays a sufficiently cationic character to be able to bind with bromocresol green (BCG), an anionic dyestuff, to form a blue-green complex.

The colour intensity of the blue-green colour is directly proportional to the albumin concentration and can be determined photometrically.

### Reagent Concentration:

<b>R1:</b>	
Succinate buffer, pH 4.1	75 mmol/l
Brij 35	7 ml/l
Detergents and stabilizers	
<b>R2:</b>	
Succinate buffer, pH 4.1	75 mmol/l
Bromocresol green	0.15 mmol/l
Brij 35	7 ml/l
Detergents and stabilizers	

### Preparation and stability:

R1: Ready for use  
R2: Ready for use  
Unopened kits: up to the expiration date at +2°C to +25°C.  
Onboard stability: R1: 28 days.  
R2: 28 days.

### Specimen:

Collect serum using standard sampling tubes.  
Heparin or EDTA plasma. Separate serum or plasma from the clot or cells within one hour and analyse immediately, or store as follows:  
Stability: < 3 days at +4°C  
6 months at -20°C

Centrifuge samples containing precipitates before performing the assay.

### Notes:

For in vitro diagnostic use.  
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.  
Exercise the normal precautions required for handling all laboratory reagents.

### Limitation interference:

Criterion: Recovery within ± 10% of initial value at an albumin concentration of 3.5g/dl.

Icterus: No significant interference up to an index I of 92 (approximate bilirubin concentration: 92 mg/dl).

Hemolysis: No significant interference up to an index H of 1100 (approximate hemoglobin concentration: 1100 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 1075 (approximate triglycerides concentration: 2150 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

Applications for automated systems are available on request.

#### Materials provided

• Working solutions as described above

#### Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

### Manual procedure:

Wavelength:	628 nm, Hg 623 (620 – 640 nm)
Temperature:	+20 to +37°C
Cuvette:	1 cm light path
Zero adjustment:	each series needs one reagent blank only

	Blank	Sample / Calibrator
Sample / Calibrator	-	10 µl
R1	1000 µl	1000 µl
R2	200 µl	200 µl

Mix and incubate 10 minutes. Read the absorbance against blank within 30 minutes.

### Calculation:

$$\frac{\Delta A \text{ sample}}{\Delta A \text{ Calibrator}} \times \text{Calibrator conc.} = \text{Albumin in g/dl}$$

### Measuring range:

0.2g/dl – 6.0 g/dl or 2 g/l – 60 g/l  
Determine samples having higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl solution (e.g. 1 + 1). Multiply the result by the appropriate dilution factor (e.g. 2).

### Reference value:

#### Expected values according to Tietz'

Adults		3.4-4.8 g/dl or 34-48 g/l
Newborn	0 - 4 days	2.8-4.4 g/dl or 28-44 g/l
Children	4 days - 14 years	3.8-5.4 g/dl or 38-54 g/l
	14 - 18 years	3.2-4.5 g/dl or 32-45 g/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, albumin results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

### Analytical sensitivity (lower detection limit):

0.2 g/dl or 2 g/l  
The lower detection limit represents the lowest measurable albumin concentration that can be distinguished from zero. It is calculated as three standard deviations of the lowest standard.

### Imprecision:

Reproducibility was determined using in an internal protocol. The following results were obtained.

Sample	Within run		
	Mean g/dl	SD g/dl	% CV
Control serum 1	3.35	0.04	1.16
Control serum 2	3.50	0.03	0.94
Control serum 3	3.71	0.03	0.83

Sample	Between Day		
	Mean g/dl	SD g/dl	% CV
Control serum 1	3.37	0.04	1.08
Control serum 2	3.52	0.04	1.17
Control serum 3	3.98	0.03	0.84

### **Method comparison:**

A comparison of the Analyticon Fluitest ALB (y) with a commercial obtainable assay

(x) gave the following result:

$$y = 0.914 x + 0.306; \quad r = 0.996$$

### **Quality Control:**

Human Control Serum:

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320
Protein Control Level 1	3 x 1 ml	#7661
Protein Control Level 2	3 x 1 ml	#7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### **Calibration:**

Standardization: The albumin method was standardized against the CRM 470 reference preparation.

S1:	0.9% NaCl		
S2:	Bio Cal®	20 x 3 ml	#1420
	Bio Cal® E	10 x 3 ml	#1430
	Bio Cal® P	3 x 1 ml	#1470

### **Calibration frequency:**

A two-point-calibration is recommended in case of:

- change of lot
- quality control requirements

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