Order information:

Catalog No.	Contents
9136	R1 6 x 100 ml R4 1 x 5 ml

Intended use:

In vitro test for the quantitative determination of albumin in human serum and plasma.


Summary:

Albumin is a carbohydrate-free protein, which constitutes 55-65% of total plasma protein. It maintains oncotic plasma pressure, is also involved in the transport and storage of a wide variety of ligands and is a source of endogenous amino acids. Albumin binds and solubilizes various compounds, e. g. bilirubin, calcium and long-chain fatty acids. Furthermore albumin is capable of binding toxic heavy metals ions as well as numerous pharmaceuticals, which is the reason why lower albumin concentrations in blood have a significant effect on pharmacokinetics. Hyperalbuminemia is of little diagnostic significance except in the case of dehydration. Hypoalbuminemia occurs during many illnesses and is caused by several factors: Compromised synthesis due either to liver disease or as a consequence of reduced protein uptake, elevated catabolism due to tissue damage (severe burns) or inflammation, malabsorption of amino acids (Crohn's disease), proteinuria as a consequence of nephrotic syndrome; protein loss via the stool (neoplastic disease). In severe cases of hypoalbuminemia, the maximum albumin concentration of plasma is 2.5g/dl. Due to the low osmotic pressure of the plasma water permeates through blood capillaries into tissue (edema). The determination of albumin allows monitoring of a controlled patient dietary supplementation and serves also as an excellent test of liver function.

Test principle:

Colorimetric assay, endpoint method

At a pH value of 4.2 albumin displays a sufficiently cationic character to be able to bind with bromocresol green (BCG), an anionic dyestuff, to form a blue-green complex.

The color intensity of the blue-green color is directly proportional to the albumin concentration and can be determined photometrically

Reagent Concentration:

R1:	
Succinate buffer, pH 4.2	75 mmol/l
Bromocresol green	0.15 mmol/l
Brij 35	7 ml/l
Detergents and stabilizers	>0.1 %
R4	
Bovine albumin concentration according to CRM 470 (IFCC)	4.5 g/dl

Preparation and stability:

All reagents are ready to use.

The reagent is stable at + 2°C to + 8°C up to the date of expiration specified.

Avoid direct sunlight.

On board stability: 28 days

Specimen:

Collect serum using standard sampling tubes.

Heparin or EDTA plasma. Separate serum or plasma from the clot or cells within one hour and analyze immediately, or store as follows:

< 3 days	at +4°C
6 months	at -20°C

Centrifuge samples containing precipitates before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitation interference:

Criterion: Recovery within ± 10% of initial value at an albumin concentration of 3.5 g/dl.

Icterus: No significant interference up to an index I of 92 (approximate bilirubin concentration: 92 mg/dl).

Hemolysis: No significant interference up to an index H of 1100 (approximate hemoglobin concentration: 1100 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 1075 (approximate triglycerides concentration: 2150 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

Manual procedure :

Wavelength:	628 nm (Hg623, 620-640nm)
Temperature:	+20 to +37°C
Cuvette:	1 cm light path
Zero adjustment:	reagent blank/each series needs one blank only

	Blank	Sample/Calibrator
Sample/ Calib.	---	5 µl
R1	1000 µl	1000 µl

Mix and incubate 10 minutes. Read absorbance against blank within 30min.

Calculation:

$$\frac{\Delta A \text{ sample}}{\Delta A \text{ Calibrator}} \times \text{Conc. of Calib./Standard} = \text{Albumin in g/dl}$$

Measuring:

0.2g/dl – 6.0 g/dl or 2 g/l – 60 g/l

Determine samples having higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl solution (e.g. 1 + 1). Multiply the result by the appropriate dilution factor (e.g. 2).

Reference value:

Expected values according to Tietz

Adults		3.4-4.8 g/dl or 34-48 g/l
Newborn	0 - 4 days	2.8-4.4 g/dl or 28-44 g/l
Children	4 days - 14 years	3.8-5.4 g/dl or 38-54 g/l
	14 - 18 years	3.2-4.5 g/dl or 32-45 g/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, albumin results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

0.2 g/dl or 2g/l

The lower detection limit represents the lowest measurable albumin concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using in an internal protocol. The following results were obtained.

Sample	Within run		
	Mean g/dl	SD g/dl	CV %
Control serum 1	2.04	0.034	1.66
Control serum 2	3.19	0.028	0.88
Control serum 3	3.67	0.042	1.14

Sample	Between Day		
	Mean g/dl	SD g/dl	CV %
Control serum 1	3.37	0.036	1.07
Control serum 2	3.52	0.039	1.11
Control serum 3	3.98	0.033	0.83

Method comparison:

A comparison of the Analyticon Fluitest ALB (y) with a commercial obtainable assay (x) gave the following result:
 $y = 1.102 x + 0.323$; $r = 0.988$

Quality Control:

Human Control Serum

Contronorm [®] Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220

Controptath [®] Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Protein Control Level 1	3 x 1 ml	#7661
Protein Control Level 2	3 x 1 ml	#7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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