Order information:

Catalog No.	Contents
H0701 Hit I / 917 (AU/ILab*)	R1 6 x 20 ml R2 6 x 5 ml

(*) Kit contains only reagent barcodes for Hitachi systems.

System information:

Hitachi 911/917: ACN 005

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of albumin in human urine and CSF (cerebrospinal fluid).

Summary:

Albumin is a non-glycosylated protein with a molecular weight of 66,000 daltons. It is synthesized in liver parenchymal cells at a rate of 14 g/day. Quantitatively, albumin is the most important protein component (> 50%) in plasma, CSF and urine. A small, but abnormal albumin excretion in urine is known as microalbuminuria. Causes of microalbuminuria can be glomerular (e.g. due to diabetic microangiopathy, hypertension, minor glomerular lesion), tubular (inhibition of reabsorption) or postrenal. Albumin is also a marker protein for various forms of proteinuria. In selective glomerular proteinuria, 100-3000 mg albumin/g creatinine are excreted. Non-selective glomerular proteinuria is characterized by elevated excretion of high-molecular weight proteins (IgG more than 10% of the albumin value). Pre-renal proteinuria is recognized by a discrepancy between albumin and total protein (albumin accounting for less than 30%). Simultaneous elevation of albumin and microproteins is found in glomerulotubular proteinuria occurring due to overloading of tubular reabsorption in glomerulopathy (e.g. nephrotic syndrome), combined glomerular tubulointerstitial nephropathy or in renal failure following diabetic nephropathy or other causes (overflow proteinuria). Albumin has two main functions in plasma: maintaining the oncotic pressure (80 % due to albumin in plasma) and transport. It is the most important transport protein for substances having low water solubility (such as free fatty acids, bilirubin, metal ions, hormones and pharmaceuticals). Depressed albumin levels are caused by hyperhydration, hepatocellular synthesis insufficiency, secretion disorders in the intravascular space, abnormal distribution between the intravascular and extravascular space, catabolism and loss of albumin, acute phase reactions and congenital analbuminemia. A variety of methods, such as radial immunodiffusion, nephelometry and turbidimetry, are available for the determination of albumin.

Test principle:

Immunoturbidimetric assay
Anti-albumin antibodies react with the antigen in the sample to form antigen/antibody complexes which, following agglutination, are measured turbidimetrically.

Reagent concentration:

R1:	
TRIS*-Buffer pH 7,8	100 mmol/l
PEG	4.5 %
NaCl	150 mmol/l
Preservative	
R2:	
Polyclonal anti-human-albumin-antibodies;	dependent on titer
TRIS*-Buffer pH 7.8	80 mmol/l
NaCl	150 mmol/l
Preservative	

*TRIS = Tris(hydroxymethyl)-aminomethan

Preparation and stability:

R1: Ready for use
R2: Ready for use
Unopened kit components are stable up to the expiration date at 2-8°C
On board stability at 2-8°C: R1 90 days
R2 90 days

Specimen:

Urine
Spontaneous, 24-hour urine or 2nd morning urine. Stability:

7 days	at + 20°C	to + 25°C
1 month	at + 2°C	to + 8°C

CSF

Stability: 1 day at + 20°C to + 25°C 6 days at + 2°C to + 8°C

Centrifuge samples containing precipitates before performing the assay.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.

Limitation interference:

Criterion: Recovery within ± 10% of initial value.

Urine

Icterus: No significant interference up to an approximate conjugated and unconjugated bilirubin concentration of 10 mg/dl.

Hemolysis: No significant interference up to an approximate hemoglobin concentration of 1000 mg/dl.

No interference by acetone <60mmol/l; ascorbic acid <1 g/l; creatinine <5g/l; glucose <20g/l; uric acid <700mg/l; urea <700mmol/l and urobilinogen <200mg/l.

A high-dose hook effect may occur at albumin concentrations above 3000mg/l.

17 frequently used pharmaceuticals were tested in vitro. No interference with the assay was found. The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

• Controls as indicated below
• 0.9% NaCl

Manual procedure :		
Wavelength:	340 nm	
Temperature:	37°C	
Cuvette:	1 cm light path	
Zero adjustment:	reagent blank	
	Blank	Sample/Calibrator
Sample/Calibrator	-	60 µl
R1	1000 µl	1000 µl
Mix, incubate for 5 minutes and read A ₁ . Then add:		
R2	200 µl	200 µl
Mix, incubate for 5 min. and read absorbance A ₂ .		
Calculation:		
$\Delta A = [(A_2 - A_1) \text{ sample or calibrator}] - [(A_2 - A_1) \text{ reagent blank}]$		
The concentration of albumine in urine or CSF has to be calculated from ΔE using mathematic function as logit/log or can be read from a graph using values of 4 levels of standards in the concentration range of 0 to 500 mg/l albumin. For zero value is recommended to use saline solution (0,9%)		

Measuring:

Urine

Measuring range: 3-500 mg/l

At higher concentrations manually dilute the samples with 0.9% NaCl (e.g. 1+1). Multiply the result by the appropriate dilution factor (e.g. 2) or determine via the rerun function.

Expected values:

Urine

2nd morning urine

Adults: < 20 mg albumin/g creatinine or < 2.26 g albumin/mol creatinine
Children (3-5 years) < 20 mg/l albumin

24-hour urine < 37 mg albumin/g creatinine < 20 mg/l < 30 mg/24 h

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges. For diagnostic purposes, albumin results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 3 mg/l

The detection limit represents the lowest measurable albumin concentration that can be distinguished from zero. It is calculated as the concentration lying three standard deviations above that of the lowest standard (zero standard + 3 SD, within-run precision, n = 21).

Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol between day (n = 20). The following results were obtained:

Between day			
Sample	Mean mg/l	SD mg/l	CV %
Control 1	18.0	0.78	4.38
Control 2	58.7	0.91	1.54
Control 3	187.7	2.90	1.54

Reproducibility was determined using human samples and controls in an internal protocol within run (n = 20). The following results were obtained:

Within run			
Sample	Mean mg/l	SD mg/l	CV %
Control 1	16.9	0.28	1.66
Control 2	58.9	0.91	1.54
Control 3	192.6	2.04	1.06

Method comparison:

A comparison of the Analyticon Turbitex[®] ALB (y) with a commercial obtainable assay (x) gave the following result (n=46):

$$y = 0.918 x - 2.823; \quad r = 0.999$$

Quality control:

Human Urine Control:

Urine control Set 8 x 5 ml #1507

The control intervals and limits must be adapted to the individual laboratory and country specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1-S6: Bio Cal ALB Set 5 x 1 ml #14070

Calibration frequency

Full calibration is recommended

- after lot change
- as required following quality control procedures

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