BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B1631	200 / 600	R1	6 x	47 ml
		R2	6 x	13 ml
B1633	300 / 600*	R1	6 x	55 ml
		R2	6 x	14 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

In vitro test for the quantitative determination of alkaline phosphatase (ALP) in human serum and plasma.

Summary:


Alkaline phosphatase in serum consists of four structural genotypes: the liver-bone-kidney type, the intestinal type, the placental type and the variant from germ cells. It occurs in osteoblasts, hepatocytes the kidneys, spleen, placenta, prostate, leukocytes and the small intestine. The liver-bone-kidney type is particularly important.

A rise in the alkaline phosphatase activity occurs with all forms of cholestasis, particularly with obstructive jaundice. It is also elevated in diseases of the skeletal system, such as Paget's disease, hyperparathyroidism, rickets and osteomalacia, as well as with fractures and malignant tumors. A considerable rise in the alkaline phosphatase activity is sometimes seen in children and juveniles. It is caused by increased osteoblast activity following accelerated bone growth. Various reference values for the purposes of clinical evaluation have been assigned to differing age groups.

In 1946, Bessey, Lowry and Brock published a method for the determination of alkaline phosphatase using p-nitrophenyl phosphate as substrate buffered with glycine/NaOH. In 1967, Hausamen et al improved the method by using diethanolamine as buffer. The assay described here meets the recommendations of the IFCC.

Test principle:

Colorimetric assay in accordance with a standardized method.

In the presence of magnesium and zinc ions, p-nitrophenylphosphate is hydrolyzed by phosphatases to form phosphate and p-nitrophenol. In this process AMP serves as transient phosphate acceptor. The release of coloured p-nitrophenol is proportional to the ALP activity and can be measured photometrically.

Reagent concentration:

R1:	
AMP* buffer, pH 10.44 (30°C)	1.0 mol/l
Magnesiumacetate	2.0 mmol/l
Zincsulfate	0.5 mmol/l
EDTA	2.0 mmol/l
Detergents and stabilizers	<0.1%
R2:	
p - Nitrophenylphosphate	70.0 mmol/l
Stabilized liquid	
*2-Amino-2-methyl-1-propanol	

Preparation and stability:

R1: Ready for use.
R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C

Onboard stability:	R1:	42 days
	R2:	42 days

Specimen:

Collect serum using standard sampling tubes.
Heparinized plasma.

Stability:	2 days	at +20°C to +25°C
	7 days	at +4°C to + 8°C
	4 weeks	at -20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 40 (approximate conjugated and unconjugated bilirubin: 40 mg/dl)

Hemolysis: No significant interference up to an index H of 300 (approximate hemoglobin concentration: 300 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 1200 (approximate triglycerides concentration: 2400 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below
• 0.9% NaCl

Measuring /reportable range:

30 - 450 U/l (0.5 - 7.5 µkat/l)

Determine samples having higher concentrations via the rerun function using 0.9% NaCl solution as diluent

Expected values:

		+37°C
U/l	Men	40-129
	Women	35-104
µkat/l	Men	0.67-2.15
	Women	0.58-1.74

For reference values for children, please refer to „Pediatric reference ranges“ 3. Edition. S.J. Soldin, C. Brugnara, I.M. Hicks, AAC Press.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the alkaline phosphatase results should always be assessed in conjunction with the patients' medical history, clinical examination and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 0.3 U/l or 0.005 µkat/l

The lower detection limit represents the lowest measurable ALP concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using human samples and controls within run (n = 20). The following results were obtained:

Sample	Within run		
	Mean U/l	SD U/l	CV %
Sample 1	82.5	4.21	5.1
Sample 2	242.7	3.61	1.5
Sample 3	235.5	4.89	2.1

Reproducibility was determined using human samples and controls in an internal protocol within (between day n = 20). The following results were obtained:

Sample	Run to run		
	Mean U/l	SD U/l	CV %
Sample 1	84.5	1.37	1.62
Sample 2	197.9	2.59	1.31
Sample 3	242.5	3.68	1.52

Method comparison:

A comparison of the Analyticon Fluitest® ALP IFCC (y) against a commercial obtainable IFCC assay (x) with 133 samples gave the following result:

$$y = 0.921 x + 2.308; \quad r = 0.9990$$

Quality Control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Fluitest[®] ALP IFCC

ALKALINE PHOSPHATASE

Calibration:

Calibration Type: Linear

S1: 0.9% NaCl

S2: Bio Cal[®] E 10 x 3 ml #1430

Calibration frequency:

Two-point calibration is recommended

- after reagent lot change
- as required following quality control procedures

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