Order information:

Catalog No.	Contents			
1722	R1	8 x	15 ml	R2 1 x 25 ml
1725	R1	4 x	50 ml	R2 4 x 10 ml
H1601 Hit I (ILab*)	R1	6 x	47 ml	R2 6 x 11 ml
H1602 Hit II (ILab*)	R1	6 x	95 ml	R2 6 x 21 ml
H1603 Hit 917 (AU*)	R1	6 x	60 ml	R2 6 x 14 ml
AU1603 AU	R1	6 x	60 ml	R2 6 x 14 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911/917: ACN 158

*For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

In vitro test for the quantitative determination of alkaline phosphatase (ALP) in human serum and plasma.

Summary:


Alkaline phosphatase in serum consists of four structural genotypes: the liver-bone-kidney type, the intestinal type, the placental type and the variant from germ cells. It occurs in osteoblasts, hepatocytes the kidneys, spleen, placenta, prostate, leukocytes and the small intestine. The liver-bone-kidney type is particularly important.

A rise in the alkaline phosphatase activity occurs with all forms of cholestasis, particularly with obstructive jaundice. It is also elevated in diseases of the skeletal system, such as Paget's disease, hyperparathyroidism, rickets and osteomalacia, as well as with fractures and malignant tumors. A considerable rise in the alkaline phosphatase activity is sometimes seen in children and juveniles. It is caused by increased osteoblast activity following accelerated bone growth. Various reference values for the purposes of clinical evaluation have been assigned to differing age groups.

In 1946, Bessey, Lowry and Brock published a method for the determination of alkaline phosphatase using p-nitrophenyl phosphate as substrate buffered with glycine/NaOH. In 1967, Hausamen et al improved upon the method by using diethanolamine as buffer. The "optimized standard method" by using diethanolamine as buffer. The assay described here meets the recommendations of the IFCC.

Test principle:

Colorimetric assay in accordance with a standardized method.

In the presence of magnesium and zinc ions, p-nitrophenylphosphate is hydrolyzed by phosphatases to form phosphate and p-nitrophenol. In this process AMP serves as transient phosphate acceptor. The release of coloured p-nitrophenol is proportional to the ALP activity and can be measured photometrically.

Reagent concentration:

R1:	
AMP* buffer, pH 10.44 (30°C)	1.0 mol/l
Magnesiumacetate	2.0 mmol/l
Zinculfate	0.5 mmol/l
EDTA	2.0 mmol/l
Detergents and stabilizers	<0.1%
R2:	
p - Nitrophenylphosphate	70.0 mmol/l
Stabilized liquid	
*2-Amino-2-methyl-1-propanol	

Preparation and stability:

Sample start:

5 Parts of R1 are mixed with one part R2. The resulting working reagent is stable:

60 days	at +2°C to +8°C
4 days	at +20°C to +25°C

Substrate start/automated analyzer:

R1: Ready for use.
R2: Ready for use.
Unopened kit components: Up to the expiration date at +2°C to +8°C
Onboard stability: R1: 42 days
R2: 42 days

Specimen:

Collect serum using standard sampling tubes.

Heparinized plasma.

Stability:	2 days	at +20°C to +25°C
	7 days	at +4°C to + 8°C
	4 weeks	at -20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 30 (approximate conjugated and unconjugated bilirubin: 30 mg/dl).

Hemolysis: No significant interference up to an index H of 200 (approximate hemoglobin concentration: 200 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 1000 (approximate triglycerides concentration: 2000 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below
• 0.9% NaCl

Manual procedure for sample start:			
Wavelength:	Hg 405 nm (400-420nm)		
Temperature:	+25°C / +30°C / +37°C		
Cuvette:	1 cm light path		
Zero adjustment:	air or distilled water		
	Macro	Semi	Micro
Working reagent	2500 µl	1000 µl	500 µl
Sample	50 µl	20 µl	10 µl
Mix and wait 30 seconds. Start stopwatch and read adsorption exactly after 1, 2, and 3 min. Record mean increase and calculate ΔA/min.			
Calculation for sample start:			
Valid with 1 cm cuvette, Hg 405 nm:			
ΔA/min x 2757 = Activity (U/l)			

Manual procedure for substrate start:			
Wavelength:	Hg 405 nm (400-420nm)		
Temperature:	+25°C / +30°C / +37°C		
Cuvette:	1 cm light path		
Zero adjustment:	air or distilled water		
	Macro	Semi	Micro
R1	2500 µl	1000 µl	500 µl
Sample	50 µl	20 µl	10 µl
R2	500 µl	200 µl	100 µl
Mix, read initial absorbance and start stopwatch simultaneously. Read again after exactly 1, 2 and 3 minutes and calculate ΔA/min.			
Calculation for substrate start:			
Formula valid with 1 cm cuvette, Hg 405 nm:			
ΔA/min x 3298 = Activity (U/l)			

Measuring /reportable range:

Up to 1200 U/l (20 µkat/l)

Determine samples having higher activities via the rerun function. On instruments without rerun function, dilute the samples with 0.9% NaCl or distilled/deionized water manually (e.g. 1 + 4). Multiply the result by the appropriate dilution factor (e.g. factor 5).

Conversion factor: U/l x 0.0167 = µkat/l

Expected values:

		+37°C
U/l	Men	40-129
	Women	35-104
µkat/l	Men	0.67-2.15
	Women	0.58-1.74

For reference values for children, please refer to „Pediatric reference ranges“ 3rd Edition. S.J. Soldin, C. Bruignara, I.M. Hicks, AACCC Press.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the alkaline phosphatase results should always be assessed in conjunction with the patients' medical history, clinical examination and other findings.

Fluitest® ALP IFCC

ALKALINE PHOSPHATASE

Analytical sensitivity (lower detection limit)

Detection limit: 2 U/l or 0.03 µkat/l

The lower detection limit represents the lowest measurable ALP concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using human samples and controls within run (n = 20). The following results were obtained:

Sample	Within run		
	Mean U/l	SD U/l	CV %
Sample 1	82.5	1.08	1.31
Sample 2	136	1.68	1.24
Sample 3	225	1.45	0.64

Reproducibility was determined using human samples and controls in an internal protocol within (between day n = 20). The following results were obtained:

Sample	Between day		
	Mean U/l	SD U/l	CV %
Sample 1	84.5	1.37	1.62
Sample 2	197.9	2.59	1.31
Sample 3	242.5	3.68	1.52

Method comparison:

A comparison of the Analyticon Fluitest® ALP IFCC (y) against the corresponding DGKC method (x) with 62 samples gave the following result:

$$y = 0.255x + 0.124; \quad r = 0.996$$

Quality Control:

Human Control Serum

Contronorm® Plus 5 x 5 ml #1205
 20 x 5 ml #1220

Contropath® Plus 5 x 5 ml #1305
 20 x 5 ml #1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl

S2: Bio Cal® E 10 x 3 ml #1430

Calibration frequency

Two-point calibration is recommended

- after reagent lot change
- as required following quality control procedures

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