### Order information:

| Catalog No. | Contents |     |       |    |     |       |
|-------------|----------|-----|-------|----|-----|-------|
| 1622        | R1       | 8 x | 15 ml | R2 | 1 x | 25 ml |
| 1625        | R1       | 4 x | 50 ml | R2 | 4 x | 10 ml |

### Intended use:

In vitro test for the quantitative determination of alkaline phosphatase (ALP) in human serum and plasma.

### Summary:


Alkaline phosphatase in serum consists of four structural genotypes: the liver-bone-kidney type, the intestinal type, the placental type and the variant from germ cells. It occurs in osteoblasts, hepatocytes the kidneys, spleen, placenta, prostate, leukocytes and the small intestine. The liver-bone-kidney type is particularly important.

A rise in the alkaline phosphatase activity occurs with all forms of cholestasis, particularly with obstructive jaundice. It is also elevated in diseases of the skeletal system, such as Paget's disease, hyperparathyroidism, rickets and osteomalacia, as well as with fractures and malignant tumours. A considerable rise in the alkaline phosphatase activity is sometimes seen in children and juveniles. It is caused by increased osteoblast activity following accelerated bone growth. Various reference values for the purposes of clinical evaluation have been assigned to differing age groups.

In 1946, Bessey, Lowry and Brock published a method for the determination of alkaline phosphatase using p-nitrophenyl phosphate as substrate buffered with glycine/NaOH. In 1967, Hausamen et al improved upon the method by using diethanolamine as buffer. The "optimized standard method" assay described here meets the 1972 recommendations of the "Deutschen Gesellschaft für Klinische Chemie (German Society of Clinical Chemistry)".

### Test principle:

Colorimetric assay in accordance with a standardized method.

In the presence of magnesium and zinc ions, p-nitrophenyl phosphate is hydrolyzed by phosphatases to form phosphate and p-nitrophenol. The p-nitrophenol released is proportional to the ALP activity and can be measured photometrically.

### Reagent Concentration:

|   |            |
|---|------------|
| <b>R1:</b>                                      |            |
| Diethanolamine buffer, pH 9.8                   | 1.0 mol/l  |
| Magnesiumsulfate                                | 0.6 mmol/l |
| Detergents and stabilizers                      | >0.1%      |
| <b>R2:</b>                                      |            |
| p - nitrophenylphosphate<br>(Stabilized liquid) | 2.0 mmol/l |

### Preparation and stability:

#### Sample start:

5 Parts of R1 are mixed with one part R2. The resulting working reagent is stable:

|        |                   |
|--------|-------------------|
| 6 days | at +2°C to +8°C   |
| 4 days | at +20°C to +25°C |

#### Substrate start:

R1: Ready for use.

R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C

Onboard stability: R1: 28 days

R2: 28 days

### Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

### Analytical sensitivity (lower detection limit)

Detection limit: 5 U/l or 0.08 µkat/l

The lower detection limit represents the lowest measurable ALP concentration that can be distinguished from zero.

### Specimen:

Collect serum using standard sampling tubes.

Heparinized plasma.

|            |         |                 |
|------------|---------|-----------------|
| Stability: | 2 days  | at +20 to +25°C |
|            | 7 days  | at +4 to +8°C   |
|            | 4 weeks | at -20°C        |

Centrifuge samples containing precipitate before performing the assay.

### Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 45 (approximate conjugated and unconjugated bilirubin: 45 mg/dl).

Hemolysis: No significant interference up to an index H of 1000 (approximate hemoglobin concentration: 1000 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 1000 (approximate triglycerides concentration: 2000 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

The effect of drugs or their metabolic products on the test is not known in all cases.

In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

#### Materials provided

• Working solutions as described above

#### Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

#### Manual procedure for sample start:

|                  |                        |  |  |
|------------------|------------------------|--|--|
| Wavelength:      | Hg 405 nm (400-420nm)  |  |  |
| Temperature:     | +25°C / +30°C / +37°C  |  |  |
| Cuvette:         | 1 cm light path        |  |  |
| Zero adjustment: | air or distilled water |  |  |

|                 | Macro   | Semi    | Micro  |
|-----------------|---------|---------|--------|
| Working reagent | 2500 µl | 1000 µl | 500 µl |
| Sample          | 50 µl   | 20 µl   | 10 µl  |

Mix and wait 30 seconds. Read absorbance and start stopwatch simultaneously. Read again absorbance after exactly 1, 2 and 3 minutes. Calculate ΔA/min.

#### Calculation for sample start:

Valid with 1 cm cuvette, Hg 405 nm:  
 $\Delta A/\text{min} \times 2757 = \text{Activity (U/l)}$

#### Manual procedure for substrate start:

|                  |                        |  |  |
|------------------|------------------------|--|--|
| Wavelength:      | Hg 405 nm (400-420nm)  |  |  |
| Temperature:     | +25°C / +30°C / +37°C  |  |  |
| Cuvette:         | 1 cm light path        |  |  |
| Zero adjustment: | air or distilled water |  |  |

|        | Macro   | Semi    | Micro  |
|--------|---------|---------|--------|
| R1     | 2500 µl | 1000 µl | 500 µl |
| Sample | 50 µl   | 20 µl   | 10 µl  |
| R2     | 500 µl  | 200 µl  | 100 µl |

Mix, read initial absorbance and start stopwatch simultaneously. Read again after exactly 1, 2 and 3 minutes and calculate ΔA/min.

#### Calculation for substrate start:

Formula valid with 1 cm cuvette, Hg 405 nm:  
 $\Delta A/\text{min} \times 3298 = \text{Activity (U/l)}$

### Measuring /reportable range:

Manual up to 687 U/l (11.5 µkat/l) at 37°C

Hitachi up to 2000 U/l (33.3 µkat/l)

Determine samples having higher activities via the rerun function. On instruments without rerun function, manually dilute the samples with 0.9% NaCl or distilled/deionized water (e.g. 1 + 4). Multiply the result by the appropriate dilution factor (e.g. factor 5).

### Expected values:

#### Reference range according to Rosalki

|        |       | 25°C  | 30°C  | 37°C  |
|--------|-------|-------|-------|-------|
| U/l    | Men   | <180  | <220  | <270  |
|        | Women | <160  | <195  | <240  |
| µkat/l | Men   | <3.00 | <3.67 | <4.50 |
|        | Women | <2.67 | <3.25 | <4.00 |

For reference values for children, please refer to „Pediatric reference ranges“ 3. Edition. S.J. Soldin, C. Brugnara, I.M. Hicks AACC Press.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the alkaline phosphatase results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

# Fluitest<sup>®</sup> ALP DGKC

ALKALINE PHOSPHATASE

## **Imprecision:**

Reproducibility was determined using controls in an internal protocol within (between day n = 20). The following results were obtained:

| Sample   | Between Day |           |      |
|----------|-------------|-----------|------|
|          | Mean<br>U/l | SD<br>U/l | % CV |
| Sample 1 | 143.1       | 2.3       | 1.59 |
| Sample 2 | 448.4       | 6.6       | 1.46 |
| Sample 3 | 449.2       | 5.9       | 1.32 |

## **Method comparison:**

A comparison of the Analyticon Fluitest<sup>®</sup> ALP DGKC (y) with a commercial obtainable assay (x) gave with 80 samples the following result:

$$y = 0.992 x + 2.64; \quad r = 0.991$$

## **Quality Control:**

Human Control Serum

|                              |           |       |
|------------------------------|-----------|-------|
| Contronorm <sup>®</sup> Plus | 5 x 5 ml  | #1205 |
|                              | 20 x 5 ml | #1220 |
| Contropath <sup>®</sup> Plus | 5 x 5 ml  | #1305 |
|                              | 20 x 5 ml | #1320 |

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

## **Calibration:**

S1: 0.9% NaCl

S2: Bio Cal<sup>®</sup> E                      10 x 3 ml                      #1430

## **Calibration frequency:**

Two-point calibration is recommended

- after reagent lot change
- as required following quality control procedures

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