### Order information:

Catalog No.	Contents
11439	R1 9 x 20 ml

### Intended use:

Enzymatic in vitro test for the quantitative determination of α-amylase in human serum, plasma and urine.

### Summary:

The α-amylases (1,4-α-D-glucanohydrolases, EC 3.2.1.1) catalyze the hydrolytic degradation of polymeric carbohydrates such as amylose, amylopectin and glycogen by cleaving 1,4-α-glucosidic bonds. In polysaccharides and oligosaccharides, several glycosidic bonds are hydrolyzed simultaneously. Maltotriose, the smallest such unit, is converted into maltose and glucose, albeit very slowly.


Two types of α-amylases can be distinguished, the pancreatic type (P-type) and the salivary type (S-type). Whereas the P-type can be attributed almost exclusively to the pancreas and is therefore organ-specific, the S-type can originate from a number of sites. As well as appearing in the salivary glands it can also be found in tears, sweat, human milk, amniotic fluid, the lungs, testes and the epithelium of the fallopian tube.

Because of the sparsity of specific clinical symptoms of pancreatic diseases, α-amylase determinations are of considerable importance in pancreatic diagnostics. They are mainly used in the diagnosis and monitoring of acute pancreatitis. Hyperamylasemia does not, however, only occur with acute pancreatitis or in the inflammatory phase of chronic pancreatitis, but also in renal failure (reduced glomerular filtration), tumors of the lungs or ovaries, pulmonary inflammation, diseases of the salivary gland, diabetic ketoacidosis, cerebral trauma, surgical interventions or in the case of macroamylasemia. To confirm pancreatic specificity, it is recommended that an additional pancreas-specific enzyme - lipase or pancreatic-α-amylase - also be determined.

Numerous methods have been described for the determination of α-amylase. These either determine the decrease in the amount of substrate viscometrically, turbidimetrically, nephelometrically and amyloclastically or measure the formation of degradation products saccharogenically or kinetically with the aid of enzyme-catalyzed subsequent reactions. The kinetic method described here is based on the cleavage of 2-chloro-4-nitrophenyl-αD-maltotrioside (CNP-G3) by α-amylase.

### Test principle:

Colorimetric test with 2-chloro-4-nitrophenyl-αD-maltotrioside (CNP-G3) as direct substrate. Colour is released directly as a result of a cleavage at the α-glycone:

(CNP = chloro-nitrophenol; G = glucose)

The increase of absorption of chloro-nitrophenol is directly proportional to the α-amylase concentration. The hydrolysis pattern in the formulation of the reagent show about less than 10 % CNP-G2 and less than 1% CNP-G4 as by products.

### Reagent concentration:

<b>R1:</b>	
MES buffer, pH 6.0	100 mmol/l
NaCl	350 mmol/l
Ca-Acetate	6 mmol/l
Potassium thiocyanate	900 mmol/l
CNP-G3	2.27 mmol/l
Stabilizers and detergents	> 0.1 %

### Preparation and stability:

R1: Ready for use	
The reagent is stable:	
Unopened:	up to expiry date at +2°C to +8°C
Opened:	14 days at +20°C to +25°C
	28 days at +2°C to +8°C.

### Specimen:

#### Serum/plasma

Collect serum using standard sampling tubes.

Heparinized or EDTA- plasma.	
Stability:	7 days at +20°C -25°C
	1 month at +2°C - 8°C

#### Urine

Collect without additives.

Stability:	2 days at +20°C - 25°C
	10 days at +2°C - +8°C

α-Amylase is unstable in acid urine. Assay promptly or adjust pH to alkaline range (about pH 7) before storage.

Centrifuge samples containing precipitate before performing the assay.

### Notes:

For in vitro diagnostic use. The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

A slight change in the yellow coloration of solution 2 does not interfere with the performance of the test.

Do not pipette by mouth, and ensure that the reagent does not come into contact with the skin. **(Saliva and sweat contain α-amylase!)**

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to 70 mg/dl bilirubin.

Hemolysis: No significant interference up to 170 mg/dl hemoglobin.

Lipemia (Intralipid): No significant interference up to 2600 mg/dl triglycerides.

An increase of the initial absorbance of the reagent to A > 0.3 (405 nm) indicates a contamination of the reagent.

### Testing procedure:

Applications for automated systems are available on request.

#### Materials provided

• Working solutions as described above

#### Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

<b>Manual procedure:</b>		
Wavelength:	Hg 405 nm (400-420nm)	
Temperature:	+37°C	
Cuvette:	1 cm light path	
Zero adjustment:	against air	
	Serum/plasma	Urine
R1	1000 µl	1000 µl
Serum/plasma	20 µl	---
Urine	---	10 µl
Mix and incubate 1 min at +37°C. Then read initial absorbance and start stopwatch simultaneously. Read again after exactly 1, 2 and 3 minutes. Determine the mean change of absorbance per minute (ΔA/min) and use this for the calculation.		
<b>Calculation:</b>		
Use absorption differences to calculate ΔA/min. Multiply with the following factors:		
	Serum/plasma	Urine
Activity; +37°C (U/l)	10256 x ΔA/min.	20118 x ΔA/min

### Measuring /reportable range:

Measuring range: Up to 1500 U/l (25.8 µka/l)

Dilute samples having higher activities with 0.9% NaCl (e.g. 1 + 9). Multiply the result by the appropriate dilution factor (e.g. factor 10).

### Expected values:

	+37°C
Serum/plasma	< 220 U/l
Spontaneously voided urine	< 1000 U/l
24 h urine	< 900 U/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the α-amylase results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

### Analytical sensitivity (lower detection limit)

Detection limit: 7 U/l or 0.12 µkat/l

### Imprecision:

Reproducibility within run was determined using human samples and controls (n = 20). The following results were obtained:

Serum	Within run		
	Mean (U/l)	SD (U/l)	CV %
Sample 1	187	1.48	0.79
Sample 2	446	2.10	0.47
Sample 3	507	2.64	0.52

Reproducibility was determined using human samples and controls between day (n = 20). The following results were obtained:

Serum	Between day		
	Mean (U/l)	SD (U/l)	CV %
Sample 1	196	2.83	1.44
Sample 2	474	6.25	1.32
Sample 3	542	6.73	1.24

# Fluitest<sup>®</sup> AMYL CNPG3

α-AMYLASE

## **Method comparison:**

A comparison of the Analyticon Fluitest<sup>®</sup> AMYL (y) with a commercial obtainable assay (x) gave with 36 samples the following result:

$$y = 0.972x + 1.282; r = 0.999$$

## **Quality Control:**

Human Control Serum

Contronorm <sup>®</sup> Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath <sup>®</sup> Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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