### Order information:

Catalog-No.	Contents					
11732	R1	8 x	15 ml	R2	1 x	25 ml
H2001 Hit I (Ilab*)	R1	6 x	47 ml	R2	6 x	11 ml
H2003 Hit 917 (AU*)	R1	6 x	60 ml	R2	6 x	14 ml
AU2003 AU	R1	6 x	60 ml	R2	6 x	14 ml

(\* ) Kit contains only reagent barcode for Hitachi system.

### System information:

Hitachi 911: ACN 028

Hitachi 917: ACN 570

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

### Intended use:

Enzymatic in vitro test for the quantitative determination of α-amylase in human serum, plasma and urine.

### Summary:

The α-amylases (1,4-α-D-glucanohydrolases, EC 3.2.1.1) catalyze the hydrolytic degradation of polymeric carbohydrates such as amylose, amylopectin and glycogen by cleaving 1,4-α-glucosidic bonds. In polysaccharides and oligosaccharides, several glycosidic bonds are hydrolyzed simultaneously. Maltotriose, the smallest such unit, is converted into maltose and glucose, albeit very slowly.

Two types of α-amylases can be distinguished, the pancreatic type (P-type) and the salivary type (S-type). Whereas the P-type can be attributed almost exclusively to the pancreas and is therefore organ-specific, the S-type can originate from a number of sites. As well as appearing in the salivary glands it can also be found in tears, sweat, human milk, amniotic fluid, the lungs, testes and the epithelium of the fallopian tube.



Because of the sparsity of specific clinical symptoms of pancreatic diseases, α-amylase determinations are of considerable importance in pancreatic diagnostics. They are mainly used in the diagnosis and monitoring of acute pancreatitis. Hyperamylasemia does not only occur with acute pancreatitis or in the inflammatory phase of chronic pancreatitis, but also in renal failure (reduced glomerular filtration), tumors of the lungs or ovaries, pulmonary inflammation, diseases of the salivary gland, diabetic ketoacidosis, cerebral trauma, surgical interventions or in the case of macroamylasemia. To confirm pancreatic specificity, it is recommended that an additional pancreas specific enzyme - lipase or pancreatic-α-amylase- also be determined.

Numerous methods have been described for the determination of α-amylase. These either determine the decrease in the amount of substrate viscometrically, turbidimetrically, nephelometrically and amyloclastically or measure the formation of degradation products saccharogenically or kinetically with the aid of enzyme-catalyzed subsequent reactions. The kinetic method described here is based on the well-proven cleavage of 4,6-ethylidene-(G<sub>7</sub>)-1,4-nitrophenyl (G<sub>1</sub>)-D-maltoheptaoside (Ethylidene Protected Substrate = EPS) by α-amylase and subsequent hydrolysis of all the degradation products to p-nitrophenol with the aid of α-glucosidase (100% chromophore liberation). The results of this method correlate with those obtained by HPLC.

### Test principle:

Enzymatic colorimetric assay according to the IFCC-method.

Defined oligosaccharides such as 4,6-ethylidene-(G<sub>7</sub>) p-nitro phenyl-(G<sub>1</sub>)-α-D-maltoheptaoside (ethylidene-G<sub>7</sub>PNP) are cleaved under the catalytic action of α-amylases. The G<sub>2</sub>PNP, G<sub>3</sub>PNP and G<sub>4</sub>PNP fragments so formed are completely hydrolyzed to p-nitrophenol and glucose by α-glucosidase.

(PNP = p-nitrophenol; G = glucose)

The color intensity of the p-nitrophenol formed is directly proportional to the α-amylase activity and is measured photometrically.

### Reagent Concentration:

#### R1:

Hepes* buffer, pH 7.15	52.5 mmol/l
NaCl	87 mmol/l
MgCl <sub>2</sub>	12.6 mmol/l
CaCl <sub>2</sub>	0.075 mmol/l
α-Glucosidase mod.	> 8 KU/l
Preservative	

#### R2:

Hepes* buffer, pH 7.15	52.5 mmol/l
4,6-ethylidene-G <sub>7</sub> PNP	22 mmol/l

\*Hepes = 2[4-(2-hydroxyethyl)-1-piperazinyl]-ethanesulfonic acid

### Preparation and stability:

#### Substrate start/Hitachi:

R1: Ready for use

R2: Ready for use

Stability: Unopened kit components: Up to the expiration date at 2-8°C

Onboard Stability: R1: 28 days

R2: 28 days

#### Sample start:

5 Parts of R1 are mixed with one part R2. The resulting working reagent is stable:

30 days at +2°C to +8°C

10 days at +20°C to +25°C

#### Specimen:

Collect serum using standard sampling tubes.

Heparinized or EDTA- plasma.

Stability: 7 days at 20-25°C

1 month at 2-8°C

#### Urine

Collect without additives.

Stability: 2 days at 20-25°C

10 days at 2-8°C

α-Amylase is unstable in acid urine. Assay promptly or adjust pH to alkaline range (about pH 7) before storage.

Centrifuge samples containing precipitate before performing the assay.

#### Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

#### Limitations - interference:

A slight change in the yellow coloration of solution 2 does not interfere with the performance of the test.

Do not pipette by mouth, and ensure that the reagent does not come into contact with the skin. (Saliva and sweat contain α-amylase!)

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 60

(approximate unconjugated bilirubin concentration: 60 mg/dl)

Hemolysis: No significant interference up to an index H of 500 (approximate haemoglobin concentration: 500 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 1000 (approximate triglycerides concentration: 2000 mg/dl). There is poor correlation between turbidity and triglycerides concentration

Highly turbid and grossly lipemic samples may cause Abs. flags. Glucose: No interference from glucose up to 2000 mg/dl. Approximately 10% lower recovery was found at glucose concentrations of 4500 mg/dl.

Ascorbic acid: No interference from ascorbic acid up to 100 mg/dl. Approximately 10% lower recovery was found at ascorbic acid concentrations of 880 mg/dl.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

#### Testing procedure:

Applications for automated systems are available on request.

#### Materials provided

• Working solutions as described above

#### Additional materials required

• Calibrators and controls as indicated below  
• 0.9% NaCl

#### Manual procedure for substrate start:

Wavelength:	405 nm (400-420nm)
Temperature:	37°C
Cuvette:	1 cm
Zero adjustment:	air or distilled water

	Serum / Plasma	Urine
R1	500 µl	500 µl
Sample	20 µl	10 µl
R2	100 µl	100 µl

Mix, read initial absorbance and start stopwatch simultaneously. Read again after exactly 1, 2 and 3 minutes.

#### Calculation:

##### Serum:

ΔA/min x 3180 = activity (U/l)

##### Urine:

ΔA/min x 6255 = activity (U/l)

### Manual Procedure for sample start:

Wavelength:	405 nm (400-420nm)
Temperature:	37°C
Cuvette:	1 cm
Zero adjustment:	air or distilled water

	Serum / Plasma	Urine
Working solution	1000 µl	1000 µl
Sample	20 µl	10 µl

Mix, read initial absorbance and start stopwatch simultaneously. Read again after exactly 1, 2 and 3 minutes.

### Calculation:

**Serum:**  
 $\Delta A/\text{min} \times 5230 = \text{activity (U/l)}$

**Urine:**  
 $\Delta A/\text{min} \times 10357 = \text{activity (U/l)}$

### Measuring /reportable range:

Measuring range: 3-1500 U/l (0.05-25.00 µkat/l)

Determine samples having higher activities via the rerun function. On instruments without rerun function, manually dilute the samples with 0.9% NaCl or distilled/deionized water (e.g. 1 + 4). Multiply the result by the appropriate dilution factor (e.g. factor 5).

### Expected values:

	U/l	µkat/l
Serum/plasma	28 - 100	0.47 – 1.67
Spontaneously voided urine	≤ 460	≤ 7.67
α-amylase/Creatinin Quotient	≤ 310 U/g	≤ 5.17 µkat/g

EDTA plasma values are approximately 8% lower than serum values.

**Note:** The results (U/l or µkat/l) can be approximated to the reference range for α-amylase EPS by multiplying by the factor 2.2. This conversion factor only applies to human samples, not control sera.

### α-Amylase/ creatinine quotient

To allow for fluctuations in the α-amylase activity in urine, it is advisable to determine the α-amylase/creatinine quotient. To do this, determine the α-amylase activity and creatinine concentration in spontaneously voided urine.

$$\text{Quotient [U/g] or } \mu\text{kat/mmol} = \frac{\alpha\text{-amylase (U/l or } \mu\text{kat/l)}}{\text{creatinine [g/l or mmol/l]}}$$

### Amylase/Creatinine Clearance Ratio (ACCR)

The ACCR is calculated from amylase activity and creatinine concentration. Both the serum and urine samples should be collected at the same time.

$$\text{ACCR [\%]} = \frac{\text{Urine amylase [U/l]} \times \text{serum creatinine [mg/l]}}{\text{Serum amylase [U/l]} \times \text{urine creatinine [mg/l]}} \times 100$$

ACCR approximately equal to 2-5%

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the α-amylase results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

### Analytical sensitivity (lower detection limit)

Detection limit: 3 U/l or 0.05 µkat/l

The lower detection limit represents the lowest α-amylase activity that can be distinguished from zero

### Method comparison:

A comparison of the Analyticon Fluitest® AMYL-IFCC (y) with a commercial obtainable assay (x) gave the following result:

$$y = 1.058 + 3.482x; \quad r = 1.000$$

### Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol (n = 20). The following results were obtained:

Probe	Within run			Between day		
	MW U/l	SD U/l	VK %	MW U/l	SD U/l	VK %
Probe 1	79.4	1.20	1.51	79.7	0.91	1.15
Probe 2	182	3.39	1.86	186	1.53	0.82
Probe 3	198	3.10	1.56	203	1.56	0.77

### Quality Control:

Human Control Serum:		
Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
ControPath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320
Urine Control:		
Urine Control Set	8 x 5 ml	#1507

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### Calibration:

S1: 0.9% NaCl		
S2: Bio Cal® E	10 x 3 ml	#1430

### Calibration frequency:

Two point calibration is recommended:

- after lot change
- as required following quality control procedures

Calibration verification: Not necessary

### Disposal:

Please note the legal regulations.

### Literature:

1. Greiling H, Gressner AM eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3<sup>rd</sup> ed. Stuttgart/New York: Schattauer Verlag, 1995.
2. Junge W, Troge B, Klein G et al. Evaluation of a New Assay for Pancreatic Amylase: Performance Characteristics and Estimation of Reference Intervals. Clin Biochem 1989;22:109-114
3. Keller H ed. Klinisch-chemische Labordiagnostik für die Praxis, 2<sup>nd</sup> ed. Stuttgart/New York: Georg Thieme Verlag, 1991;354-361.
4. Kurse-Jarres JD, Hafkenschied JCM, Hohenwallner W et al. Evaluation of a New α-Amylase Assay Using 4.6-Ethylidene-(G<sub>7</sub>)-1-4-nitrophenyl-(G<sub>1</sub>)-α-D-maltoheptaoside as Substrate. J Clin Chem Clin Biochem 1989;27:103-113.
5. Kurrle-Weitenhiller A, Hölzel W, Engel D et al. Method for the determination of total and pancreatic α-amylase based on 100% Cleavage of the protected substrate ethylidene-4-nitrophenyl-maltoheptaoside. Clin Chem 1996; 42:98
6. Lott J. Inflammatory diseases of the pancreas CRC Critical Reviews in Clinical Laboratory Sciences 1982;17:201--228.
7. Rizzotti P, Klein G. Evaluation of a Specific immunoinhibition Method for the Determination of Pancreatic α-Amylase. Eur J Clin Chem Clin Biochem 1994;32:97-106
8. Tietz NW, Burlina A, Gerhardt W et al. Multicenter Evaluation of a Specific Pancreatic Isoamylase Assay Based on a Double Monoclonal Antibody Technique. Clin Chem 1988;34:2096-2102.
9. Tietz NW, Huang WY, Rauh DF et al. Laboratory tests in the differential diagnosis of hyperamylasemia. Clin Chem 1986; 32: 301-307.
10. Gerber
11. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790

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