Order information:

Catalog No.	Contents			
11333	R1	1 x	65 ml	R2 20 x for 3 ml
11330	R1	1 x	110 ml	R2 10 x for 10 ml
11332	R1	1 x	110 ml	R2 5 x for 20 ml

Intended use:

Enzymatic in vitro test for the quantitative determination of α-amylase in human serum, plasma and urine.

Summary:

The α-amylases (1,4-α-D-glucanohydrolases, EC 3.2.1.1) catalyze the hydrolytic degradation of polymeric carbohydrates such as amylose, amylopectin and glycogen by cleaving 1,4-α-glucosidic bonds. In polysaccharides and oligosaccharides, several glycosidic bonds are hydrolyzed simultaneously. Maltotriose, the smallest such unit, is converted into maltose and glucose, albeit very slowly.

Two types of α-amylases can be distinguished, the pancreatic type (P-type) and the salivary type (S-type). Whereas the P-type can be attributed almost exclusively to the pancreas and is therefore organ-specific, the S-type can originate from a number of sites. As well as appearing in the salivary glands it can also be found in tears, sweat, human milk, amniotic fluid, the lungs, testes and the epithelium of the fallopian tube.



Because of the sparsity of specific clinical symptoms of pancreatic diseases, α-amylase determinations are of considerable importance in pancreatic diagnostics. They are mainly used in the diagnosis and monitoring of acute pancreatitis. Hyperamylasemia does not, however, only occur with acute pancreatitis or in the inflammatory phase of chronic pancreatitis, but also in renal failure (reduced glomerular filtration), tumors of the lungs or ovaries, pulmonary inflammation, diseases of the salivary gland, diabetic ketoacidosis, cerebral trauma, surgical interventions or in the case of macroamylasemia. To confirm pancreatic specificity, it is recommended that an additional pancreas-specific enzyme - lipase or pancreatic-α-amylase- also be determined.

Numerous methods have been described for the determination of α-amylase. These either determine the decrease in the amount of substrate viscometrically, turbidimetrically, nephelometrically and amyloclastically or measure the formation of degradation products saccharogenically or kinetically with the aid of enzyme-catalyzed subsequent reactions. The kinetic method described here is based on the well-proven cleavage of 4,6-ethylidene-(G₇)-1,4-nitrophenyl-(G₁)-, α₁D-maltoheptaoside (Ethylidene Protected Substrate = EPS) by α-amylase and subsequent hydrolysis of all the degradation products to p-nitrophenol with the aid of α-glucosidase (100% chromophore liberation).

Test principle:

Enzymatic colorimetric assay. The assay meets the recommendations of the "Deutschen Gesellschaft für Klinische Chemie (German Society of Clinical Chemistry)".

Defined oligosaccharides such as 4,6-ethylidene-(G₇) p-nitrophenyl-(G₁)-α₁D-maltoheptaoside (ethylidene-G₇PNP) are cleaved under the catalytic action of α-amylases. The G₂PNP G₃PNP and G₄PNP fragments so formed are completely hydrolyzed to p-nitrophenol and glucose by α-glucosidase.

(PNP = p-nitrophenol; G = glucose)

The color intensity of the p-nitrophenol formed is directly proportional to the α-amylase activity and is measured photometrically.

Reagent Concentration:

R1:	
Hepes* buffer, pH 7,1	100 mmol/l
NaCl	50 mmol/l
MgCl ₂	10 mmol/l
R2:	
α-Glucosidase	>23 U/l
4.6-Ethyliden-G ₇ PNP	3 mmol/l

*Hepes = 2[4-(2-hydroxyethyl)-1-piperazinyl]-ethanesulfonic acid

Preparation and stability:

Unopened kit components: Up to the expiration date at +2°C to +8°C.
Dilute volume of enzyme/substrate reagent / R2 with the corresponding volume of buffer /R1.

This working reagent is stable:

5 days	at +20°C to +25°C
28 days	at +2°C to +8°C.

Specimen:

Collect serum using standard sampling tubes.

Heparinized or EDTA- plasma.

Stability: 7 days at +20°C to +25°C
1 month at +2°C to +8°C

Urine Collect without additives.

Stability: 2 days at +20°C to +25°C
10 days at +2°C to +8°C

α-Amylase is unstable in acid urine. Assay promptly or adjust pH to alkaline range (about pH 7) before storage.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

A slight change in the yellow coloration of solution 2 does not interfere with the performance of the test.

Do not pipette by mouth, and ensure that the reagent does not come into contact with the skin. (**Saliva and sweat** contain α-amylase!)

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 28

(approximate bilirubin concentration: 28 mg/dl)

Hemolysis: No significant interference up to an index H of 270 (approximate hemoglobin concentration: 270 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 270 (approximate triglycerides concentration: 540 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Manual Procedure:				
Wavelength:	Hg 405 nm (400-420nm)			
Temperature:	+25 / +30°C / +37°C			
Cuvette:	1 cm			
Zero adjustment:	against air/ aqua dest.			
+25°C/+30°C	Macro	Semi-micro		
	serum/plasma	urine	serum/plasma	urine
Working solution	2500 µl	2500 µl	1000 µl	1000 µl
	sample	100 µl	50 µl	40 µl
+37°C	Working solution	2500 µl	1000 µl	1000 µl
	sample	50 µl	20 µl	20 µl
Mix and incubate 3 minutes at +25°C to + 30°C or 2 minutes at +37°C. Then read absorbance and start stopwatch at the same time. Repeat readings after exactly 1, 2 and 3 minutes. Determine the mean absorbance change per minute (ΔA/min), and use this for the calculation.				
Calculation:				
+25 /+30°C	Serum/plasma	urine		
	7893 x ΔA/min	15408 x ΔA/min		
+37°C	15408 x ΔA/min	37957 x ΔA/min		

Measuring /reportable range:

Serum/plasma: up to 2000 U/l (33.3 µkat/l)

Urine: up to 5000 U/l (83.3 µkat/l)

Dilute samples having higher activities with 0.9% NaCl or distilled/deionized water (e.g. 1+4). Multiply the result by the appropriate dilution factor (e.g. factor 5).

Expected values:

	25°C	30°C	37°C
Serum / plasma	< 120 U/l	< 160 U/l	< 220 U/l
Spontaneously voided urine	600 U/l	800 U/l	1000 U/l
24 h urine	450 U/l	650 U/l	900 U/l

α-Amylase/ creatinine quotient

To allow for fluctuations in the α-amylase activity in urine, it is advisable to determine the α-amylase/creatinine quotient. To do this, determine the α-amylase activity and creatinine concentration in spontaneously voided urine.

$$\text{Quotient [U/g or } \mu\text{kat/mmol]} = \frac{\alpha\text{-amylase (U/l or } \mu\text{kat/l)}}{\text{creatinine [g/l or mmol/l]}}$$

Amylase/Creatinine Clearance Ratio (ACCR)

The ACCR is calculated from amylase activity and creatinine concentration. Both the serum and urine samples should be collected at the same time.

$$\text{ACCR [\%]} = \frac{\text{Urine amylase [U/l]} \times \text{serum creatinine [mg/l]}}{\text{Serum amylase [U/l]} \times \text{urine creatinine [mg/l]}} \times 100$$

ACCR approximately equal to 2-5%

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the α-amylase results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 7 U/l or 0.12 μkat/l

The lower detection limit represents the lowest α-amylase activity that can be distinguished from zero

Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol (n = 20). The following results were obtained:

Between day			
Sample	Mean U/l	SD U/l	%CV
Control serum 1	191	4.30	2.24
Control serum 2	440	11.17	2.54
Control serum 3	498	14.88	2.99
Within run			
Sample	Mean U/l	SD U/l	%CV
Control serum 1	189	1.72	0.91
Control serum 2	433	2.97	0.68
Control serum 3	496	2.71	0.55

Method comparison:

A comparison of the Analyticon AMYL (y) with a commercial obtainable assay (x) gave the following result:

$$y = 1.030x - 3.331; \quad r = 1.000$$

Quality control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
ControPath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl		
S2: Bio Cal® E	10 x 3 ml	#1430

Disposal:

Please note the legal regulations.

Literature:

- Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474
- Greiling H, Gressner AM ed. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3. Auflage. Stuttgart/New York: Schattauer Verlag, 1995.
- Junge W, Troge B, Klein G et al. Evaluation of a New Assay for Pancreatic Amylase: Performance Characteristics and Estimation of Reference Intervals. Clin Biochem 1989;22:109-114
- Junge W, Waldenström J, Bouman A et al. Evaluation of the Assay for Total and Pancreatic Alpha-Amylase based on 100% Cleavage of Et-G7-PNP at 6 European Clinical Centres (Poster Medlab 97). Basel, Switzerland: 12. IFCC European Congress of Clinical Chemistry, August 17-22, 1997
- Keller H ed. Klinisch-chemische Labordiagnostik für die Praxis, 2. Auflage. Stuttgart/New York: Georg Thieme Verlag, 1991;354-361.
- Kurrle-Jarres JD, Hafkenschied JCM, Hohenwallner W et al. Evaluation of a New α-Amylase Assay Using 4,6-Ethylidene-(G₇)-1,4-nitrophenyl-(G₁)-α-D-maltoheptaoside as Substrate. J Clin Chem Clin Biochem 1989;27:103-113.
- Kurrle-Weitenhiller A, Hölzel W, Engel D et al. Method for the Determination of Total and Pancreatic α-Amylase based on 100% Cleavage on the protected substrate Ethylidene-4-nitrophenyl- Maltoheptaoside. Clin Chem 1996; 42:98
- Lorentz K. Approved Recommendation on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes. Part 9. IFCC-Method for α-Amylase (1,4-α-D-Gluca 4-Glucano- Hydrolase, EC 3.2.1.1.). Clin Chem Lab Med 1998;38:195-203
- Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720
- Rauscher E et al. Fresenius Z Analyt Chem 1986; 324:304.
- Salt WB II, Schenker S. Amylase-its clinical significance: a review of the literature [Review]. Medicine 1976; 55:269-281 .
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis [Review]. Ann Intern Med 1985; 102:576 - 580.
- Tietz NW ed: Clinical Guide to Laboratory Tests. 3. Auflage. Philadelphia, PA: W B Saunders Company; 1995:46-51
- Tietz NW, Huang WY, Rauh DF et al. Laboratory tests in the differential diagnosis of hyperamylasemia. Clin Chem 1986; 32: 301-307.

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Literatur:

- Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474
- Greiling H, Gressner AM (Hrsg.). Lehrbuch der Klinischen Chemie und Pathobiochemie, 3. Auflage. Stuttgart/New York: Schattauer Verlag, 1995.
- Junge W, Troge B, Klein G et al. Evaluation of a New Assay for Pancreatic Amylase: Performance Characteristics and Estimation of Reference Intervals. Clin Biochem 1989;22:109-114
- Junge W, Waldenström J, Bouman A et al. Evaluation of the Assay for Total and Pancreatic Alpha-Amylase based on 100% Cleavage of Et-G7-PNP at 6 European Clinical Centres (Poster Medlab 97). Basel, Switzerland: 12. IFCC European Congress of Clinical Chemistry, August 17-22, 1997
- Keller H (Hrsg.). Klinisch-chemische Labordiagnostik für die Praxis, 2. Auflage. Stuttgart/New York: Georg Thieme Verlag, 1991;354-361.
- Kruse-Jarres JD, Hafkenscheid JCM, Hohenwallner W et al. Evaluation of a New α-Amylase Assay Using 4.6-Ethylidene-(G₇)-1-4-nitrophenyl-(G₁)-α-D-maltoheptaoside as Substrate. J Clin Chem Clin Biochem 1989;27:103-113.
- Kurrle-Weitenhiller A, Hölzel W, Engel D et al. Method for the Determination of Total and Pancreatic α-Amylase based on 100% Cleavage on the protected substrate Ethylidene-4-nitrophenyl-Maltoheptaoside. Clin Chem 1996; 42:98
- Lorentz K. Approved Recommendation on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes. Part 9. IFCC-Method for α-Amylase (1,4-α-D-Gluca 4-Glucono- Hydrolase, EC 3.2.1.1.). Clin Chem Lab Med 1998;38:195-203
- Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720
- Rauscher E et al. Fresenius Z Analyt Chem 1986; 324:304.
- Salt WB II, Schenker S. Amylase-its clinical significance: a review of the literature [Review]. Medicine 1976; 55:269-281 .
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis [Review]. Ann Intern Med 1985; 102:576 - 580.
- Tietz NW (Hrsg.) : Clinical Guide to Laboratory Tests. 3. Auflage. Philadelphia, PA: W B Saunders Company; 1995:46-51
- Tietz NW, Huang WY, Rauh DF et al. Laboratory tests in the differential diagnosis of hyperamylasemia. Clin Chem 1986; 32: 301-307.