Turbitex® APO A1

APOLIPOPROTEIN A1

BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B0831	200 / 600	R1	6 x	20 ml
		R2	6 x	8 ml
B0833	300 / 600*	R1	6 x	20 ml
		R2	6 x	8 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of Apolipoprotein A-1 in human serum and plasma.

Summary:

Apolipoprotein A-1 is the main component of the high-density lipoprotein class. It accounts for approximately 30 % of the HDL-particle. Apolipoprotein A-1 originates in equal amounts from the small intestine and the liver. It has a physiological role due to its cofactor activity in the lecithin-cholesterol-acryl-transferase (LCAT) reaction and its ability to take up free cholesterol from cells. These processes are important for the reverse cholesterol transport to the liver.

As published in the scientific literature, Apolipoprotein determinations are primarily used for the early detection and assessment of the individual coronary risk and for the diagnosis of the relatively rare (but only serious) apolipoproteinopathy, which is generally inherited. The determination of apolipoprotein A-1 is used in the detection of Tangier disease, an analphalipoproteinemia (Apo A-1 low, cholesterol low). A positive angiography finding, i.e. the demonstration of stenosis, is associated with pathological Apo A-1/Apo B findings. The Apolipoprotein determinations supplement modern HDL, LDL and Lp(a) diagnostics.

The determination of apolipoprotein A-1 is made by immunoassay (RIA; ELISA), electroimmunodiffusion (EID), radial immunodiffusion, nephelometry or turbidimetry.

Test principle:

Immunoturbidimetric assay.

Anti-Apolipoprotein A-antibodies react with the antigen in the sample to form antigen/antibody complexes which, following agglutination, can be measured turbidimetrically.

Reagent concentration:

R1:

Tris buffer pH 7.8	50 mmol/l
PEG	5 %
NaCl	150 mmol/l

Preservative/ Detergent

R2:

Polyclonal anti-human APO A-1	
Antibody (goat)	dependent on titer
Tris buffer pH 7.8	150 mmol/l
NaCl	150 mmol/l
Preservative	

Preparation and stability:

R1: Ready for use.

R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C

Onboard stability at +2°C to +8°C:	R1:	42 days
	R2:	42 days

Specimen:

Collect serum using standard sampling tubes

Na-heparin, or EDTA plasma

Stability:	1 day	at +15°C - +25°C
	3 days	at + 2°C - +8°C
	2 months	at - 20°C (only freeze once)

Centrifuge samples containing precipitates before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 20 (approximate conjugated and unconjugated bilirubin concentration: 20 mg/dl).

Hemolysis: No significant interference up to an index H of 1200 (approximate haemoglobin concentration: 1200 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 212.5 (approximate triglycerides concentration: 425 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

No high-dose hook-effect up to an apolipoprotein A-1 concentration of 600 mg/dl.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

- Controls and calibrators as indicated below
- 0.9% NaCl

Expected values:

Men:	104 – 202 mg/dl (1.04 – 2.02 g/l)
Women:	108 – 225 mg/dl (1.08 – 2.25 g/l)

Each laboratory should investigate the transferability of the expected values to its own patients' population and if necessary determine its own reference range. For diagnostic purposes the APO A-1 results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

40 mg/dl (0.4 g/l)

The lower detection limit represents the lowest measurable apolipoprotein A-1 concentration that can be distinguished from zero.

Measuring:

80-750 mg/dl or 0.8-7.5 g/l

Determine samples having higher concentrations via the rerun function using 0,9 % NaCl solution as diluents.

Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol (n = 20). The following results were obtained:

Sample	Within run		
	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	303.5	10.1	3.3
Sample 2	264	8.9	3.4
Sample 3	168	5.1	3.0

Sample	Run to run		
	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	163.7	13.9	8.5
Sample 2	235	11.3	4.8
Sample 3	254	13.4	5.3

Method comparison:

A comparison of the Analyticon APO A1 (y) with a commercial obtainable assay (x) gave the following result:

$$y = 0.981 x + 4.81; r = 0.991$$

Quality Control:

Contronorm® L	5 x 2 ml	#1302
Controptath® L	5 x 2 ml	#1303

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The apolipoprotein A-1 method was standardized using the WHO/IFCC SP1-01 reference standard (WHO-IRP 1992).

Calibration Type: Logit-Log

The kalibrator has one level, the dilutions will be done automatically by the instrument.

S1 - S5: Bio Cal® L	1 x 2 ml	#1401
	5 x 2 ml	#1402

Calibration frequency:

Full calibration is recommended:

- after lot change
- as required following quality control procedures

Disposal:

Please note the legal regulations.

Literature:

1. Evans K, Mitcheson J, Laker M. Effect of Storage at 4°C and -20°C on Lipid, Lipoprotein, and Apolipoprotein Concentrations. Clin Chem 1995;41:392-396.
2. Greiling H, Gressner AM (Hrsg. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer, 1995:325-326.
3. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag, 1996.

Turbitex[®] APO A1

APOLIPOPROTEIN A1

4. Karl J, Engel WD. Determination of Apolipoprotein A 1 and B without sample dilution. Poster presented at the 57th Meeting of the European Atherosclerosis Society, Lisbon and the IX European Congress of Clinical Chemistry, Cracow, 1991
5. Marcovina SM et al. Clin Chem 1994;40:586-592.
6. Occupational Safety and Health Standards: bloodborne pathogens (29 CFR 1910 1030). Federal Register. July 1 1998;6: 267 280.
7. Richtlinie des Rates (90/679/EWG). Amtsblatt d. Europ. Gem. Nr. L374 vom 31.12.1990:1-12.

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