Order information:

Catalog-No.	Contents
H0801 Hit 1 / 917 (Ilab* / AU*)	R1 6 x 19 ml R2 6 x 5 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911/917: ACN 168
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of Apolipoprotein A-1 in human serum and plasma.

Summary:

Apolipoprotein A-1 is the main component of the high-density lipoprotein class. It accounts for approximately 30 % of the HDL-particle. Apolipoprotein A-1 originates in equal amounts from the small intestine and the liver. It has a physiological role due to its cofactor activity in the lecithin-cholesterol-acryl-transferase (LCAT) reaction and its ability to take up free cholesterol from cells. These processes are important for the reverse cholesterol transport to the liver.

As published in the scientific literature, Apolipoprotein determinations are primarily used for the early detection and assessment of the individual coronary risk and for the diagnosis of the relatively rare (but only serious) apolipoproteinopathy, which is generally inherited. The determination of apolipoprotein A-1 is used in the detection of Tangier disease, an alpha lipoproteinemia (Apo A-1 low, cholesterol low). A positive angiography finding, i.e. the demonstration of stenosis, is associated with pathological Apo A-1/Apo B findings. The Apolipoprotein determinations supplement modern HDL, LDL and Lp(a) diagnostics.

The determination of apolipoprotein A-1 is made by immunoassay (RIA; ELISA), electroimmunodiffusion (EID), radial immunodiffusion, nephelometry or turbidimetry.

Test principle:

Immunoturbidimetric assay.
• Sample and addition of R1 (buffer)
• Addition of R2 (anti-apolipoprotein A-1-antibody) and start of reaction:

Anti-Apolipoprotein A-antibodies react with the antigen in the sample to form antigen/antibody complexes which, following agglutination, can be measured turbidimetrically.

Reagent concentration:

R1:	
Tris buffer pH 7.8	50 mmol/l
PEG	5 %
NaCl	150 mmol/l
Preservative/ Detergent	
R2:	
Polyclonal anti-human APO A-1 Antibody (goat)	dependent on titer
Tris buffer pH 7.8	150 mmol/l
NaCl	150 mmol/l
Preservative	

Preparation and stability:

R1: Ready for use.
R2: Ready for use.
Unopened kit components: Up to the expiration date at +2°C to +8°C
Onboard stability at +2°C to +8°C: R1: 42 days
R2: 42 days

Specimen:

Collect serum using standard sampling tubes
Na-heparin, or EDTA plasma
Stability: 1 day at +15°C - +25°C
3 days at + 2°C - +8°C
2 months at - 20°C (only freeze once)

Centrifuge samples containing precipitates before performing the assay.

For analyzers with automatic sample predilution, the samples are prediluted with 0.9% NaCl solution by the instrument.

For analyzers without automatic sample predilution, manually predilute samples 1+20 with 0.9% NaCl solution.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.
Icterus: No significant interference up to an index I of 30 (approximate conjugated and unconjugated bilirubin concentration: 30 mg/dl).
Hemolysis: No significant interference up to an index H of 1000 (approximate haemoglobin concentration: 1000 mg/dl).
Lipemia (Intralipid): No significant interference up to an index L of 1000 (approximate triglycerides concentration: 2000 mg/dl). There is poor correlation between turbidity and triglycerides concentration.
No high-dose hook-effect up to an apolipoprotein A-1 concentration of 600 mg/dl.
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been

discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above
- Additional materials required
- Controls and calibrators as indicated below
- 0.9% NaCl

Manual procedure:

Wavelength: 376 nm
Temperature: +37°C
Cuvette: 1cm
Zero adjustment: against blank

	Blank	Sample/ Calibrator
Sample/ Calibrator	---	12 µl
R1	500 µl	500 µl

Mix, incubate 5 min. and read A₁. Then add:

R2	100 µl	100 µl
----	--------	--------

Mix, incubate 5 min. and read A₂.

Calculation:

$$\Delta A = [(A_2 - A_1) \text{ sample or calibrator}] - [(A_2 - A_1) \text{ blank}]$$

The concentration of Apo A-1 in patient sera has to be calculated from ΔA using mathematic function as logit/log or can be read from a graph using values of 4 levels of standards. For zero value is recommended to use saline solution (0.9%).

Measuring/reportable range:

Roche/Hitachi systems automatically calculate the apolipoprotein A-1 concentration of each sample.

Conversion factors: mg/dl x 0.01 = g/l
g/l x 100 = mg/dl

- Roche/Hitachi 902 analyzer:
20 - 400 mg/dl (0.20 - 4.00 g/l)

If the APO A-1 concentration of the sample is below the measuring range, manually dilute the samples 1+10 (instead of 1+20) with 0.9% NaCl solution. Multiply the result by the factor 0.523.

- Roche/Hitachi 904/911/912/917 analyzer:
20 - 400 mg/dl (0.20 - 4.00 g/l)*
Extended measuring range with rerun: 10 - 400 mg/dl (0.10 - 4.0 g/l).

Expected values:

Men: 104 - 202 mg/dl (1.04 - 2.02 g/l)
Women: 108 - 225 mg/dl (1.08 - 2.25 g/l)

Each laboratory should investigate the transferability of the expected values to its own patients' population and if necessary determine its own reference range. For diagnostic purposes the APO A-1 results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

20.0 mg/dl (0.20 g/l)
The lower detection limit represents the lowest measurable apolipoprotein A-1 concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol (n = 12). The following results were obtained:

	within run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	86.8	1.21	1.39
Sample 2	153	2.42	1.58
Sample 3	172	1.92	1.12

Method comparison:

A comparison of the Analyticon APO A1 (y) with a commercial obtainable assay (x) gave with 50 samples the following result:
 $y = 0.981 x + 4.81$; $r = 0.991$

Turbitex® APO A1

APOLIPOPROTEIN A1

Quality Control:

Contropath® L 5 x 2 ml #1303
Contronorm® L 5 x 2 ml #1302

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The apolipoprotein A-1 method was standardized using the WHO/IFCC SP1-01 reference standard (WHO-IRP 1992).

S1: 0,9% NaCl
S2 - S5: Bio Cal® L 1 x 2 ml #1401
5 x 2 ml #1402

• Hitachi 902 analyzer:

Dilute the Bio Cal L in the following manner to prepare S2 – S6:

No.	Dilution Sample + NaCl (0.9%)	Conversion factor
S2	1 + 100	0.208
S3	1 + 50	0.412
S4	1 + 20	1.000
S5	1 + 15	1.313
S6	1 + 9.5	2.000

Multiply the lot specific Bio Cal L value by the factors above to determine the standard concentrations for the 6-point calibration curve.

• Hitachi 904/911/912/917 analyzers:

Multiply the lot-specific Bio Cal L values by the factors given below to determine the standard concentrations for the 6-point calibration curve.

S2: 0.198 S5: 1.349
S3: 0.399 S6: 2.080
S4: 1.000

Calibration frequency

Full calibration is recommended:

- after lot change
- as required following quality control procedures

Disposal:

Please note the legal regulations.

Literature:

1. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
2. Evans K, Mitcheson J, Laker M. Effect of Storage at 4°C and -20°C on Lipid, Lipoprotein, and Apolipoprotein Concentrations. Clin Chem 1995;41:392-396.
3. Greiling H, Gressner AM (Hrsg. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer, 1995:325-326.
4. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag, 1996.
5. Karl J, Engel WD. Determination of Apolipoprotein A 1 and B without sample dilution. Poster presented at the 57th Meeting of the European Atherosclerosis Society, Lisbon and the IX European Congress of Clinical Chemistry, Cracow, 1991
6. Marcovina SM et al. Clin Chem 1994;40:586-592.
7. Occupational Safety and Health Standards: bloodborne pathogens (29 CFR 1910 1030). Federal Register. July 1 1998;6: 267 280.
8. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720.
9. Richtlinie des Rates (90/679/EWG). Amtsblatt d. Europ. Gem. Nr. L374 vom 31.12.1990:1-12.

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Literatur:

1. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
2. Evans K, Mitcheson J, Laker M. Effect of Storage at 4°C and -20°C on Lipid, Lipoprotein, and Apolipoprotein Concentrations. Clin Chem 1995;41:392-396.
3. Greiling H, Gressner AM (Hrsg. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3. Auflage. Stuttgart/New York: Schattauer, 1995:325-326.
4. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes Preamerical Variables. Broschüre in: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag, 1996.
5. Karl J, Engel WD. Determination of Apotipoprotein A I and B without sample dilution. Poster presented at the 57th Meeting of the European Atherosclerosis Society, Lisbon and the IX European Congress of Clinical Chemistry, Cracow, 1991
6. Occupational Safety and Health Standards: bloodborne pathogens (29 CFR 1910 1030). Federal Register. July 1 1998;6: 267 280.
7. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720.
8. Richtlinie des Rates (90/679/EWG). Amtsblatt d. Europ. Gem. Nr. L374 vom 31.12.1990:1-12.