### BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B0931	200 / 600	R1	6 x	20 ml
		R2	6 x	6 ml
B0933	300 / 600*	R1	6 x	20 ml
		R2	6 x	6 ml

\*only for instruments with a reagent tray for 20 / 70ml bottles

### Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of Apolipoprotein B in human serum and plasma.

### Summary:

Apolipoprotein B is the main component of the LDL and IDL, and it is also an important component of VLDL and chylomicrons. It occurs in two forms: apo B-100 in the VLDL/LDL and apo B-48 in the chylomicrons.

Apo B-100 originates in the liver, whereas apo B-48 (the N-terminal half of apo B-100) is mainly formed in the small intestine. Apolipoprotein B is important for the transport of cholesterol and triglycerides as well as for their receptor-mediated cellular uptake.

As published in the scientific literature, Apolipoprotein determinations are primarily used in the early detection and assessment of the individual coronary risk and for the diagnosis of the relatively rare (but very serious) apolipoproteinopathy, which is generally inherited. The determination of apolipoprotein B is used in the detection of hyperapobetalipoproteinemia (apoB elevated, LDL-cholesterol normal), which is found in almost half of patients with coronary diseases. A positive angiography finding, i.e. the demonstration of stenosis, is associated with pathological apo A-1/apo B findings. The Apolipoprotein determinations supplement modern HDL, LDL and Lp(a) diagnostics.

The determination of apolipoprotein B is made by immunoassay (RIA; ELISA), electroimmunodiffusion (EID), radial immunodiffusion, nephelometry or turbidimetry.

### Test principle:

Immunoturbidimetric assay.

Anti-apolipoprotein B-antibodies react with the antigen in the sample to form antigen/antibody complexes which following agglutination, can be measured turbidimetrically.

### Reagent concentration:

**R1:**  
Tris\* buffer pH 7.8 100 mmol/l  
PEG 4.5 %  
NaCl 150 mmol/l  
Preservative; Detergent

**R2:**  
Polyclonal anti-human apolipoprotein B dependent on titer  
Antibody (goat)  
Tris\* buffer pH 7.8 80 mmol/l  
NaCl 150 mmol/l  
Preservative

\*TRIS= Tris (hydroxymethyl)-aminomethane

### Preparation and stability

R1: Ready for use.  
R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C  
Onboard stability at +2°C to + 8°C: R1: 42 days  
R2: 42 days

### Specimen

Collect serum using standard sampling tubes  
Na-heparin, or Na-EDTA plasma

Stability: 1 day at +20°C to +25°C  
3 days at + 4°C to + 8°C  
2 months at - 20°C (only freeze once!)

Centrifuge samples containing precipitate before performing the assay.

### Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to 100 mg/dl).

Hemolysis: No significant interference up to an index H of 350 (approximate haemoglobin concentration: 350 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 1250 (approximate triglycerides concentration: 2500 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

No high-dose hook-effect up to an apolipoprotein B concentration of 1000 mg/dl.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

*Materials provided*

• Working solutions and calibrators as described above

*Additional materials required*

• Controls as indicated below

• 0.9% NaCl

### Measuring/ reportable range:

30-320 mg/dl (0.3-3.2 g/l)

Determine samples having higher concentrations via the rerun function using 0.9% NaCl solution as diluent.

### Expected values:

Men: 66 – 133 mg/dl (0.66 – 1.33 g/l)

Women: 60 – 117 mg/dl (0.60 – 1.17 g/l)

Each laboratory should investigate the transferability of the expected values to its own patients' population and if necessary determine its own reference range. For diagnostic purposes the APO B results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

### Analytical sensitivity (lower detection limit):

22 mg/dl (0.22 g/l)

The lower detection limit represents the lowest measurable apolipoprotein B concentration that can be distinguished from zero.

### Imprecision:

Reproducibility was determined within run using controls in an internal protocol (n = 20). The following results were obtained:

<b>Within run</b>			
Sample	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	136.3	2.7	2.0
Sample 2	77.8	1.6	2.1
Sample 3	97.2	2.8	2.9

<b>Run to run</b>			
Sample	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	73.1	4.5	6.1
Sample 2	101.9	3.5	3.4
Sample 3	186.2	4.1	2.2

### Method comparison:

A comparison of the Analyticon Turbitex APO B (y) with a commercial obtainable assay (x) gave the following result (mg/dl):

$$y = 0.9438x - 5.4391 \quad r = 0.9845$$

### Quality Control:

Contronorm® L 5 x 2 ml #1302  
Controptath® L 5 x 2 ml #1303

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### Calibration:

Standardization: The Apolipoprotein B method was standardized using the WHO/IFCC SP3-07 reference standard.

Calibration Type: Logit-Log

The kalibrator has one level, the dilutions will be done automatically by the instrument.

S1 – S5: Bio Cal® L 1 x 2 ml #1401  
5 x 2 ml #1402

### **Preparation of calibrator for multi-point calibration:**

Prepare a fresh dilution of Bio Cal<sup>®</sup> L with 0.9 % NaCl solution:

S1: 0.9% NaCl	
S2: 1:20	e.g. 10 µl Bio Cal <sup>®</sup> L + 190 µl 0.9% NaCl
S3: 1:10	e.g. 20 µl Bio Cal <sup>®</sup> L + 180 µl 0.9% NaCl
S4: 1:5	e.g. 40 µl Bio Cal <sup>®</sup> L + 160 µl 0.9% NaCl
S5: 1:2	e.g. 100 µl Bio Cal <sup>®</sup> L + 100 µl 0.9% NaCl
S5: Bio Cal <sup>®</sup> L undiluted	

Calculate calibration values from the package insert of the calibrator (e.g. for 1:20 dilution, divide the given value by 20.)

### **Calibration frequency:**

Full calibration is recommended:

- after lot change
- as required following quality control procedures

### **Disposal:**

Please note the legal regulations.

Text passages with grey background were changed in the latest revision of this package insert.

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\*TRIS= Tris (hydroxymethyl)-aminomethan-hydrochlorid

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