Order information:

Catalog-No.	Contents
H0901 Hit 1 / 917 (Ilab* / AU*)	R1 6 x 20 ml R2 6 x 4.5 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911/917: ACN 151
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of Apolipoprotein B in human serum and plasma.

Summary:

Apolipoprotein B is the main component of the LDL and IDL, and it is also an important component of VLDL and chylomicrons. It occurs in two forms: apo B-100 in the VLDL/LDL and apo B-48 in the chylomicrons.

Apo B-100 originates in the liver, whereas apo B-48 (the N-terminal half of apo B-100) is mainly formed in the small intestine. Apolipoprotein B is important for the transport of cholesterol and triglycerides as well as for their receptor-mediated cellular uptake.

As published in the scientific literature, Apolipoprotein determinations are primarily used in the early detection and assessment of the individual coronary risk and for the diagnosis of the relatively rare (but very serious) apolipoproteinopathy, which is generally inherited. The determination of apolipoprotein B is used in the detection of hyperapobetalipoproteinemia (apoB elevated, LDL-cholesterol normal), which is found in almost half of patients with coronary diseases. A positive angiography finding, i.e. the demonstration of stenosis, is associated with pathological apo A-1/apo B findings. The Apolipoprotein determinations supplement modern HDL, LDL and Lp(a) diagnostics.

The determination of apolipoprotein B is made by immunoassay (RIA; ELISA), electroimmunodiffusion (EID), radial immunodiffusion, nephelometry or turbidimetry.

Test principle:

Immunoturbidimetric assay.
• Sample and addition of R1 (buffer)
• Addition of R2 (anti-apolipoprotein B-antibody/buffer) and start of reaction.

Anti-apolipoprotein B-antibodies react with the antigen in the sample to form antigen/antibody complexes which following agglutination, can be measured turbidimetrically.

Reagent concentration:

R1:	
Tris* buffer pH 7.8	100 mmol/l
PEG	4.5 %
NaCl	150 mmol/l
Preservative; Detergent	

R2:	
Polyclonal anti-human apolipoprotein B	dependent on titer
Antibody (goat)	
Tris* buffer pH 7.8	80 mmol/l
NaCl	150 mmol/l
Preservative	

*TRIS= Tris (hydroxymethyl)-aminomethane

Preparation and stability

R1: Ready for use.
R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C
Onboard stability at +2°C to + 8°C: R1: 42 days
R2: 42 days

Specimen

Collect serum using standard sampling tubes

Na-heparin, or Na-EDTA plasma

Stability:	1 day	at +20°C to +25°C
	3 days	at + 4°C to + 8°C
	2 months	at - 20°C (only freeze once!)

Centrifuge samples containing precipitate before performing the assay.
For analyzers with automatic sample predilution, the samples are prediluted with 0.9% NaCl solution by the instrument.
For analyzers without automatic sample predilution, manual predilution of the samples 1 + 20 with 0.9% NaCl is necessary.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.
Icterus: No significant interference up to 80 mg/dl).
Hemolysis: No significant interference up to an index H of 1000 (approximate haemoglobin concentration: 1000 mg/dl).
Lipemia (Intralipid): No significant interference up to an index L of 600 (approximate triglycerides concentration: 1200 mg/dl). There is poor correlation between turbidity and triglycerides concentration.
No high-dose hook-effect up to an apolipoprotein B concentration of 1000 mg/dl.
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions and calibrators as described above

Additional materials required

- Controls as indicated below
- 0.9% NaCl
- General laboratory equipment

Manual procedure:		
Wavelength:	340 nm	
Temperature:	+37°C	
Cuvette:	1cm	
Zero adjustment:	against reagent blank	
	Blank	Sample/ Calibrator
R1	1000 µl	1000 µl
R2	200 µl	200 µl
Mix, incubate 2 min. and read A ₁ reagent blank. Then add:		
Sample/calibrator	---	20 µl
Mix, incubate 5 min. and read A ₂ .		
Calculation:		
$\Delta A = [(A_2) \text{ sample/calibrator}] - [(A_1) \text{ blank}]$		
The concentration of Apo B in patient sera can be read from a graph using values of 6 levels of standards in the concentration range of 0 to 460 mg/dl APO B. For zero value is recommended to use saline solution (0.9%).		

Measuring/reportable range:

Roche/Hitachi systems automatically calculate the apolipoprotein B concentration of each sample.

Conversion factors: mg/dl x 0.01 = g/l
g/l x 100 = mg/dl

• Roche/Hitachi 911/917 analyzer:
20 - 400 mg/dl (0.20-4.00 g/l)*
Extended measuring with rerun: 10 - 400 mg/dl (0.10 - 4.0 g/l)**

* Maximum reportable range is dependent on the highest standard concentration.
**Maximum extended reportable range is an approximate value on the calculated value of the highest standard.

Expected values:

Men: 66 - 133 mg/dl (0.66 - 1.33 g/l)
Women: 60 - 117 mg/dl (0.60 - 1.17 g/l)

Each laboratory should investigate the transferability of the expected values to its own patients' population and if necessary determine its own reference range. For diagnostic purposes the APO B results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

10.0 mg/dl (0.10 g/l)

The lower detection limit represents the lowest measurable apolipoprotein B concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined within run using controls in an internal protocol (n = 20). The following results were obtained:

Within run			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	61.0	0.47	0.77
Sample 2	104	0.78	0.75
Sample 3	185	1.68	0.91

Reproducibility was determined using human samples and controls in an internal protocol (n = 18). The following results were obtained:

Between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	74.7	3.20	4.28
Sample 2	180	5.87	3.27

Method comparison:

A comparison of the Analyticon TurbiteX APO B (y) with a commercial obtainable assay (x) gave with 45 samples the following result (mg/dl):

$$y = 1.08 x + 3.358 \quad r = 0.984$$

Quality Control:

ControPath® L	5 x 2 ml	#1303
Contronorm® L	5 x 2 ml	#1302

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The apolipoprotein B method was standardized using the WHO/IFCC SP3-07 reference standard.

S1: 0.9% NaCl		
S2 - S5: Bio Cal® L	1 x 2 ml	#1401
	5 x 2 ml	#1402

Preparation of calibrator for multi-point calibration:

• Manual testing procedure:

Prepare a fresh dilution of Bio Cal® L with 0.9 % NaCl solution:

S1: 0.9% NaCl
S2: 1:20 e.g. 10 µl Bio Cal® L + 190 µl 0.9% NaCl
S3: 1:10 e.g. 20 µl Bio Cal® L + 180 µl 0.9% NaCl
S4: 1:5 e.g. 40 µl Bio Cal® L + 160 µl 0.9% NaCl
S5: 1:2 e.g. 100 µl Bio Cal® L + 100 µl 0.9% NaCl
S5: Bio Cal® L undiluted

Calculate calibration values from the package insert of the calibrator (e.g. for 1:20 dilution, divide the given value by 20.)

• Hitachi 911 analyzers:

Multiply the lot-specific Bio Cal L values by the factors given below to determine the standard concentrations for the 6-point calibration curve.

S2: 0.252	S5: 2.084
S3: 0.496	S6: 3.504
S4: 1.527	

• Hitachi 917 analyzers:

Multiply the lot-specific Bio Cal L values by the factors given below to determine the standard concentrations for the 6-point calibration curve.

S2: 0.318	S5: 2.100
S3: 0.600	S6: 3.500
S4: 1.615	

Calibration frequency:

Full calibration is recommended:

- after lot change
- as required following quality control procedures

Disposal:

Please note the legal regulations.

Text passages with grey background were changed in the latest revision of this package insert.

Literature:

1. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
2. Council Directive (90/679/EEG) Official Journal of the European Communities No. L374 form Dec. 31 1990
3. Department of Labour, Occupational Safety and Health Administration, 29 CFR 1910 1030. Occupational Exposure to Bloodborne Pathogens; Final Rule. Federal Register 1991; 56:64 175-182
4. Greiling H, Gressner AM. (Hrsg. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer, 1995:325-326.
5. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag, 1996.
6. Karl J, Engel WD. Determination of Apolipoprotein A I and B without sample dilution. Poster presented at the 57th Meeting of the European Atherosclerosis Society, Lisbon and the IX European Congress of Clinical Chemistry, Cracow, 1991
7. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720.

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*TRIS= Tris (hydroxymethyl)-aminomethan-hydrochlorid

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Literatur:

1. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
2. Department of Labor, Occupational Safety and Health Administration, 29 CFR 1910.1030. Occupational Exposure to Bloodborne Pathogens; Final Rule. Federal Register 1991; 56:64 175-182
3. Greiling H, Gressner AM (Hrsg. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3. Auflage. Stuttgart/New York: Schattauer, 1995:325-326.
4. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes Preenalytical Variables. Broschüre in: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag, 1996.
5. Karl J, Engel WD. Determination of Apolipoprotein A I and B without sample dilution. Poster presented at the 57th Meeting of the European Atherosclerosis Society, Lisbon and the IX European Congress of Clinical Chemistry, Cracow, 1991
6. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720.
7. Richtlinie des Rates (90/679/EWG). Amtsblatt d. Europ. Gem. Nr. L374 vom 31.12.1990:1-12.