### BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B8601	200 / 600	R1	4 x	12 ml
		R2	4 x	20 ml
B8603	300 / 600*	R1	4 x	14 ml
		R2	4 x	20 ml

\*only for instruments with a reagent tray for 20 / 70ml bottles

### Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of antistreptolysin-O in human serum and plasma.

### Summary:

Immunological testing for specific antibodies to streptococcal metabolic products yields important information about previous streptococcus infections. Antibodies are produced against the pathogen and its metabolic products. One example is the antibody to streptolysin O, an enzyme produced by Lancefield group A β-haemolytic streptococci. Determination of antistreptolysin O is performed when toxic and sensitizing associated illnesses occur, such as rheumatic fever (major symptoms: carditis, polyarthritis, chorea minor, subcutaneous nodules, erythema annulare) and poststreptococcal acute glomerulonephritis.

Various methods are available for assaying antistreptolysin O, such latex agglutination and hemolysis inhibition. These ASLO assay is based on the principle of immunological agglutination test principle using latex particles as reaction enhance.

### Test principle:

Immunoturbidimetric assay. Streptolysin O coated on latex reacts with antibodies in the sample to form an antibody/antigen complex. Following agglutination this is measured turbidimetrically.

### Reagent concentration:

#### R1:

TRIS buffer\* pH 8.2 170 mmol/l

#### R2:

Borate buffer pH 8.2 10 mmol/l  
Latex particles coated with streptolysin O 1.7 g/l

\*TRIS= Tris (hydroxymethyl)-aminomethane

### Preparation and stability:

R1: Ready for use.

R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C

Onboard Stability R1: 90 days  
R2: 90 days

**Mix reagent R2 well before first use and once weekly!**

### Specimen:

Collect serum using standard sampling tubes

Heparinized and EDTA-Plasma

Stability: 2 days at +20°C to +25°C  
2 days at +4°C to +8°C  
6 months at - 20°C

Centrifuge samples containing precipitate before performing the assay.

### Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 60 (approximate conjugated and unconjugated bilirubin concentration: 60 mg/dl).

Hemolysis: No significant interference up to an index H of 1200 (approximate haemoglobin concentration: 1200 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 1725 (approximate triglycerides concentration: 3450 mg/dl). There is poor correlation between turbidly and triglycerides concentration.

Rheumatoid factors < 180 IU/l do not interfere.

A high-dose hook-effect may occur at antistreptolysin O concentrations above 4000 IU/ml.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

#### Materials provided

• Working solutions as described above

#### Additional materials required

• Calibrators and controls as indicated below  
• 0.9% NaCl

### Measuring/reportable range:

50 - 300 IU/ml

Determine samples having higher concentrations via the rerun function using 0.9% NaCl solution as diluents.

### Expected values:

Adults: Up to 200 IU/ml

Children: Up to 150 IU/ml

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.

For diagnostic purposes the ASL O results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

### Analytical sensitivity (lower detection limit):

Detection limit: 18 IU/ml

The lower detection limit represents the lowest measurable ASL-O concentration that can be distinguished from zero. It is calculated as three standard deviations of 20 replicates of the lowest standard.

### Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol (n = 20). The following results were obtained:

<u>Within run</u>			
Sample	Mean (IU/ml)	SD (IU/ml)	CV %
Human serum	154.9	4.74	3.1
Protein Control 1	307.3	4.5	1.5
Protein Control 2	454.7	6.36	1.4

<u>Run to run</u>			
Sample	Mean (IU/ml)	SD (IU/ml)	CV %
Human serum	508.5	18.71	3.7
Protein Control 1	121.0	5.97	4.9
Protein Control 2	303.5	12.24	4.0

### Method comparison:

A comparison of the Analyticon Turbitex® ASL O (y) with a commercial obtainable assay (x) gave the following result (IU/ml).

$y = 1.3449x - 39.695$ ;  $r = 0.9840$

### Quality Control:

Human Control Serum:

Protein Control Level 1 3 x 1 ml #7661

Protein Control Level 2 3 x 1 ml #7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### Calibration:

Standardization: This ASL-O method was calibrated against an international standard defined by the WHO for ASL-O.

Calibration Type: Logit-Log

S1: 0.9 % NaCl

S2-S6: Bio Cal® ASL Calibration Set 5 x 1 ml #14150

### Calibration frequency:

Calibration is recommended

- after lot change
- as required following quality control procedures

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\*TRIS= Tris (hydroxymethyl)-aminomethan-hydrochlorid

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