Order information:

Catalog No.	Contents				
7381	R1	1 x	12 ml	R2 1 x 20 ml	
H8601	Hit 911 / 917 (ILab*/ AU*)	R1	4 x	12 ml	R2 4 x 20 ml

(*) Kit contains only reagent barcode for Hitachi system

System information:

Hitachi 911: ACN: 361-400 (user defined method)
 Hitachi 917: ACN/BCN: 901-905 (user defined method)
 For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Immuno-turbidimetric assay for the in vitro quantitative determination of antistreptolysin-O in human serum and plasma on automated clinical chemistry analyzers.

Summary:

Immunological testing for specific antibodies to streptococcal metabolic products yields important information about previous streptococcus infections. Antibodies are produced against the pathogen and its metabolic products. One example is the antibody to streptolysin O, an enzyme produced by Lancefield group A β-haemolytic streptococci. Determination of antistreptolysin O is performed when toxic and sensitising associated illnesses occur, such as rheumatic fever (major symptoms: carditis, polyarthritides, chorea minor, subcutaneous nodules, erythema annulare) and poststreptococcal acute glomerulonephritis.

Various methods are available for assaying antistreptolysin O, such as latex agglutination and hemolysis inhibition. These ASLO assay is based on the principle of immunological agglutination test principle using latex particles as reaction enhance.

Test principle:

Immuno-turbidimetric assay
 • Sample and addition of R1 (buffer)
 • Addition of R2 (streptolysin O latex/buffer) and start of reaction

Streptolysin O coated on latex reacts with antibodies in the sample to form an antibody/antigen complex. Following agglutination this is measured turbidimetrically.

Reagent concentration:

R1:
 TRIS buffer* pH 8.2 170 mmol/l
R2:
 Borate buffer pH 8.2 10 mmol/l
 Latex particles coated with streptolysin O 1.7 g/l

*TRIS= Tris (hydroxymethyl)-aminomethane

Preparation and stability:

R1: Ready for use.
 R2: Ready for use.
Mix reagent R2 well before first use and once weekly!
 Unopened kit components: Up to the expiration date at +2°C to +8°C
 R1: 90 days opened and refrigerated on the analyser
 R2: 90 days opened and refrigerated on the analyser

Specimen:

Collect serum using standard sampling tubes
 Heparinized and EDTA-Plasma
 Stability: 2 days at +20°C to +25°C
 2 days at +4°C to +8°C
 6 months at -20°C
 Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
 The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
 Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.
 Icterus: No significant interference up to an index I of 26 (approximate conjugated and unconjugated bilirubin concentration: 26 mg/dl).
 Hemolysis: No significant interference up to an index H of 1000 (approximate haemoglobin concentration: 1000 mg/dl).
 Lipemia (Intralipid): No significant interference up to an index L of 1000 (approximate triglycerides concentration: 2000 mg/dl). There is poor correlation between turbidity and triglycerides concentration.
 Rheumatoid factors < 180 IU/l do not interfere.
 A high-dose hook-effect may occur at antistreptolysin O concentrations above 4000 IU/ml.
 The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above
- Additional materials required
- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure:

Wavelength: 570 nm
 Temperature: +37°C
 Cuvette: 1 cm
 Zero adjustment: against reagent blank

	Blank	Sample/ Calibrator
Sample/ Calibrator	-	15 µl
R1	500 µl	500 µl
Mix, incubate. 5 min. Then add.		
R2	850 µl	850 µl

Mix, read absorbance A₁ after 30 sec.
 Incubate 5 min. and read absorbance A₂.

Calculation:

$$\Delta A = [(A_2 - A_1) \text{ sample or calibrator}] - [(A_2 - A_1) \text{ blank}]$$

The concentration of ASO in patient sera has to be calculated from ΔA using mathematic function as logit/log or can be read from a graph using values of 6 levels of standards in the concentration range of 0 to 560 IU/ml ASL-O. For zero value is recommended to use saline solution (0.9%)

Measuring/reportable range:

20-600 IU/ml
 At higher concentrations, dilute the sample with 0.9% NaCl (e.g. 1 + 1). Multiply the result by the appropriate factor (e.g. 2).

Expected values:

Adults: Up to 200 IU/ml
 Children: Up to 150 IU/ml

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.
 For diagnostic purposes the ASL O results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 20 IU/ml
 The lower detection limit represents the lowest measurable ASL-O concentration that can be distinguished from zero. It is calculated as three standard deviations of 21 replicates of the lowest standard.

Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol (n = 21). The following results were obtained:

Sample	Within run		
	Mean IU/ml	SD IU/ml	CV %
Human serum	126	1.6	1.3
Protein Control 1	150	1.7	1.1
Protein Control 2	258	2.6	1.0

Sample	Between day		
	Mean IU/ml	SD IU/ml	CV %
Human serum	147	5.0	3.4
Protein Control 1	180	4.6	2.5
Protein Control 2	310	6.6	2.1

Method comparison:

A comparison of the Analyticon Turbitex® ASL O (y) with a commercial obtainable assay (x) gave the following result (IU/ml).
 $y = 1.07x + 2.01$; $r = 0.992$

Quality Control:

Human Control Serum:

Protein Control Level 1	3 x 1 ml	#7661
Protein Control Level 2	3 x 1 ml	#7662
Contronorm® ARC	2 x 1 ml	#7562
Contropath® ARC	2 x 1 ml	#7561

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: This ASL-O method was calibrated against an international standard defined by the WHO for ASL-O.

S1-S5: Bio Cal® ASL 5 x 1 ml #14150

Calibration frequency:

Calibration is recommended

- after lot change
- as required following quality control procedures

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