Fluitest® BIL direct DCA

BILIRUBIN, DCA METHOD

Order information:

| Catalog No. | Contents | | | | | |
|---------------------|----------|-----|-------|----|-----|-------|
| 3348 | R1 | 5 x | 20 ml | R2 | 1 x | 20 ml |
| H3351 Hit I (ILab*) | R1 | 6 x | 47 ml | R2 | 6 x | 11 ml |

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911: ACN 269
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Quantitative determination of direct bilirubin in serum, heparinized plasma or EDTA plasma by DCA method.

Summary:

80-85% of bilirubin originates on degradation of hemoglobin with the other 15-20% being derived from cytochrome, myoglobin and catalases. Unconjugated bilirubin, which binds to plasma albumin, is produced in the course of degradation in the reticuloendothelial system, liver Kupffer cells, spleen and bone marrow. Unconjugated (primary, indirect, water-insoluble) bilirubin is soluble in lipids and toxic. With the aid of the glucuronyl transferase enzyme, bilirubin is conjugated primarily by glucuronic acid in the microsomes of hepatic parenchymal cells. In contrast to unconjugated bilirubin, conjugated (secondary, direct) bilirubin is soluble in water, and is excreted via the kidneys. Bilirubin assays are suitable for evaluating the degree of severity of icteric clinical symptoms as well as for monitoring and objectively assessing these symptoms. Distinguishing between direct and indirect bilirubin is a valuable aid in the differential diagnosis of different forms of jaundice. A direct bilirubin value of < 20 % total bilirubin is an indicator of jaundice of pre-hepatic origin. This value can increase to > 50 % in hepatic and post-hepatic jaundice.

Test principle:

Colorimetric test with stabilized dichloraniline (DCA)

- Sample and addition of R1
- Addition of R2 (DCA) and start of the reaction: Stabilized diazonium salt (DCA) reacts in acidic medium with bilirubin to form an azobilirubin (red dye). In the direct bilirubin assay, only conjugated bilirubin reacts.

The color intensity of the red azo dye formed is directly proportional to the direct (conjugated) bilirubin concentration and can be determined photometrically.

Reagents – contents and concentrations:

| | |
|----------------------|-------------|
| R1: | |
| EDTA-Na ₂ | 0,07 mmol/l |
| Sodium chloride | 6,6 g/l |
| Sulfamic acid | 70 mmol/l |

| | |
|----------------------------------|-------------|
| R2: | |
| 2,4-Dichlorophenyl-Diazoniumsalt | 0,09 mmol/l |
| HCl | 130 mmol/l |
| EDTA-Na ₂ | 0,02 mmol/l |

Preparation and stability:

R1 ready for use.
R2 ready for use.
(Un)opened kit components: Up to the expiration date at +2 to +8°C
Avoid contamination. Do not freeze! Protect reagent 2 from light

Specimen:

Fresh serum, heparinized plasma or EDTA plasma.
Keep out of light and protect the sample from the effects of sunlight.
Perform the assay immediately.
Stability: 2 days at +18 to +25°C
7 days at +2 to +8°C
3 month at < -20°C freeze immediately and only once!
Centrifuge samples containing precipitate before performing the assay.

Note:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Lipemia (Intralipid): No significant interference up to an index L of 500 (approximate triglycerides concentration: 1000 mg/dl). There is poor correlation between turbidity and triglycerides concentration.
Hemolysis: No significant interference up to an index H of 50 (approximate hemoglobin concentration: 50 mg/dl).
Ascorbic acid: No significant interference up to an index A of 30 (approximate ascorbic acid concentration: 30 mg/dl).
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above
- Additional materials required
- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure:

Wavelength: 546 nm (540-560 nm)
Temperature: +25 / +37°C
Cuvette: 1 cm light path
Zero adjustment: against reagent blank

| | Blank | Sample / Calibr. |
|---------------------|---------|------------------|
| Sample / Calibrator | --- | 100 µl |
| Aqua dest. | 100 µl | --- |
| R1 | 1000 µl | 1000 µl |

Mix, incubate for 3–5 min. at assay temperature, read initial absorbance A₁, then add:

| | | |
|----|--------|--------|
| R2 | 200 µl | 200 µl |
|----|--------|--------|

Mix and incubate exactly 5 minutes at 37°C or 10min at 20–25°C. Then read absorbance A₂.

Calculation:

$$\Delta A = (A_2 - A_1) \text{Sample/Calibrator} - (A_2 - A_1) \text{Blank}$$

$$\frac{\Delta A \text{ sample}}{\Delta A \text{ calibrator}} \times \text{calibrator conc.} = \text{direct bilirubin in mg/dl}$$

Measuring /reportable range:

Up to 10 mg/dl (171 µmol/l)
Determine samples having higher concentrations via the rerun function. On instruments without rerun function dilute manually with 0.9% NaCl or distilled/deionized water (e.g. 1 + 1). Multiply the result by the appropriate dilution factor (e.g. factor 2).

Expected values in serum:

Direct bilirubin: < 0.2 mg/dl (3.4 µmol/l)

Analytical sensitivity (lower detection limit):

0.1 mg/dl (1.71 µmol/l)
The lower detection limit represents the lowest measurable Bilirubin concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls in an internal protocol within (between day n = 20). The following results were obtained:

| Sample | Between day | | |
|----------|--------------|------------|------|
| | Mean (mg/dl) | SD (mg/dl) | CV % |
| Sample 1 | 0.35 | 0.01 | 3.34 |
| Sample 2 | 0.75 | 0.01 | 1.00 |
| Sample 3 | 2.13 | 0.02 | 0.71 |

| Sample | Within run | | |
|----------|--------------|------------|------|
| | Mean (mg/dl) | SD (mg/dl) | CV % |
| Sample 1 | 0.36 | 0.01 | 3.12 |
| Sample 2 | 0.76 | 0.01 | 1.46 |
| Sample 3 | 2.07 | 0.03 | 1.30 |

Method comparison:

A comparison of the Analyticon Fluitest BIL-Direct DCA (y) with a commercial obtainable assay (x) gave the following result:

$$y = 0.95x + 0,04; \quad r = 0.995$$

Quality Control:

| | | |
|---------------------|-----------|-------|
| Human Control Serum | | |
| Contronorm® Plus | 5 x 5 ml | #1205 |
| | 20 x 5 ml | #1220 |
| Controptath® Plus | 5 x 5 ml | #1305 |
| | 20 x 5 ml | #1320 |

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Fluitest[®] BIL direct DCA

BILIRUBIN, DCA METHOD

Calibration:

Multicalibrator

| | | |
|------------|-----------|-------|
| Bio Cal® | 20 x 3 ml | #1420 |
| Bio Cal® E | 10 x 3 ml | #1430 |

Text passages with grey background were changed in the latest revision of this package insert.

Literature:

1. Thomas L ed. Clinical Laboratory Diagnostics, 1st ed. Frankfurt: TH-Books Verlagsges, 1998:192-202.
2. Tolman KG, Rej R. Liver function, in: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999.p.1125-77.
3. Rand RN, di Pasqua A. A new diazo method for the determination of bilirubin. Clin Chem 1962;6:570-8

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