# Fluitest® BIL direct

## BILIRUBIN JENDRASSIK/GRÓF METHOD

### BioLyzer® Order information:

Catalog No.	BioLyzer	Contents
B3301	200 / 600	R1 6 x 100 ml
		R2 4 x 40 ml
		R3 4 x 8 ml
B3333	300 / 600*	R1 5 x 47 ml
		R2 5 x 12 ml
		R3 5 x 3 ml

\*only for instruments with a reagent tray for 20 / 70ml bottles

### Intended use:

Quantitative determination of direct bilirubin in serum, heparinized plasma or EDTA plasma by the Jendrassik-Gróf method.

### Summary:

80-85% of bilirubin originates on degradation of hemoglobin with the other 15-20% being derived from cytochrome, myoglobin and catalases. Unconjugated bilirubin, which binds to plasma albumin, is produced in the course of degradation in the reticuloendothelial system, liver Kupffer cells, spleen and bone marrow. Unconjugated (primary, indirect, water-insoluble) bilirubin is soluble in lipids and toxic. With the aid of the glucuronyl transferase enzyme, bilirubin is conjugated primarily by glucuronic acid in the microsomes of hepatic parenchymal cells. In contrast to unconjugated bilirubin, conjugated (secondary, direct) bilirubin is soluble in water, and is excreted via the kidneys. Bilirubin assays are suitable for evaluating the degree of severity of icteric clinical symptoms as well as for monitoring and objectively assessing these symptoms. Distinguishing between direct and indirect bilirubin is a valuable aid in the differential diagnosis of different forms of jaundice. A direct bilirubin value of < 20 % total bilirubin is an indicator of jaundice of pre-hepatic origin. This value can increase to > 50 % in hepatic and post-hepatic jaundice.

### Test principle:

Jendrassik-Gróf method

In the presence of caffeine accelerator, total bilirubin couples with sulfanilic acid to form a red azobilirubin dye. The color intensity is proportional to the bilirubin concentration. Determination of direct bilirubin is performed without caffeine additive.

### Reagents – contents and concentrations:

<b>R1:</b>	
Sodium chloride	150 mmol/l
Preservative	
<b>R2:</b>	
Sulfanilic acid	29 mmol/l
HCl	0.17 N
<b>R3:</b>	
Sodium nitrite	25 mmol/l

### Preparation and stability:

R1 is ready for use.

Stability of R1: 99 days at 2-8°C Protect from light.

Working reagent:

Prepare working reagent by pouring the content of bottle R3 into R2. Mix, the resulting reagent is ready to use.

Stability of working reagent: 5 days at 2-8°C Protect from light.

Unopened kit components: Up to the expiration date at 2-8°C

Close bottles after use.

### Specimen:

Fresh serum, heparinized plasma or EDTA plasma. Hemolysis interferes with the test.

Keep out of light and protect the sample from the effects of sunlight. Perform the assay immediately.

Centrifuge samples containing precipitate before performing the assay.

### Note:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Hemolysis: No significant interference up to an H index of 13 (approximate haemoglobin concentration: 13 mg/dl).

Lipemia (Intralipid): No significant interference up to an L index of 100 (approximate triglycerides concentration: 200 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Protect samples from light to avoid false low values.

Hepatotoxic drugs which cause cholestasis and haemolysis produce elevated recoveries.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

### Measuring /reportable range:

0.1 – 6.0 mg/dl (1.7 – 102 µmol/l)

Determine samples having higher concentrations via the rerun function using 0.9% NaCl solution as diluents.

### Expected values in serum:

Direct bilirubin: < 0.25 mg/dl (4.3 µmol/l)

### Analytical sensitivity (lower detection limit):

0.01 mg/dl (0.17 µmol/l)

The lower detection limit represents the lowest measurable Bilirubin concentration that can be distinguished from zero.

### Imprecision:

Reproducibility was determined using controls in an internal protocol within (between day n = 20). The following results were obtained:

Sample	Run to run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	1.02	0.039	3.8
Sample 2	13.0	0.322	2.47

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	1.00	0.009	0.9
Sample 2	2.60	0.025	1.0
Sample 3	2.85	0.025	0.9

### Method comparison:

A comparison of the Analyticon Fluitest® BIL-Direct (y) with a commercial obtainable assay (x) gave the following result:

$$y = 0.9433x - 0.0984; \quad r = 0.9934$$

### Quality Control:

Human Control Serum

Contronorm® Plus 5 x 5 ml #1205

20 x 5 ml #1220

ControPath® Plus 5 x 5 ml #1305

20 x 5 ml #1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### Calibration:

Calibration Type: Linear

S1: 0.9% NaCl

S2: Bio Cal® E 10 x 3 ml #1430

### Calibration frequency:

Two point calibration is recommended;

- Daily of working reagent
- after lot change
- as required following quality control procedures

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## How to handle extra bottle labels for R2

For the Bilirubin direct test a working solution has to be prepared prior to use by mixing R2 and R3. This working solution has a stability of 5 days.

To ensure that the total volume of reagent can be used even at low sample throughput, the volume of R2 can be divided into half and additional bottle labels for R2 are provided with the kit.

Mix half of the volume of R2 and R3 in an extra system bottle and put one R2 label on the bottle. Prepare the next working solution by pouring the remaining R3 into the remaining bottle of R2 (this bottle already has the appropriate barcode).

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