Fluitest® BIL direct

BILIRUBIN JENDRASSIK/GRÓF METHOD

Order information:

Catalog No.		Contents					
H3301	Hit I (ILab*)	R1	6 x	100 ml	R2	4 x	40 ml
		R3	4 x	8 ml			
H3303	Hit 917 (AU*)	R1	5 x	60 ml	R2	5 x	12 ml
		R3	5 x	3 ml			

(*) Kit contains only reagent barcodes for Hitachi systems

System information:

Hitachi 911: ACN 202

Hitachi 917: ACN 294

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Quantitative determination of direct bilirubin in serum, heparinized plasma or EDTA plasma by the Jendrassik-Gróf method.

Summary:

80-85% of bilirubin originates on degradation of hemoglobin with the other 15-20% being derived from cytochrome, myoglobin and catalases. Unconjugated bilirubin, which binds to plasma albumin, is produced in the course of degradation in the reticuloendothelial system, liver Kupffer cells, spleen and bone marrow. Unconjugated (primary, indirect, water-insoluble) bilirubin is soluble in lipids and toxic. With the aid of the glucuronyl transferase enzyme, bilirubin is conjugated primarily by glucuronic acid in the microsomes of hepatic parenchymal cells. In contrast to unconjugated bilirubin, conjugated (secondary, direct) bilirubin is soluble in water, and is excreted via the kidneys. Bilirubin assays are suitable for evaluating the degree of severity of icteric clinical symptoms as well as for monitoring and objectively assessing these symptoms. Distinguishing between direct and indirect bilirubin is a valuable aid in the differential diagnosis of different forms of jaundice. A direct bilirubin value of < 20 % total bilirubin is an indicator of jaundice of pre-hepatic origin. This value can increase to > 50 % in hepatic and post-hepatic jaundice.

Test principle:

Jendrassik-Gróf method

In the presence of caffeine accelerator, total bilirubin couples with sulfanilic acid to form a red azobilirubin dye. The colour intensity is proportional to the bilirubin concentration. Determination of direct bilirubin is performed without caffeine additive.

Reagents – contents and concentrations:

R1:

Sodium chloride 150 mmol/l

Preservative

R2:

Sulfanilic acid 29 mmol/l

HCl 0.17 N

R3:

Sodium nitrite 25 mmol/l

Preparation and stability:

R1 ready for use.

Working reagent:

Prepare working reagent by pouring the content of bottle R3 into R2. Mix, the resulting reagent is ready to use.

Stability of working reagent: 5 days at 2-8°C, protect from light

Unopened kit components: Up to the expiration date at +2 to +8°C

Close bottles after use.

Specimen:

Fresh serum, heparinized plasma or EDTA plasma.

Haemolysis interferes with the test. Don't use lipemic sera.

Keep out of light and protect the sample from the effects of sunlight.

Perform the assay immediately.

Centrifuge samples containing precipitate before performing the assay.

Note:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Lipemia (Intralipid): Triglycerides may interfere. There is poor correlation between turbidity and triglycerides concentration.

Haemolysis interferes with the test. Usually too low values are obtained.

Protect samples from light to avoid false low values.

Lipemia causes falsely high values.

Hepatotoxic drugs which cause cholestasis and haemolysis produce elevated recoveries.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

Manual procedure:

Performance of a manual test procedure or similar adaption to a (semi)automatic analyzer needs reagent R4 (tartrate/NaOH) and is therefore only possible with the Fluitest Bili-T/D article #1965.

Measuring /reportable range:

Up to 25 mg/dl (425 µmol/l)

Determine samples having higher concentrations via the rerun function. On instruments without rerun function dilute manually with 0.9% NaCl or distilled/deionized water (e.g. 1 + 3). Multiply the result by the appropriate dilution factor (e.g. factor 4).

Expected values in serum:

Direct bilirubin: < 0.25 mg/dl (4.3 µmol/l)

Analytical sensitivity (lower detection limit):

0.03 mg/dl (0.58 µmol/l)

The lower detection limit represents the lowest measurable Bilirubin concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls in an internal protocol within (between day n = 20). The following results were obtained:

Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	0.73	0.02	2.90
Sample 2	2.15	0.04	1.63
Sample 3	0.90	0.03	3.28

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	0.33	0.01	2.38
Sample 2	0.95	0.01	0.80
Sample 3	2.16	0.01	0.62

Method comparison:

A comparison of the Analyticon Fluitest® BIL-Direct (y) with a commercial obtainable assay (x) gave the following result:

$$y = 1,001 x + 0,0245; \quad r = 0.996$$

Quality Control:

Human Control Serum

Contronorm® Plus 5 x 5 ml #1205

20 x 5 ml #1220

Contropath® Plus 5 x 5 ml #1305

20 x 5 ml #1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl

S2: Bio Cal® E 10 x 3 ml #1430

Bio-Cal® 20 x 3 ml #1420

Calibration frequency:

Two point calibration is recommended;

- Daily of working reagent
- after lot change
- as required following quality control procedures

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