# Fluitest® BIL total

## BILIRUBIN JENDRASSIK/GRÓF METHOD

### Order information:

Catalog No.		Contents					
H3201	Hit I (ILab*)	R1	6 x	100 ml	R2	4 x	40 ml
		R3	4 x	8 ml			
H3203	Hit 917 (AU*)	R1	5 x	60 ml	R2	5 x	12 ml
		R3	5 x	3 ml			

(\* ) Kit contains only reagent barcode for Hitachi system.

### System information:

Hitachi 911: ACN 201  
 Hitachi 917: ACN 018  
 For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

### Intended use:

Quantitative determination of total bilirubin in serum, heparinized plasma or EDTA plasma by the Jendrassik-Gróf method.

### Summary:

80-85% of bilirubin originates on degradation of hemoglobin with the other 15-20% being derived from cytochrome, myoglobin and catalases. Unconjugated bilirubin, which binds to plasma albumin, is produced in the course of degradation in the reticuloendothelial system, liver Kupffer cells, spleen and bone marrow. Unconjugated (primary, indirect, water-insoluble) bilirubin is soluble in lipids and toxic. With the aid of the enzyme glucuronyl transferase, bilirubin is conjugated primarily by glucuronic acid in the microsomes of hepatic parenchymal cells. In contrast to unconjugated bilirubin, conjugated (secondary, direct) bilirubin is soluble in water and excreted via the kidneys. Bilirubin assays are suitable for evaluating the degree of severity of icteric clinical symptoms as well as for monitoring and objectively assessing these symptoms. Distinguishing between direct and indirect bilirubin is a valuable aid in the differential diagnosis of different forms of jaundice. A direct bilirubin value of < 20 % total bilirubin is an indicator of jaundice of pre-hepatic origin. This value can increase to > 50 % in hepatic and post-hepatic jaundice.

### Test principle:

Jendrassik-Gróf method: In the presence of caffeine accelerator, total bilirubin couples with sulphanic acid to form a red azobilirubin dye, the color intensity of which is proportional to the bilirubin concentration. Determination of direct bilirubin is performed without caffeine additive.

### Reagents concentration:

<b>R1:</b>	
Caffeine	0.26 mol/l
Sodium benzoate	0.52 mol/l
<b>R2:</b>	
Sulphanilic acid	29 mmol/l
HCl	0.17 N
<b>R3:</b>	
Sodium nitrite	25 mmol/l

### Preparation and stability:

R1 ready for use.  
 Working reagent: Prepare working reagent by pouring the content of bottle R3 into R2. Mix, the resulting reagent is ready to use.  
 Stability of working reagent: 5 days at 2-8°C, protect from light  
 Unopened kit components: Up to the expiration date at +2 to +8°C  
 Close bottles after use.

### Specimen:

Fresh serum, heparinized plasma or EDTA plasma.  
 Hemolysis interferes with the test. Don't use lipemic sera.  
 Keep out of light and protect the sample from the effects of sunlight.  
 Perform the assay immediately.  
 Centrifuge samples containing precipitate before performing the assay.

### Note:

For in vitro diagnostic use.  
 The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.  
 Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Criterion: Recovery within ±10% of initial value.  
 Lipemia (Intralipid): No significant interference up to an L index of 38 (approximate triglycerides concentration: 72 mg/dl). There is poor correlation between turbidity and triglycerides concentration.  
 Hemolysis: No significant interference up to an H index of 40 (approximate hemoglobin concentration: 40 mg/dl).  
 Hemolysis interferes with the test. Usually too low values are obtained.  
 Protect samples from light to avoid false low values.  
 A propranolol metabolite produces false elevated values. Lipemia causes false high values. Hepatotoxic drugs which cause cholestasis and hemolysis produce elevated recoveries.  
 The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

Applications for automated systems are available on request.

#### Materials provided

- Working solutions as described above
- Additional materials required*
- Calibrators and controls as indicated below
- 0.9% NaCl

### Manual procedure:

Performance of a manual test procedure or similar adaption to a (semi)automatic analyzer needs reagent R4 (tartrate/NaOH) and is therefore only possible with the Fluitest Bili-T/D article #1965.

### Measuring /reportable range:

Up to 25 mg/dl (425 µmol/l)  
 Determine samples having higher concentrations via the rerun function. On instruments without rerun function dilute manually with 0.9% NaCl or distilled/deionized water (e.g. 1 + 3). Multiply the result by the appropriate dilution factor (e.g. factor 4).

### Expected values in serum:

Total bilirubin: < 1 mg/dl (17 µmol/l)

### Analytical sensitivity (lower detection limit):

(0.03 mg/dl 0.58 µmol/l)  
 The lower detection limit represents the lowest measurable Bilirubin concentration that can be distinguished from zero.

### Imprecision:

Reproducibility was determined using controls in an internal protocol. The following results were obtained:

Between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	1.00	0.01	1.40
Sample 2	1.84	0.02	1.00
Sample 3	4.79	0.05	1.07

Within run			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	1.12	0.017	1.52
Sample 2	5.33	0.034	0.64
Sample 3	7.12	0.030	0.42

### Method comparison:

A comparison of the Analyticon Fluitest® BILI T/D (y) with a commercial obtainable assay (x) gave the following result:  
 $y = 1,001 x + 0,0245; \quad r = 0,996$

### Quality Control:

Human Control Serum		
Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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