# Fluitest<sup>®</sup> BIL total

## BILIRUBIN WALTERS/GERARDE METHOD

### Order information:

Catalog No.	Contents		
804	R1 2 x	100 ml	R2 1 x 10 ml

### Intended use:

In vitro test for the quantitative determination of total bilirubin in human serum and plasma by the Walters/Gerarde method.

### Summary:

80-85% of bilirubin originates on degradation of hemoglobin with the other 15-20% being derived from cytochrome, myoglobin and catalases. Unconjugated bilirubin, which binds to plasma albumin, is produced in the course of degradation in the reticuloendothelial system, liver Kupffer cells, spleen and bone marrow. Unconjugated (primary, indirect, water-insoluble) bilirubin is soluble in lipids and toxic. With the aid of the glucuronyl transferase enzyme, bilirubin is conjugated primarily by glucuronic acid in the microsomes of hepatic parenchymal cells. In contrast to unconjugated bilirubin, conjugated (secondary, direct) bilirubin is soluble in water, and is excreted via the kidneys. Bilirubin assays are suitable for evaluating the degree of severity of icteric clinical symptoms as well as for monitoring and objectively assessing these symptoms. Distinguishing between direct and indirect bilirubin is a valuable aid in the differential diagnosis of different forms of jaundice. A direct bilirubin value of < 20 % total bilirubin is an indicator of jaundice of pre-hepatic origin. This value can increase to > 50 % in hepatic and post-hepatic jaundice.

### Test principle:

Walters/Gerarde Method

The azobilirubin produced by the reaction between bilirubin and the diazonium salt of sulfanilic acid shows maximum absorption at 555 nm in an acid medium. The intensity of the colour produced is proportional to the quantity of bilirubin which has reacted. In the absence of an accelerator, only conjugated bilirubin reacts. In the presence of an accelerator, the dimethylsulphoxide (DMSO), the non-conjugated bilirubin also participates to the reaction, thus to determine the level of total bilirubin.

### Reagents – contents and concentrations:

#### R1:

Sulfanilic acid	29 mmol/l
HCl	0.17 N
DMSO	7 mol/l

#### R2:

Sodium nitrite	17.1 mmol/l
----------------	-------------

### Preparation and stability:

R1: ready for use.

R2: ready for use.

Stability: at +2°C to +8°C up to the expiration date

Always keep bottles tightly closed and protected from light.

### Specimen:

Fresh serum, heparinized plasma or EDTA plasma.

Hemolysis interferes with the test. Don't use lipaemic sera.

Keep out of light and protect the sample from light.

Perform the test immediately.

Centrifuge samples containing precipitate before performing the assay.

### Note:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Hemolysis: Elevated levels of haemoglobin may interfere.

Lipemia (Intralipid): Elevated levels of triglycerides may interfere. There is poor correlation between turbidity and triglycerides concentration.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Measuring /reportable range:

Up to 10 mg/dl (171 µmol/l)

In case of higher results, dilute sample 1:2 with NaCl solution and repeat test. Multiply result by 2.

### Normal values:

Serum: Up to 1.00 mg/dl (17.1 µmol/l):

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.

### Analytical sensitivity (lower detection limit):

0.01mg/dl or 0.17 µmol/l

The lower detection limit represents the lowest measurable Bilirubin concentration that can be distinguished from zero.

### Testing procedure:

Applications for automated systems are available on request.

#### Materials provided

• Working solutions as described above

#### Additional materials required

• Calibrators and controls as indicated below  
• 0.9% NaCl

### Manual procedure:

Wavelength:	Hg 555 nm (530 – 580 nm)
Temperature:	+25 or +37°C
Cuvette:	1 cm light path
Zero adjustment:	against reagent blank

	Sample blank	sample
R1	1000 µl	1000 µl
R2	-----	50 µl
Sample	100 µl	100 µl

Mix and incubate at room temperature for 10 minutes or 5 minutes at +37°C. Read absorbance of sample against sample blank

### Calculation:

$$\frac{\Delta A \text{ sample}}{\Delta A \text{ calibrator}} \times \text{calibrator conc.} = \text{bilirubin conc.}$$

### Imprecision:

Reproducibility was determined using controls in an internal protocol within (n = 20).

The following results were obtained:

between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	1.06	0.048	4.57
Sample 2	4.49	0.075	1.67
Sample 3	4.70	0.059	1.25

within run			
Sample	Mean mg/dl	SD mg/dl	CV %
Sample 1	1.05	0.018	1.69
Sample 2	2.46	0.027	1.11
Sample 3	4.57	0.028	0.60

### Method comparison:

A comparison of the Analyticon Fluitest BIL-T (y) with a commercial obtainable assay (x) gave the following result:

$$y = 0.944 x + 0.194; \quad r = 0.999$$

### Quality control:

Human Control Serum

Contronorm <sup>®</sup> Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath <sup>®</sup> Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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