### Order information:

Catalog-No.	Contents					
1965	R1	1 x	50 ml	R2	1 x	10 ml
	R3	1 x	100 ml	R4	1 x	100 ml

### Intended use:

In vitro test for the quantitative determination of direct bilirubin and total bilirubin in serum, heparinized plasma or EDTA plasma by the Jendrassik-Gróf method.



### Summary:

80-85% of bilirubin originates on degradation of hemoglobin with the other 15-20% being derived from cytochrome, myoglobin and catalases. Unconjugated bilirubin, which binds to plasma albumin, is produced in the course of degradation in the reticuloendothelial system, liver Kupffer cells, spleen and bone marrow. Unconjugated (primary, indirect, water-insoluble) bilirubin is soluble in lipids and toxic. With the aid of the glucuronyl transferase enzyme, bilirubin is conjugated primarily by glucuronic acid in the microsomes of hepatic parenchymal cells. In contrast to unconjugated bilirubin, conjugated (secondary, direct) bilirubin is soluble in water, and is excreted via the kidneys. Bilirubin assays are suitable for evaluating the degree of severity of icteric clinical symptoms as well as for monitoring and objectively assessing these symptoms. Distinguishing between direct and indirect bilirubin is a valuable aid in the differential diagnosis of different forms of jaundice. A direct bilirubin value of <20% total bilirubin is an indicator of jaundice of pre-hepatic origin. This value can increase to >50% in hepatic and post-hepatic jaundice.

### Test principle:

Jendrassik-Gróf method

In the presence of caffeine accelerator, total bilirubin couples with sulfanilic acid to form a red azobilirubin dye, the color intensity which is proportional to the bilirubin concentration. Determination of direct bilirubin is performed without caffeine additive. The addition of alkaline tartrate causes a transformation from the red azobilirubin dye to a blue dye and the absorbance maximum from 546nm to 578nm.

### Reagents – contents and concentrations:

<b>R1:</b>		
Sulfanilic acid		29 mmol/l
HCl		0.17 mol/l
<b>R2:</b>		
Sodium nitrite		25 mmol/l
<b>R3:</b>		
Caffeine		0.26 mol/l
Sodium benzoate		0.52 mol/l
<b>R4:</b>		
Tartrate		0.93 mol/l
NaOH		1.9 mol/l

### Preparation and stability:

R1 ready for use.

Unopened kit components: Up to the expiration date at + 2 to + 8°C.  
Close bottles after use.

### Specimen:

Fresh serum, heparinized plasma or EDTA plasma.  
Hemolysis interferes with the test. Don't use lipemic sera.  
Keep out of light and protect the sample from the effects of sunlight.  
Centrifuge samples containing precipitate before performing the assay.

Stability:	Total	Direct	
	1 day	2 days	at +20 to +25°C
	7 days	7 days	at +4 to +8°C
	6 months	6 months	at - 20°C

### Note:

For in vitro diagnostic use.  
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.  
Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Criterion: Recovery within ± 10% of initial value.  
Hemolysis: Elevated levels of haemoglobin may interfere.  
Lipemia (Intralipid): Elevated levels of triglycerides may interfere. There is poor correlation between turbidity and triglycerides concentration.  
Protect samples from light to avoid false low values.  
The effect of drugs or their metabolic products on the test is not known in all cases.  
In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

Applications for automated systems are available on request.

#### Materials provided

- Working solutions as described above
- Additional materials required
- Calibrators and controls as indicated below
- NaCl - sol. 0.9%

#### Procedure Bilirubin Total:

Wavelength: Hg 578 nm (560 – 600 nm)  
Temperature: +20 to +25°C  
Cuvette: 1 cm  
Zero adjustment: against sample blank

	Sample blank	sample
R1	200 µl	200 µl
R2	-----	50 µl
R3	1000 µl	1000 µl
sample	200 µl	200 µl

Mix and incubate at room temperature for 10 – 60 minutes. Add:

R4	1000 µl	1000 µl
----	---------	---------

Mix and incubate at room temperature for 5 – 30 minutes. Read absorbance of sample against sample blank.

#### Calculation:

##### with factors

Concentration in mg/dl  $10.8 \times \Delta A_{\text{Bilirubin Total}}$   
Concentration in µmol/l  $185 \times \Delta A_{\text{Bilirubin Total}}$

##### with calibrator

$\frac{\Delta E \text{ sample}}{\Delta E \text{ calibrator}} \times \text{calibrator conc.} = \text{bilirubin conc.}$

#### Procedure Bilirubin Direct:

Wavelength: Hg 546 nm (530 – 555 nm)  
Temperature: +20°C to +25°C  
Cuvette: 1 cm  
Zero adjustment: against sample blank

	Sample blank	sample
R1	200 µl	200 µl
R2	-----	50 µl
NaCl - sol. 0.9%	2000 µl	2000 µl
Sample	200 µl	200 µl

Mix and incubate at room temperature for 5 minutes. Read absorbance of sample against sample blank.

#### Calculation:

##### with factors

Concentration in mg/dl  $13.7 \times \Delta A_{\text{Bilirubin Direct}}$   
Concentration in µmol/l  $234 \times \Delta A_{\text{Bilirubin Direct}}$

##### with calibrator

$\frac{\Delta E \text{ sample}}{\Delta E \text{ calibrator}} \times \text{calibrator conc.} = \text{bilirubin conc.}$

### Measuring /reportable range:

Bilirubin direct: Up to 6.0 mg/dl (102 µmol/l)

Bilirubin total: Up to 14 mg/dl (238 µmol/l)

Dilute determine samples having higher concentrations with 0.9% NaCl or distilled/deionized water (e.g. 1 + 3) and repeat test. Multiply the result by the appropriate dilution factor (e.g. factor 4).

### Expected values in serum:

Total bilirubin up to 1.0 mg/dl or 17.1 µmol/l

Direct bilirubin up to 0.25 mg/dl or 4.3 µmol/l

### Analytical sensitivity (lower detection limit):

0.01 mg/dl (0.17 µmol/l)

The lower detection limit represents the lowest measurable Bilirubin concentration that can be distinguished from zero.

# Fluitest® BIL T/D

## BILIRUBIN JENDRASSIK/GRÓF METHOD

### Imprecision:

Reproducibility was determined using controls in an internal protocol within (n = 20).  
The following results were obtained:

#### BILIRUBIN TOTAL

Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	1.00	0.014	1.40
Sample 2	1.84	0.018	1.00
Sample 3	4.79	0.051	1.07

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	1.14	0.017	1.47
Sample 2	1.82	0.021	1.18
Sample 3	5.14	0.020	0.39

#### BILIRUBIN DIRECT

Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	0.73	0.021	2.90
Sample 2	0.90	0.029	3.28
Sample 3	2.15	0.035	1.63

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample 1	0.33	0.008	2.38
Sample 2	0.95	0.008	0.80
Sample 3	2.16	0.013	0.62

### Method comparison:

A comparison of the Analyticon Fluitest BILI T/D (y) with a commercial obtainable assay (x) gave the following result (mg/dl):

$$y = 1.001 x + 0.0245; \quad r = 0.996$$

### Quality Control:

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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