# Turbitex® C 3c

COMPLEMENT FACTOR C3c

## βBiolzyer® Order information:

Catalog No.	Biolzyer	Contents		
B3731	200 / 600	R1	2 x	20 ml
		R2	2 x	7 ml
B3733	300 / 600*	R1	2 x	11 ml
		R2	2 x	4 ml

\*only for instruments with a reagent tray for 20 / 70ml bottles

## Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of human C3c in human serum and plasma.

## Summary:

Activation of the complement system takes place via a classical and an alternative route. The two pathways come together in a joint terminal branch. As complement factor C3c is a factor common to both pathways, the C3c concentration can be evaluated as a parameter for activation of the complement system. Activation is indicated by lowered values. Additional differentiation can be made by determining C4. If the C4 level is normal, activation of the alternative route is likely. Depressed values are observed in a number of inflammatory and infectious diseases. Primary causes are systemic lupus erythematosus (SLE), rheumatoid arthritis, subacute bacterial endocarditis, viremia, parasitic infections or bacterial sepsis. A considerable decrease in C3c can be found in patients with partial liposystrophy or membranoproliferative glomerulonephritis when the C3c-nephritis factor is present. C3c measurements are used in the diagnosis and treatment of various inflammatory disorders including glomerulonephritis, rheumatoid arthritis and bacterial infections. As an acute phase protein, C3c is produced to an increased extent during inflammatory processes. It is elevated in systemic infections, non-infectious chronic inflammatory conditions (primarily chronic polyarthritis) and physiological states (pregnancy). The elevation rarely exceeds twice the normal value and can mask a reduction in the current consumption. A variety of methods, such as nephelometry, radial immunodiffusion and turbidimetry are available for the determination of complement factor C3c.

## Test principle:

Immunoturbidimetric assay  
Anti-C3c antibodies react with antigen in the sample to form an antigen/antibody complexes which, following agglutination, are measured turbidimetrically.

## Reagent concentration:

**R1:**  
TRIS\* buffer pH 7.8 100 mmol/l  
Polyethylene glycol 4.5%  
NaCl 150 mmol/l  
Preservative

**R2:**  
Anti-human C3c antibody (rabbit) dependent on titer;  
TRIS\* buffer pH 7.8 80 mmol/l  
NaCl 150 mmol/l  
Preservative

\*TRIS= Tris (hydroxymethyl)-aminomethane

## Preparation and stability:

R1: Ready for use.  
R2: Ready for use.

Unopened kit components: Up to the expiration date at +2 to +8°C

Onboard stability R1 84 days  
R2 84 days

## Specimen:

Collect serum using standard sampling tubes  
Na-heparin-, Li-heparin plasma or EDTA-Plasma  
Stability: at +20°C - +25°C 4 days  
at +4°C - +8°C 8 days

Centrifuge samples containing precipitate before performing the assay.

C3 is converted into two fragments C3c and C3dg, which are functionally inactive and cross react with the antibody uses in the test. C3 should be determined in EDTA-plasma and C3c in serum. Two ensure that all of the C3c present has been converted to C3c, the determination should be performed only after 24-48 h or after incubation for 1 h at +37°C.

## Notes:

For in vitro diagnostic use.  
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.  
Exercise the normal precautions required for handling all laboratory reagents.

## Testing procedure:

### Materials provided

• Working solutions as described above

### Additional materials required

• Calibrators and controls as indicated below  
• 0.9% NaCl solution

## Limitation interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 85 (approximate conjugated and unconjugated bilirubin: 85 mg/dl)

Hemolysis: No significant interference up to an index H of 1200 (approximate hemoglobin concentration: 1200 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 1050 (approximate triglycerides concentration: 2100 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Rheumatoid factors up to 130 IU/ml do not interfere.

No high-dose hook-effect is seen up to a C3c concentration of 500 mg/dl.

Monoclonal gammopathy sera of the IgA or IgM type can interfere with the C3c determination.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

## Measuring/reportable range:

3 - 400 mg/dl (0.03 - 4 g/l)

Determine samples having higher concentrations via the rerun function using 0.9% NaCl solution as diluents.

## Expected values:

90 - 180 mg/dl (0.9 - 1.8 g/l)\*

\* Reference values according to CRM 470 Protein Standardization

Each laboratory should investigate the transferability of the expected values to its own patients' population and if necessary determine its own reference range.

For diagnostic purposes the C3c results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

## Analytical sensitivity (lower detection limit)

Detection limit: 0.8 mg/dl (0.008 g/l)

The lower detection limit represents the lowest measurable C 3c concentration that can be distinguished from zero.

## Imprecision:

Reproducibility was determined using controls (n = 20). The following results were obtained:

Sample	Within run		
	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	57.76	1.296	2.2
Sample 2	127.10	2.593	2.0
Sample 3	185.74	7.381	4.0

Sample	Run to run		
	Mean (mg/dl)	SD (mg/dl)	CV %
Human serum 1	182.39	6.143	3.4
Human serum 2	111.63	2.706	2.4
Human serum 3	215.24	10.783	5.0

## Method comparison:

A comparison of the Analyticon Turbitex® C3c (y) with a commercial obtainable assay (x) gave the following result.

$$y = 0.9191x + 8.0516; \quad r = 0.9535$$

## Quality control:

Protein Control Level 1 3 x 1 ml #7661  
Protein Control Level 2 3 x 1 ml #7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

## Calibration:

Standardization: The C3c method was calibrated against the Reference Preparation CRM 470. (RPPHS - Reference Preparation for Proteins in Human serum)

Calibration Type: Logit-Log

S1: NaCl (0.9%)

S2-S6: Bio Cal® P Calibration Set 5 x 1 ml #1475

## Calibration frequency:

Full recalibration is recommended

• after lot change  
• as required following quality control procedures

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