### Order information:

| Catalog No.                     | Contents                   |
|---------------------------------|----------------------------|
| H3701 Hit I / 917 (AU* / ILab*) | R1 2 x 18 ml   R2 2 x 4 ml |

(\* ) Kit contains only reagent barcodes for Hitachi systems.

### System information:

Hitachi 911/917: ACN 036  
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

### Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of human C3c in human serum and plasma.

### Summary:

Activation of the complement system takes place via a classical and an alternative route. The two pathways come together in a joint terminal branch. As complement factor C3c is a factor common to both pathways, the C3c concentration can be evaluated as a parameter for activation of the complement system. Activation is indicated by lowered values. Additional differentiation can be made by determining C4. If the C4 level is normal, activation of the alternative route is likely. Depressed values are observed in a number of inflammatory and infectious diseases. Primary causes are systemic lupus erythematosus (SLE), rheumatoid arthritis, subacute bacterial endocarditis, viremia, parasitic infections or bacterial sepsis. A considerable decrease in C3c can be found in patients with partial liposystrophy or membranoproliferative glomerulonephritis when the C3c-nephritis factor is present. C3c measurements are used in the diagnosis and treatment of various inflammatory disorders including glomerulonephritis, rheumatoid arthritis and bacterial infections. As an acute phase protein, C3c is produced to an increased extent during inflammatory processes. It is elevated in systemic infections, non-infectious chronic inflammatory conditions (primarily chronic polyarthritis) and physiological states (pregnancy). The elevation rarely exceeds twice the normal value and can mask a reduction in the current consumption.

A variety of methods, such as nephelometry, radial immunodiffusion and turbidimetry are available for the determination of complement factor C3c.

### Test principle:

Immunoturbidimetric assay  
• Sample and addition of R1 (buffer)  
• Addition of R2 (anti-C3c antibody/buffer) and start of reaction

Anti-C3c antibodies react with antigen in the sample to form an antigen/antibody complexes which, following agglutination, are measured turbidimetrically.

### Reagent concentration:

|                                   |                     |
|-----------------------------------|---------------------|
| <b>R1:</b>                        |                     |
| TRIS* buffer pH 7.8               | 100 mmol/l          |
| Polyethylene glycol               | 4.5%                |
| NaCl                              | 150 mmol/l          |
| Preservative                      |                     |
| <b>R2:</b>                        |                     |
| Anti-human C3c antibody (rabbit), | dependent on titer; |
| TRIS* buffer pH 7.8               | 80 mmol/l           |
| NaCl                              | 150 mmol/l          |
| Preservative                      |                     |

\*TRIS= Tris (hydroxymethyl)-aminomethane

### Preparation and stability:

R1: Ready for use.  
R2: Ready for use.  
Unopened kit components: Up to the expiration date at +2 to +8°C  
Onboard stability at +2 to +8°C: R1 84 days  
R2 84 days

### Specimen:

Collect serum using standard sampling tubes  
Na-heparin-, Li-heparin plasma or EDTA-Plasma  
Stability: at +20°C - +25°C 4 days  
at +4°C - +8°C 8 days

Centrifuge samples containing precipitate before performing the assay.

C3 is converted into two fragments C3c and C3dg, which are functionally inactive and cross react with the antibody uses in the test. C3 should be determined in EDTA-plasma and C3c in serum. To ensure that all of the C3c present has been converted to C3c, the determination should be performed only after 24-48 h or after incubation for 1 h at +37°C.

### Notes:

For in vitro diagnostic use.  
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.  
Exercise the normal precautions required for handling all laboratory reagents.

### Limitation interference:

Criterion: Recovery within ±10% of initial value.  
Icterus: No significant interference up to a bilirubin concentration of 80 mg/dl.  
Hemolysis: No significant interference up to a haemoglobin concentration of 1100 mg/dl.  
Lipemia (Intralipid): Elevated levels of triglycerides may interfere. There is poor correlation between turbidity and triglycerides concentration.  
Rheumatoid factors up to 1200 IU/ml do not interfere.  
No high-dose hook-effect is seen up to a C3c concentration of 500 mg/dl.  
Monoclonal gammopathy sera of the IgA or IgM type can interfere with the C3 determination.  
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

Applications for automated systems are available on request.

#### Materials provided

• Working solutions as described above

#### Additional materials required

• Calibrators and controls as indicated below  
• 0,9% NaCl solution

| <b>Manual procedure:</b>  |               |                   |
|---|---------------|-------------------|
| Wavelength:   | 340 nm        |                   |
| Temperature:  | +37°C         |                   |
| Cuvette:  | 1 cm          |                   |
| Zero adjustment:  | against blank |                   |
| Sample/Calibrator   | Blank         | Sample/Calibrator |
| R1  | ---           | 9 µl              |
|   | 750 µl        | 750 µl            |
| Mix, incubate 5 min. and read A <sub>1</sub> . Then add:  |               |                   |
| R2  | 150 µl        | 150 µl            |
| Mix, incubate 5 min. and read A <sub>2</sub> .  |               |                   |
| <b>Calculation:</b>   |               |                   |
| $\Delta A = [(A_2 - A_1) \text{ sample or calibrator}] - [(A_2 - A_1) \text{ blank}]$   |               |                   |
| The concentration of C3c in patient sera has to be calculated from $\Delta A$ using mathematic function as logit/log or can be read from a graph using values of 6 levels of standards in the concentration range of 0 to 300 mg/dl C3c. For zero value is recommended to use saline solution (0.9%). |               |                   |

### Measuring/reportable range:

4 - 300 mg/dl (0.04 - 3.0 g/l)\*  
At higher concentrations, dilute the sample with 0.9% NaCl (e.g. 1+1). Multiply the result by the appropriate factor (e.g. 2).  
\*Maximum reportable range is dependent on the highest standard concentration.

### Expected values:

90 - 180 mg/dl (0.9 - 1.8 g/l)\*\*

\*\* Reference values according to CRM 470 Protein Standardization

Each laboratory should investigate the transferability of the expected values to its own patients' population and if necessary determine its own reference range. For diagnostic purposes the C3c results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

### Imprecision:

Reproducibility within run was determined using controls (n = 20). The following results were obtained:

| Sample   | Within run |          |      |
|----------|------------|----------|------|
|          | Mean mg/dl | SD mg/dl | CV % |
| Sample 1 | 75.9       | 1.03     | 1.35 |
| Sample 2 | 101.0      | 0.82     | 0.81 |
| Sample 3 | 152.3      | 1.98     | 1.30 |

Reproducibility between day was determined using controls (n = 20). The following results were obtained:

| Sample        | Between day |          |      |
|---------------|-------------|----------|------|
|               | Mean mg/dl  | SD mg/dl | CV % |
| Human serum 1 | 81.8        | 1.84     | 2.25 |
| Human serum 2 | 162.5       | 3.78     | 2.33 |
| Human serum 3 | 244.6       | 5.79     | 2.37 |

### **Analytical sensitivity (lower detection limit)**

Detection limit: 4 mg/dl (0.04 g/l)

The lower detection limit represents the lowest measurable C 3c concentration that can be distinguished from zero. It is calculated as the concentration lying three standard deviations above that of the lowest standard (zero standard + 3 SD, within run precision, n = 21).

### **Method comparison:**

A comparison of the Analyticon Turbitex<sup>®</sup> C3c (y) with a commercial obtainable assay (x) gave with 40 samples the following result.

$$y = 1.00 x + 2.75; \quad r = 0.997$$

### **Quality control:**

|                         |          |       |
|-------------------------|----------|-------|
| Protein Control Level 1 | 3 x 1 ml | #7661 |
| Protein Control Level 2 | 3 x 1 ml | #7662 |

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### **Calibration:**

Standardization: The C3c method was calibrated against the Reference Preparation CRM 470. (RPPHS - Reference Preparation for Proteins in Human serum)

#### **• Roche/Hitachi 911/917**

|                               |          |       |
|-------------------------------|----------|-------|
| S1: NaCl (0.9%)               |          |       |
| S2-S6: Bio Cal <sup>®</sup> P | 3 x 1 ml | #1470 |

Multiply the lot specific calibrator value by the factors below to determine the standard concentrations for the six-point calibration curve:

|           |           |
|-----------|-----------|
| S2: 0,143 | S5: 1,000 |
| S3: 0,286 | S6: 1,333 |
| S4: 0,350 |           |

### **Calibration frequency:**

Full recalibration is recommended

- after lot change
- as required following quality control procedures

### **Disposal:**

Please note the legal regulations.

Text passages with grey background were changed in the latest revision of this package insert.

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