Order information:

Catalog No.	Contents
H3801 Hit 1 / 917 (AU* / ILab*)	R1 2 x 18 ml R2 2 x 4 ml

(*) Kit contains only reagent barcodes for Hitachi systems.

System information:

Hitachi 911/917: ACN 032
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of human C4 in human serum and plasma.

Summary:

The complement system can be activated via the classical and an alternative route. Complement factor C4 participates in activation by the classical route. A decrease in C4 is common, but complete absence is rare. A lowered concentration or the complete absence of C4 occurs in immunocomplex diseases, systemic lupus erythematosus (SLE), autoimmune thyroiditis and juvenile dermatomyositis. The commencement of SLE in patients with C4-deficiencies can often be detected at a very early stage, and the course of the disease is milder than in patients with normal complement levels. Infections such as bacterial and viral meningitis, streptococcal and staphylococcal sepsis and pneumonia are associated with a decrease in C4. Additional differentiation can be obtained by the determination of C4 when the level of complement factor C3 is low. If in such cases the protein concentration of C4 is normal, then an activation of the alternative route is likely. The main use of C4 determinations is in assessing the course of hypocomplement conditions. As an acute phase protein, C4 is produced to an increased extent during inflammatory processes. It is elevated in systemic infections, non-infectious chronic inflammatory conditions (primarily chronic polyarthritis) and physiological states (pregnancy). The elevation rarely exceeds twice the normal value and can mask a reduction in the current consumption. A variety of methods, such as nephelometry, radial immunodiffusion and turbidimetry are available for the determination of complement factor C 4.

Test principle:

Immunoturbidimetric assay
• Sample and addition of R1 (buffer)
• Addition of R2 (anti-C4 antibody/buffer) and start of reaction

Anti-C 4 antibodies react with antigen in the sample to form an antigen/antibody complexes which, following agglutination, are measured turbidimetrically.

Reagent concentration:

R1:	
TRIS* buffer pH 7.8	100 mmol/l
Polyethylene glycol	4.5%
NaCl	150 mmol/l
Preservative	
R2:	
Anti-human C 4 antibody (rabbit);	dependent on titer
TRIS* buffer pH 7.8	80 mmol/l
NaCl	150 mmol/l
Preservative	

*TRIS= Tris (hydroxymethyl)-aminomethane

Preparation and stability:

R1: Ready for use.
R2: Ready for use.
Unopened kit components: Up to the expiration date at +2°C to +8°C
Onboard stability at +2°C to +8°C: R1 84 days
R2 84 days

Specimen:

Collect serum using standard sampling tubes
Na heparin-, Li-heparin plasma or EDTA-Plasma
Stability: 2 days at +20°C - +25°C
2 days at +4°C - +8°C

Centrifuge samples containing precipitate before performing the assay.

Limitation interference:

Criterion: Recovery within ±10% of initial value.
Icterus: No significant interference up to a bilirubin concentration of 80 mg/dl.
Hemolysis: No significant interference up to a haemoglobin concentration of 1100mg/dl.
Lipemia (Intralipid): Elevated levels of triglycerides may interfere. There is poor correlation between turbidity and triglycerides concentration.
Rheumatoid factors up to 1200 IU/ml do not interfere.
No high-dose hook-effect is seen up to a C 4 concentration of 500 mg/dl.
Monoclonal gammopathy sera of the IgA or IgM type can interfere with the C4 determination.
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been

discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above
- Additional materials required
- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure:		
Wavelength:	340 nm	
Temperature:	+37°C	
Cuvette:	1 cm	
Zero adjustment:	against blank	
sample/calibrator ---	Blank 12 µl	Sample/Calibrator 750 µl
R1	750 µl	750 µl
Mix, incubate 5 min. and read A ₁ . Then add:		
R2	150 µl	150 µl
Mix, incubate 5 min. and read A ₂ .		
Calculation:		
$\Delta A = [(A_2 - A_1) \text{ sample or calibrator}] - [(A_2 - A_1) \text{ blank}]$		
The concentration of C4 in patient sera has to be calculated from ΔA using mathematic function as logit/log or can be read from a graph using values of 6 levels of standards in the concentration range of 0 to 100 mg/dl C4. For zero value is recommended to use saline solution (0.9%)		

Measuring/reportable range:

1.5 - 100 mg/dl (0.015 - 1.0 g/l)*
At higher concentrations, dilute the sample with 0.9% NaCl (e.g. 1 + 1). Multiply the result by the appropriate factor (e.g. 2).

*Maximum reportable range is dependent on the highest standard concentration.

Expected values:

10 - 40 mg/dl (0.1 - 0.4 g/l)**

** Reference values according to CRM 470 Protein Standardization

Each laboratory should investigate the transferability of the expected values to its own patients' population and if necessary determine its own reference range. For diagnostic purposes the C 4 results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

1.5 mg/dl (0.015 g/l)

The lower detection limit represents the lowest measurable C 4 concentration that can be distinguished from zero. It is calculated as the concentration lying three standard deviations above that of the lowest standard (zero standard + 3 SD, within run precision, n = 21).

Imprecision:

Reproducibility within run was determined using controls (n = 20). The following results were obtained:

	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample			
Sample 1	19.45	0.22	1.11
Sample 2	25.52	0.27	1.07
Sample 3	38.13	0.41	1.08

Reproducibility between day was determined using controls (n = 20). The following results were obtained:

	Between day		
	Mean mg/dl	SD mg/dl	CV %
Sample			
Human serum 1	19.8	0.44	2.22
Human serum 2	39.5	0.83	2.09
Human serum 3	59.1	1.17	1.99

Method comparison:

A comparison of the Analyticon Turbitex C4 (y) with a commercial obtainable assay (x) gave with 60 samples the following result.

$$y = 1.021 x + 1.374; r = 0.974$$

Quality control:

Protein Control Level 1	3 x 1 ml	#7661
Protein Control Level 2	3 x 1 ml	#7662

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The C4 method was calibrated against the Reference Preparation CRM 470 (RPPHS - Reference Preparation for Proteins in Human serum).

• Roche/Hitachi 904/911/912/917

S1:	NaCl (0.9%)	
S2-S6:	Bio Cal [®] P	3 x 1 ml #1470

Multiply the lot specific calibrator value by the factors below to determine the standard concentrations for the six-point calibration curve:

S2: 0.204 S5: 1.750

S3: 0.407 S6: 2.500

S4: 1.00

• Roche/Hitachi 704/717/902

S1:	NaCl (0.9%)	
S2-S6:	Bio Cal [®] P Cal set	5 x 1 ml #1475

Calibration frequency

Full recalibration is recommended

- after 90 days
- after lot change
- as required following quality control procedures

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