BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B3931	200 / 600	R1	5 x	25 ml
		R2	5 x	5 ml
B3933	300 / 600*	R1	5 x	25 ml
		R2	5 x	7 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

In vitro assay for the quantitative determination of cholinesterase (CHE, EC 3.1.1.8) in human serum and plasma.

Summary:

Serum cholinesterase (pseudocholinesterase, cholinesterase II or PCHE) is found in the liver, pancreas, heart, serum and in the white matter of the brain. This serum enzyme should not be confused with acetylcholinesterase from erythrocytes (EC 3.1.1.7), which is also referred to as cholinesterase I.


The biological function of cholinesterase is unknown. Clinically, it serves as an indicator of possible insecticide poisoning, and is measured as an index of liver function. Pre-operative screening of cholinesterase is used to detect patients with atypical forms of enzyme and hence avoid prolonged apnea caused by slow elimination of muscle relaxants.

Depressed cholinesterase levels are found in cases of intoxication with organophosphorus compounds and in hepatitis, cirrhosis, myocardial infarction, acute infections and atypical phenotypes of the enzyme.

This assay is based on the DGKC '94 method.

Test principle:

The activity of the cholinesterase is determined according to the following reaction:

Reagent Concentration:

R1:		
Pyrophosphate buffer, pH 7.6	75 mmol/l	
Potassiumhexacyanoferrat (III)	2.0 mmol/l	
R2:		
s-Butyrylthiocholine iodide	15 mmol/l	

Preparation and stability:

All reagents are ready to use.

The reagent is stable at + 2°C to + 8°C up to the date of expiration specified. Avoid direct sunlight. Do not leave bottle 2 open.

On board stability: R1: 28 days
R2: 28 days

Specimen:

Collect serum using standard sampling tubes. Heparinized or EDTA plasma.

Do not use fluoride as anticoagulants.

Stability: 2 days at +20°C to +25°C
7 days at + 4°C to + 8°C
4 weeks at -20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Do not leave bottle 2 open.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 100 (approximate bilirubin concentration 100 mg/dl).

Hemolysis: No significant interference up to an index H of 1200 (approximate hemoglobin concentration: 1200 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 3000 (approximate triglycerides concentration: 6000 mg/dl) There is poor correlation between turbidity and triglycerides concentration.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

Measuring / reportable range:

200 – 16000 U/l or 3.3 – 266.6 µkat/l

Determine samples having higher concentrations via the rerun function using 0.9% NaCl solution as diluent.

Reference value:

Children, Men, Women over 40 years 5.1 to 11.7 KU/l
Women, 16 up to 39 years 4.0 to 12.6 KU/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the CHE results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 75 U/l or 1.25 µkat/l

The lower detection limit represents the lowest measurable cholinesterase activity that can be distinguished from zero.

Imprecision:

Reproducibility was determined using samples in an internal protocol. The following results were obtained.

Sample	Within run		
	Mean U/l	SD U/l	CV %
Control serum 1	5587.6	47.9	0.9
Control serum 2	3154.3	36.0	1.1
Control serum 3	4527.2	49.2	1.1

Sample	Run to run		
	Mean U/l	SD U/l	CV %
Control serum 1	5398.3	98.1	1.8
Control serum 2	2953.9	103.9	3.5

Method comparison:

A comparison of the Analyticon CHE (y) with a commercial obtainable assay (x) gave the following correlation (U/l):

$$y = 0.752x + 173.7; \quad r = 0.9852$$

Quality Control:

Human Control Serum		
Contronorm®Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath®Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Calibration Type: Linear

S1: 0.9% NaCl

S2: Bio Cal® E 10 x 3 ml #1430

Calibration frequency:

Two point calibration is recommended

- after reagent lot change
- as required following quality control procedure

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Text passages with grey background were changed in the latest revision of this package insert.

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