Order information:

Catalog No.	Contents			
6202	R1	8 x	15 ml	R2 1 x 25 ml
H3901 Hit I (ILab*)	R1	5 x	25 ml	R2 5 x 5 ml

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911: ACN 510
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

In vitro assay for the quantitative determination of cholinesterase (CHE, EC 3.1.1.8) in human serum and plasma.

Summary:


Serum cholinesterase (pseudocholinesterase, cholinesterase II or PCHE) is found in the liver, pancreas, heart, serum and in the white matter of the brain. This serum enzyme should not be confused with acetylcholinesterase from erythrocytes (EC 3.1.1.7), which is also referred to as cholinesterase I.

The biological function of cholinesterase is unknown. Clinically, it serves as an indicator of possible insecticide poisoning, and is measured as an index of liver function. Pre-operative screening of cholinesterase is used to detect patients with atypical forms of enzyme and hence avoid prolonged apnea caused by slow elimination of muscle relaxants.

Depressed cholinesterase levels are found in cases of intoxication with organophosphorus compounds and in hepatitis, cirrhosis, myocardial infarction, acute infections and atypical phenotypes of the enzyme. This assay is based on the DGKC '94 method.

Test principle:

The activity of the cholinesterase is determined according to the following reaction:

Reagent Concentration:

R1:	
Pyrophosphate buffer, pH 7.6	75 mmol/l
Potassiumhexacyanoferrat (III)	2.0 mmol/l
R2:	
s-Butyrylthiocholine iodide	15 mmol/l

Preparation and stability:

The reagents are ready for use.
The stability of the reagents is indicated on the label.
Keep at +2°C to +8°C protected from light.
It is not recommended to use a monoreagent, however it is possible to mixture of 5 volumes buffer R1 and 1 volume of substrat reagent 2 in case of manual procedures or analyzers with limited space for reagents. The working reagent is stable for:

24 hours	at +20°C to +25°C
5 days	at + 2°C to + 8°C

Specimen:

Collect serum using standard sampling tubes.
Heparinized or EDTA plasma.
Do not use fluoride as anticoagulants.

Stability:	2 days	at +20°C to +25°C
	7 days	at + 4°C to + 8°C
	4 weeks	at -20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.
Do not leave bottle 2 open.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.
Icterus: No significant interference up to an index I of 45 (approximate bilirubin concentration 45 mg/dl).
Hemolysis: No significant interference up to an index H of 1000 (approximate hemoglobin concentration: 1000 mg/dl).
Lipemia (Intralipid): No significant interference up to an index L of 750 (approximate triglycerides concentration: 1500 mg/dl) There is poor correlation between turbidity and triglycerides concentration.
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

- Working solutions as described above
- Additional materials required
- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure for substrate start:

Wavelength:	Hg 405nm (400-440 nm)
Temperature:	+37°C
Cuvette:	1 cm
Zero adjustment:	air or distilled water

Semi-Micro

R1	1000 µl
Sample	20 µl

Mix and incubate for 5min. Then add:

R2	200 µl
----	--------

Mix again and read initial absorbance after 90 seconds and start stopwatch at the same time. Repeat reading after exactly 30, 60 and 90 seconds and calculate ΔA/min.

Calculation:

$$\Delta A/\text{min} \times 65.8 = \text{Activity (KU/l)}$$

Manual procedure for sample start:

Wavelength:	Hg 405nm (400-440 nm)
Temperature:	+37°C
Cuvette:	1 cm
Zero adjustment:	air or distilled water

Semi-Micro

Working reagent	1000 µl
Sample	20 µl

Mix and read initial absorbance after 90 seconds and start stopwatch at the same time. Repeat reading after exactly 30, 60 and 90 seconds and calculate ΔA/min.

Calculation:

$$\Delta A/\text{min} \times 55.0 = \text{Activity (KU/l)}$$

Measuring / reportable range:

35 – 25000 U/l or 0.58 – 414 µkat/l

Substrate start:

If the decrease of absorbance within one minute is higher than 0.380 at 405 nm or 25 KU, dilute sample 1:5 with saline solution and repeat test. Multiply result by 5.

Sample start:

If the decrease of absorbance within one minute is higher than 0.190 at 405 nm or 12.5 KU, dilute sample 1:5 with saline solution and repeat test. Multiply result by 5.

Expected values:

	+37°C
Children, Men, Women over 40 years	5.1 to 11.7 KU/l
Women, 16 up to 39 years	4.0 to 12.6 KU/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the CHE results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 35 U/l or 0.58 µkat/l

The lower detection limit represents the lowest measurable cholinesterase activity that can be distinguished from zero.

Imprecision:

Reproducibility was determined using samples in an internal protocol. The following results were obtained.

With in run			
Sample	Mean U/l	SD U/l	CV %
Control serum 1	4252	26.7	0.63
Control serum 2	5890	41.3	0.70
Control serum 3	6810	33.4	0.49

Between day			
Sample	Mean U/l	SD U/l	CV %
Control serum 1	4246	163	3.85
Control serum 2	4848	212	4.38
Control serum 3	4591	193	4.2

Method comparison:

A comparison of the Analyticon CHE (y) with a commercial obtainable assay (x) gave the following correlation (U/l):
 $y = 0.763 x + 0.352$; $r = 0.996$

Quality Control:

Human Control Serum

Contronorm®Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath®Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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