CHOL

CHOLESTEROL

Order information:

Catalog No.	Contents	
4046	R1 1 x 110 ml	R2 5 x for 20 ml
		R4 1 x 5 ml
4010	R1 4 x 100 ml	R2 4 x for 100 ml
		R4 1 x 5 ml

Intended use:

Enzymatic in vitro test for the quantitative determination of cholesterol in human serum and plasma.

Summary:

Cholesterol is a steroid with a secondary hydroxyl group in the C₃ position. It is synthesized in many types of tissue, but particularly in the liver and intestinal wall. Approximately three quarters of cholesterol is newly synthesized and one quarter originates from dietary intake. Cholesterol assays are used for screening for atherosclerotic risk and in the diagnosis and treatment of disorders involving elevated cholesterol levels as well as lipid and lipoprotein metabolic disorders.


Cholesterol analysis was first reported by Liebermann in 1885 followed by Burchard in 1889. In the Liebermann-Burchard reaction, cholesterol forms a blue-green dye from polymeric unsaturated carbohydrates in an acetic acid/acetic anhydride/concentrated sulfuric acid medium. The Abell and Kendall method is specific for cholesterol, but is technically complex and requires the use of corrosive reagents. In 1974, Roeschlau and Allain described the first fully enzymatic method.

This method is based on the determination of Δ⁴-cholestenone after enzymatic cleavage of the cholesterol ester by cholesterol esterase. Conversion of cholesterol by cholesterol oxidase, and subsequent measurement by the Trinder reaction of the hydrogen peroxide formed. Optimization of ester cleavage (>99.5%) allows standardization using primary and secondary standards and a direct comparison with the CDC and NIST reference methods.



Test principle:

Enzymatic colorimetric test.

Cholesterol is determined enzymatically using cholesterol esterase and cholesterol oxidase.

Cholesterol esters are cleaved by the action of cholesterol esterase to yield free cholesterol and fatty acids.

Cholesterol is converted by oxygen with the aid of cholesterol oxidase to cholest-4-en 3-one and hydrogen peroxide.

The created hydrogen peroxide forms a red dyestuff by reacting with 4-aminoantipyrine and phenol under the catalytic action of peroxidase. The color intensity is directly proportional to the concentration of cholesterol and can be determined photometrically.

Reagent concentration:

R1:	
Pipes buffer, pH 6.9	90 mmol/l
Phenol	26 mmol/l

R2:	
Cholesterol oxidase	200 U/l
Cholesterol esterase	300 U/l
Peroxidase	1250 U/l
4-Aminoantipyrine	0.4 mmol/l

R4:	
Cholesterol	200 mg/dl (5.17 mmol/l)

Preparation and stability:

Dissolve contents of enzyme reagent R2 with the corresponding volume of buffer/R1.

This working solution is stable for:

48 days	at + 2°C to + 8°C
10 days	at + 20°C to + 25°C

Specimen:

Collect serum using standard sampling tubes.

Heparinized or EDTA plasma. Do not use citrate, oxalate or fluoride plasma.

Stability: 5-7 days	at +2 to +8°C
3 months	at -20°C

Fasting and no fasting samples can be used.

Centrifuge samples containing precipitate before performing assay.

Blood samples should only be drawn prior to the administration of Metaminazole.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 20

(approximate concentration: 20 mg/dl bilirubin)

Hemolysis: No significant interference up to an index H of 350 (approximate hemoglobin concentration: 350 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 85 (approximate triglycerides concentration: 170 mg/dl) There is poor correlation between turbidity and triglycerides concentration.

False low results can occur in patients taking Acetaminophen (Paracetamol), N-Acetyl cysteine (NAC) or Metaminazole (Dipyron). The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- NaCl (0.9%)

Manual procedure:

Wavelength	Hg 546nm (500-550nm)
Temperature	+25/+30/+37°C
Cuvette	1 cm light path
Zero adjustment	one reagent blank per series only

	Blank	Sample/ Calibrator
Sample / Calibrator	---	10 µl
Working reagent	1000 µl	1000 µl

Mix. Incubate 5 min. at +37°C or 10 min. at +20°C to +25°C.
Read absorbance of calibrator and sample against reagent blank within 60 min.

Calculation:

by calibrator:

$$\frac{\Delta A \text{ Sample}}{\Delta A \text{ Calibrator}} \times \text{Calibrator conc.} = \text{Cholesterol conc.}$$

by factor:

Wavelength	c (mg/dl)	c (mmol/l)
Hg 546nm	860 x ΔA	22.2 x ΔA

Expected values:

Clinical Interpretation according to the recommendations of the European atherosclerosis Society.

Lipid metabolic disturbance

Cholesterol	below 200 mg/dl	No
Triglycerides		
Cholesterol	200 – 300 mg/dl	yes if HDL-Cholesterol below 35 mg/dl
Cholesterol	above 300 mg/dl	Yes
Triglycerides	above 300 mg/dl	

Recommendations of the NCEP Adult Treatment Panel for the following risk-cut-off thresholds for the US American population

Desirable cholesterol level:	<200 mg/dl
Borderline high cholesterol:	200-239 mg/dl
High cholesterol:	>240 mg/dl

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the cholesterol results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

At least two measurements of cholesterol on separate occasions should be made before any medical decision is made, since a single point total cholesterol measurement may not represent a patient's usual cholesterol concentration.

CHOL

CHOLESTEROL

Measuring/reportable range:

3-800 mg/dl (0.08-20.7 mmol/l)

To determine samples having higher activities, manually dilute the samples with 0.9% NaCl-solution (e.g. 1 + 2). Multiply the result by the appropriate dilution factor (e.g. factor 3).

Analytical sensitivity (lower detection limit):

Detection limit: 3 mg/dl (0.08 mmol/l)

The lower detection limit represents the lowest measurable cholesterol concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls. The following results were obtained:

Sample	Within run		
	Mean mg/dl	SD mg/dl	% CV
Control serum 1	108	1.52	1.41
Control serum 2	193	4.22	2.19
Control serum 3	282	4.78	1.70

Sample	Between day		
	Mean mg/dl	SD mg/dl	% CV
Control serum 1	102	3.11	3.05
Control serum 2	195	4.41	2.26
Control serum 3			

Method comparison:

A comparison of the Analyticon CHOL (y) with a commercial obtainable assay (x) gave following result:

$$y = 1.012x + 0.026; \quad r = 1.000$$

Quality control:

Human Control Serum:

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Human Lipid Control Serum:

Controptath® L	5 x 2 ml	#1303
Contronorm® L	5 x 2 ml	#1302

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl

S2: Bio Cal® E	10 x 3 ml	#1430
Bio Cal®	20 x 3 ml	#1420
Bio Cal® L	1 x 2 ml	#1401
	5 x 2 ml	#1402

R4: Calibrator provided in kit

Disposal:

Please note the legal regulations.

Literature:

1. Abell L. et al. Standard Methods in Clinical Chemistry 1958; 26:2.
2. Allain C.C. et al. Clin Chem 1974;20:470
3. Burchard H. Beiträge zur Kenntnis der Cholesterine. Dissertation Rostock 1889.
4. Cohn J.S., McNamara J.R., Schaefer E.J.. Lipoprotein Cholesterol Concentrations in the plasma of Human Subjects as Measured in the fed and Fasted States. Clin. Chem 1988;34:2456-2459
5. Glick M.R., Ryder K.W., Jackson S.A.. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474
6. Greinling H., Gressner A.M. eds. Lehrbuch der Klinischen Chemie und Pathobiochemie. 3rd ed. Stuttgart/New York: Schattauer, 1995
7. Liebermann C. Ber Dtsch chem Ges 1885; 18:1803
8. Pisani T., Gebiski C.P., Leary E.T. et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch. Pathol Lab Med 1995;119:1127
9. Recommendations for Improving Cholesterol Measurement: A Report from the Laboratory Standardization Panel of the National Cholesterol Education Program. NIH Publication No. 90-2964, February 1990
10. Roeschlau P. et al. Z Klein Chem Klein Biochem 1974;12:226
11. Second Report of the Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults. NIH Publikation No. 93-3096, September 1993
12. Siedel J. Hägele E.O., Ziegenhorn J. et al. Clin Chem 1983;29:1075.
13. Study Group, European Atherosclerosis Society. Strategies for the prevention of coronary heart disease: A policy statement of the European Atherosclerosis Society. European Heart Journal 1987;8:77
14. Thomas, L., Labor und Diagnose, 5th. ed. (1998)
15. Tietz N.W. eds.. Clinical Guide to Laboratory Tests, 3. Auflage. Philadelphia, Pa: W.B. Saunders Company, 1995:130-131
16. Wiebe D.A., Bernert J.T.. Clin Chem. 1984;30:352

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Literatur:

1. Abell L. et al Standard Methods in Clinical Chemistry 1958; 26:2.
2. Allain C.C. et al. Clin Chem 1974;20:470
3. Burchard H. Beiträge zur Kenntnis der Cholesterine. Dissertation Rostock 1889.
4. Cohn J.S., McNamara J.R., Schaefer E.J. Lipoprotein Cholesterol Concentrations in the plasma of Human Subjects as Measured in the fed and Fasted States. Clin. Chem 1988;34:2456-2459
5. Glick M.R., Ryder K.W., Jackson S.A.. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474
6. Greinling H., Gressner A.M. (Hrsg.) Lehrbuch der Klinischen Chemie und Pathobiochemie. 3rd. Stuttgart/New York: Schattauer, 1995
7. Liebermann C. Ber Dtsch chem Ges 1885; 18:1803
8. Pisani T. , Geebsk C.P., Leary E.T. et al Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch. Pathol Lab Med 1995;119:1127
9. Recommendations for Improving Cholesterol Measurement: A Report from the Laboratory Standardization Panel of the National Cholesterol Education Program. NIH Publication No. 90-2964, February 1990
10. Roeschlau P. et al. Z Klein Chem Klein Biochem 1974;12:226
11. Second Report of the Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults. NIH Publikation No. 93-3096, September 1993
12. Siedel J., Hägele E.O., Ziegenhorn J. et al Clin Chem 1983;29:1075.
13. Study Group, European Atherosclerosis Society. Strategies for the prevention of coronary heart disease: A policy statement of the European Atherosclerosis Society. European Heart Journal 1987;8:77
14. Thomas L., Labor und Diagnose, 5th. ed. (1998)
15. Tietz N.W. (Hrsg.). Clinical Guide to Laboratory Tests, 3.Auflage. Philadelphia, Pa: WB Saunders Company, 1995:130-131
16. Wiebe D.A., Bernert J.T. Clin Chem. 1984;30:352