Order information:

Catalog No.		Contents		
H4501	Hit I (ILab*)	R1 6 x 47 ml	R2 6 x 11 ml	
H4503	Hit 917 (AU*)	R1 6 x 60 ml	R2 6 x 15 ml	
AU4503	AU	R1 6 x 60 ml	R2 6 x 15 ml	

(*) Kit contains only reagent barcode for Hitachi system.

System information:

Hitachi 911/917: ACN 057
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:


In vitro test for the quantitative determination of Creatine Kinase (CK) in human serum and plasma.

Summary:

Creatine kinase (CK) is a dimeric enzyme occurring in four different forms: a mitochondrial isoenzyme and the cytosolic isoenzymes CK-MM (muscle type), CK-BB (brain type) and CK-MB (myocardial type). The determination of CK and CK-isoenzyme activities is utilized in the diagnosis and monitoring of myocardial infarction and myopathies such as the progressive Duchenne muscular dystrophy. Following injury to the myocardium, such as occurs with acute myocardial infarction, CK is released from the damaged myocardial cells. In early cases, a rise in the CK activity can be found just 4 hours after an infarction, the CK-activities reaches a maximum after 12-24 hours and then falls back to the normal range after 3-4 days. Myocardial damage is very likely when the total CK activity is above 190 U/l, the CK-MB activity is above 24 U/l (+37°C) and the CK-MB activity fraction exceeds 6% of the total. The assay method using creatine phosphate and ADP was first described by Oliver, modified by Rosalki and further improved for optimal test conditions by Szasz. CK is rapidly inactivated by oxidation of the sulfhydryl groups in the active center. The enzyme can be reactivated by the addition of acetylcysteine (NAC). Interference by adenylate kinase is prevented by the addition of diadenosine pentaphosphate and AMP. Standardized methods for the determination of CK using the "reverse reaction" and activation by NAC were recommended by the German Society for Clinical Chemistry (DGKC) and the International Federation of Clinical Chemistry (IFCC) in 1977 and 1990 respectively. This assay meets the recommendations of the IFCC and DGKC (Standard method 94).

Test principle:

- UV Test
• Sample and addition of R1
• Addition of R2 and start of reaction:

Equimolar quantities of NADPH and creatine are formed at the same rate. The photometrically measured rate of formation of NADPH is proportional to the CK activity.

Reagent Concentration:

R1:	
Imidazole Buffer, pH 6.7	110 mmol/l
Glucose	21 mmol/l
Mg-Acetate	11 mmol/l
EDTA	2,1 mmol/l
NADP	2,4 mmol/l
N-Acetylcysteine	24 mmol/l
Hexokinase (HK)	≥ 2,5 U/l
Preservatives/Stabilizers	<1 %
R2:	
Tris Buffer, pH 9.1	50 mmol/l
ADP	2,4 mmol/l
AMP	6 mmol/l
Diadenosinpentaphosphate	12 µmol/l
G-6-P-DH	≥ 1,7 U/l
Creatinphosphate	186 mmol/l

Preparation and stability:

Unopened kid components: Up to the expiry date at +2°C to +8°C.

Serum start:

Mix 5 volumes of R1 with 1 volume of R2.

This solution is stable up to

10 days	at + 2°C to + 8°C
1 day	at +20°C to +25°C

Substrate start:

R1: Ready for use; R2: Ready for use

Onboard stability: R1: 28 days
R2: 28 days

Specimen:

Collect serum using standard sampling tubes. Heparinized- or EDTA-plasma

Stability: 7 days at + 4°C to + 8°C
2 days at +20°C to +25°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

Warning! R1 contains hazardous components.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions for handling all laboratory reagents.

Limitation interference:

Criterion: Recovery within ± 10% of initial values.

Icterus: No significant interference up to an index I of 46 (approximate bilirubin concentration: 46 mg/dl).

Hemolysis: No significant interference up to an index H of 1100 (approximate haemoglobin concentration: 1100 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 625 (approximate triglycerides concentration: 1250 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Blood samples should only be drawn prior to the administration of sulfasalazine or sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure for serum start:

Wavelength:	340 nm (Hg 334 nm, Hg 365 nm)		
Temperature:	+37°C		
Cuvette:	1 cm		
Zero adjustment:	air or distilled water		

	Macro	Semi-Micro	Micro
Working reagent	2000 µl	1000 µl	500 µl

Incubate for 30 sec. at 37°C to pre-warm reagent. Then add:

Sample	80 µl	40 µl	20 µl
--------	-------	-------	-------

Mix and incubate for 3 minutes. Measure the absorbance increase per minute for another 3 minutes.

Calculation:

	Macro	Semi-Micro	Micro
340 nm ΔE/min. x	4126	4126	4126
Hg 334 nm ΔE/min. x	4207	4207	4207
Hg 365 nm ΔE/min. x	7429	7429	7429

Manual procedure for substrate start:

Wavelength:	340 nm (Hg 334 nm, Hg 365 nm)		
Temperature:	+37°C		
Cuvette:	1 cm		
Zero adjustment:	air or distilled water		

	Macro	Semi-Micro	Micro
R1	2000 µl	1000 µl	500 µl
Sample	80 µl	40 µl	20 µl

Mix and incubate for 30 sec. at 37°C. Then add:

R2	400 µl	200 µl	100 µl
----	--------	--------	--------

Mix and incubate for 3 minutes. Measure the absorbance increase per minute for another 3 minutes.

Calculation:

	Macro	Semi-Micro	Micro
340 nm ΔE/min. x	4920	4920	4920
Hg 334 nm ΔE/min. x	5016	5016	5016
Hg 365 nm ΔE/min. x	8861	8861	8861

Measuring / reportable range:

3 – 2300 U/l (0.053 – 38.41 µkat/l)

Determine samples with higher activities via the rerun function. On instruments without rerun function, manually dilute these samples with 0.9% NaCl or distilled/deionized water (e.g. 1 + 10). Multiply the result by the appropriate dilution factor (e.g. 11).

Expected values:

On the basis of the optimized standard method of the DGKC (1977) and the IFCC method.

Myocardial infarction: There is high probability of myocardial damage when the following three conditions are fulfilled:

		+25°C	+30°C	+37°C
1	CK Men	>80 U/l (1.33 µkat/l)	>130 U/l (2.17 µkat/l)	>190 U/l* (3.20 µkat/l)
	CK Women	>70 U/l (1.17 µkat/l)	>110 U/l (1.83 µkat/l)	>170 U/l* (2.85 µkat/l)
2	CK-MB	>10 U/l (0.17 µkat/l)	>15 U/l (0.25 µkat/l)	>25 U/l* (0.42 µkat/l)
3	The CK-MB-activity accounts for 6-25% of the total CK activity.			

* Calculated values

If MI is suspected but the values obtained are below the specified limits, a fresh infarct may have occurred. In this case, tests should be repeated after 4 hours.

The following factors were used for converting the reference values from 25°C: 1.53 (+30°C) and 2.38 (+37°C).

CK varies with physical activity level and race in healthy individuals.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the CK results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

3 U/l (0.03 µkat/l)

The lower detection limit represents the lowest measurable CK concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls. The following results were obtained:

Sample	Within run		
	Mean U/l	SD U/l	CV %
Control serum 1	165.08	1.77	1.07
Control serum 2	323.40	5.00	1.55
Control serum 3	481.93	4.80	1.00

Sample	Between day		
	Mean U/l	SD U/l	CV %
Control serum 1	167.53	3.25	1.94
Control serum 2	320.63	10.34	1.67
Control serum 3	493.93	8.23	3.35

Method comparison:

A comparison of the Analyticon Fluitest® CK-NAc (y) with a commercial obtainable assay (x) gave the following result:

$$y = 0.9724x - 0.6476; \quad r = 0.9997$$

Quality Control:

Human Control Serum

Contronorm® Plus 5 x 5 ml #1205

20 x 5 ml #1220

Controptath® Plus 5 x 5 ml #1305

20 x 5 ml #1320

Serum Control CK, CK-MB 2 x 2 ml #1541

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl

S2: Bio Cal® E 10 x 3 ml #1430

BioCal® Specific Enzymes 1 x 2 ml #14575

Calibration frequency:

Two point calibration is recommended

- after reagent lot change
- as required following quality control procedure

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Text passages with grey background were changed in the latest revision of this package insert.

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