CK MB

CREATINE-KINASE MB

Order information:

Catalog No.	Contents			
153	R1	1 x 65 ml	R2	20 x for 3 ml
157	R1	1 x 110 ml	R2	5 x for 20 ml

Intended use:

Immuno-inhibition assay for the quantitative in vitro determination of the CK-MB isoenzyme of creatine kinase in human serum and plasma.

Summary:




Creatine kinase (CK) is a dimeric enzyme occurring in four different forms: a mitochondrial isoenzyme and the cytosolic isoenzymes CKMM (muscle type), CK-BB (brain type) and CK-MB (myocardial type).

The determination of CK-MB is an important element in the diagnosis of myocardial ischemia, e.g. in acute myocardial infarction or myocarditis. CK-MB is detectable in the blood about 3–8 hours after the onset of cardiac symptoms and can be detected over a lengthy period of time. CK-MB may also appear in other clinical conditions such as rhabdomyolysis and stroke. Within the scope of laboratory diagnostics, the determination of total CK, myoglobin and troponin T can contribute to the differentiation of these clinical pictures.

The sensitivity of a CK-MB determination is dependent upon the time at which the sample was taken. Follow-up assays are therefore meaningful.

Test principle:

Human CK-MB is composed of two subunits, CK-M and CK-B which both have an active site. With the aid of polyclonal antibody to CK-M, the catalytic activity of CK-M subunits in the sample is inhibited to 99.6% without affecting the CK-B subunits. The remaining CK-B activity, corresponding to half the CK-MB activity, is determined by the CK-NAC method analogously to total CK. As the CK-BB isoenzyme only rarely appears in serum and the catalytic activity of the CK-M and CK-B sub-units hardly differ, the catalytic activity of the CK-MB isoenzyme can be calculated from the measured CK-B activity by multiplying the result by 2.

Reagent concentration:

R1:	
Imidazole Buffer, pH 6.7	100 mmol/l
Glucose	20 mmol/l
Mg-Acetate	10 mmol/l
Brij 35P	< 5 %
Preservative/ stabilizer	< 0,1 %

R2:	
ADP	5 mmol/l
AMP	5 mmol/l
Diadenosine-P-5-P	10 µmol/l
NADP	2 mmol/l
Creatine phosphate	30 mmol/l
Hexokinase (HK)	2500 U/l
G-6-P-DH	1500 U/l
N-Acetylcysteine	20 mmol/l
Antibody to CK-M	inhibition up to 2000 U/l CK-M

Preparation and stability:

Dissolve contents of enzyme reagent/R2 with the corresponding volume of buffer/R1.

Gently swirl until completely dissolved. Do not shake!

This working reagent is stable for:

10 days	at +2°C to +8°C
3 days	at +20°C to +25°C.

Unopened kit components: Up to the expiration date at +2°C to +8°C

Specimen:

Collect serum using standard sampling tubes.

Heparinized- or EDTA-plasma

Stability:	24 hours	at +4°C:	< 10 %
	1 hour	at +25°C:	< 10 %

Prior to the CK-MB assay, the total CK activity should be determined by the CK NAC method. The antibody is capable of inhibiting up to 2000 U/l CK-M subunit (37°C). Accordingly, CK-MM activities up to 1000 U/l (37°C) are completely inhibited. Therefore, samples with total CK activities above 1000 U/l (37°C) require dilution because complete inhibition is no longer assured.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

Warning! R1 contains hazardous material.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitation interference:

Criterion: Recovery within ± 10% of initial values.

Hemolysis interferes.

Icterus: No significant interference up to an index I of 80 (approximately 80 mg/dl bilirubin).

Lipemia (Intralipid): no significant interference up to 150 (approximate triglyceride concentration: 300 mg/dl)

There is poor correlation between turbidity and triglycerides concentration.

In patients with a disposition to macro-CK formation, implausibly high CK-MB values may be measured in relation to the total CK, since the macroforms mainly consist of CK-B subunits. As these patients have generally not suffered a myocardial infarction, additional diagnostic measures are necessary.

Blood samples should only be drawn prior to the administration of sulfasalazine and sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.

The effect of drugs or their metabolic products on the test is not known in all cases.

In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Manual Procedure:

Determine CK-activity before CK-MB measurement.

Wavelength:	340 nm, Hg 334 or Hg 365 nm
Temperature:	+37°C
Cuvette:	1 cm light path
Zero adjustment:	air or distilled water

Cuvette	Macro	Semi	Micro
Working reagent	2500 µl	1000 µl	500 µl

Incubate for 2min. at 37°C to pre-warm reagent. Then add:

Sample	100 µl	40 µl	20 µl
--------	--------	-------	-------

Mix and incubate for 7 minutes. Read the absorbance A1. Incubate for 5 more minutes and read the absorbance A2.

Calculation:

$\Delta A = A2 - A1$ (for 5 minutes)

	CK-B [U/L]	CK-MB [U/L]
340 nm $\Delta A \times$	825	1651
Hg 334 nm $\Delta A \times$	841	1683
Hg 365 nm $\Delta A \times$	1486	2972

Measuring /reportable range:

5 – 1000 U/l (0.08 – 16.7 µkat/l)

In case of higher activities, dilute samples 1:4 with 0.9% NaCl. Multiply the result by the appropriate dilution factor (e.g. 5).

Expected values:

On the basis on the optimized standard method of the DGKC (1977) and the IFCC method.

Myocardial infarction: There is high probability of myocardial damage when the following three conditions are fulfilled:

		+25°C	+30°C	+37°C
A.	.CK Men	>80 U/l (1.33 µkat/l)	>130 U/l (2.17 µkat/l)	>174 U/l* (2.90 µkat/l)
	CK Women	>70 U/l (1.17 µkat/l)	>110 U/l (1.83 µkat/l)	>140 U/l* (2.33 µkat/l)
B.	CK-MB	>10 U/l (0.17 µkat/l)	>16 U/l (0.25 µkat/l)	>24 U/l* (0.40 µkat/l)
C.	The CK-MB-activity between 6 % and 25% of the total CK activity.			

* Calculated values

If MI is suspected but the values obtained are below the specified limits, a fresh infarct may have occurred. In this case, tests should be repeated after 4 hours.

The following factors were used for converting the reference values from 25°C : 1.53 (30°C) and 2.38 (37°C).

CK varies with physical activity level and race in healthy individuals.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the CK results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

5 U/l bzw. 0.08 µkat/l

The detection limit represents the lowest measurable CK-MB concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls. The following results were obtained:

Between day			
Sample	Mean U/l	SD U/l	% CV
Control serum 1	94	2.39	2.54
Control serum 2	-	-	-
Control serum 3	-	-	-

Within run			
Sample	Mean U/l	SD U/l	% CV
Control serum 1	62.1	2.03	3.27
Control serum 2	249	6.86	2.76
Control serum 3	392	4.43	1.13

Method comparison:

A comparison of the Analyticon CK-MB (y) with a commercial obtainable assay (x)

gave the following result (U/l):

$$y = 1.013 x + 1.726; r = 1.000$$

Quality control:

Human Control Serum

Serum Control CK, CK-MB 2 x 2 ml #1541

For quality control, use a suitable control material. The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl

S2: BioCal® Specific Enzymes 1 x 2 ml #14575

Disposal:

Please note the legal regulations.

Literature:

1. Adams JE, Abendschein DR, Jaffe AS. Biochemical markers of myocardial injury: Is MB Creatine kinase the choice for the 1990? *Circulation* 1993;88:750-763.
2. Apple FS. Diagnostic markers for detection of acute myocardial infarction and reperfusion. *Laboratory Medicine* 1992;23:297-322.
3. Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-474. Passing H.
4. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes Preanalytical Variables, Broschüre in Samples: From The Patient to the Laboratory. Darmstadt: GIT Verlag 1996
5. Jockers-Wretou E, Pileider G. *Clin Chim Acta* 1975;58:223.
6. Neumeier D, Prellwitz W, Würzburg U et al. Determination of creatine kinase isoenzyme MB activity in serum using immunological inhibition of creatine kinase M subunit activity - Activity kinetics and diagnostic significance in myocardial infarction. *Clin Chim Acta* 1976;73:445-451.
7. Remaley AT, Wilding P. Macroenzymes: Biochemical Characterization, Clinical Significance, and Laboratory Detection. *Clin Chem* 1989;35:2261-2270.
8. Rozenman Y, Gotsman MS. The earliest diagnosis of acute myocardial infarction. *Annu Rev Med* 1994;45:31-44.
9. Stein W. *Med. Welt* 1985;36:572
10. Szasz G Busch EW. Third European Congress of Clinical Chemistry Brighton, England 3.-8. Juni 1979 (Abstract)
11. Würzburg U et al. *Klin Wschr* 1976;54:357.

Text passages with grey background were changed in the latest revision of this package insert.

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Literatur:

- 1 Adams JE, Abendschein DR, Jaffe AS. Biochemical markers of myocardial injury: Is MB Creatine kinase the choice for the 1990? Circulation 1993;88:750-763.
- 2 Apple FS. Diagnostic markers for detection of acute myocardial infarction and reperfusion. Laboratory Medicine 1992;23:297-322.
- 3 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474. Passing H.
- 4 Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes Preanalytica Variables, Broschüre in Samples: From The Patient to the Laboratory. Darmstadt: GIT Verlag 1996
- 5 Jockers-Wretou E, Pileider G. Clin Chim Acta 1975;58:223.
- 6 Neumeier D, Prellwitz W, Würzburg U et al. Determination of creatine kinase isoenzyme MB activity in serum using immunological inhibition of creatine kinase M subunit activity - Activity kinetics and diagnostic significance in myocardial infarction. Clin Chim Acta 1976;73:445-451.
- 7 Remaley AT, Wilding P. Macroenzymes: Biochemical Characterization, Clinical Significance, and Laboratory Detection. Clin Chem 1989;35:2261-2270.
- 8 Rozenman Y, Gotsman MS. The earliest diagnosis of acute myocardial infarction. Annu Rev Med 1994;45:31-44.
- 9 Stein W. Med. Welt 1985;36:572
- 10 Szasz G Busch EW. Third European Congress of Clinical Chemistry Brighton, England 3.-8. Juni 1979 (Abstract)
- 11 Würzburg U et al. Klin Wschr 1976;54:357.

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