Order information:

Catalog No.	Contents
H5001 Hit I (ILab*)	R1 6 x 47 ml R2 6 x 11 ml
H5003 Hit 917 (AU*)	R1 6 x 60 ml R2 6 x 15 ml

(*) Kit contains only reagents barcode for Hitachi system

System information:

Hitachi 911: ACN 288
ACN 169 (Blank)

Hitachi 917: ACN 690

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Kinetic in vitro assay using rate-blanking and compensation for the quantitative determination of creatinine in human serum, plasma and urine.


Summary:

In muscle metabolism, creatinine is synthesized endogeneously from creatine and creatine phosphate. Under conditions of normal renal function, creatinine is excreted by glomerular filtration. Creatinine determinations are performed for the diagnosis and monitoring of acute and chronic renal disease as well as for the monitoring of renal dialysis. Creatinine concentrations in urine can be used as reference values for the excretion of certain analytes (albumin, a-amylase).

Test principle:

Kinetic colorimetric assay

- Sample and addition of R1 (sodium hydroxide)
- Addition of R2 (picric acid) and start of reaction:

In alkaline solution, creatinine forms a yellow-orange complex with picrate. The color intensity is directly proportional to the creatinine concentration and can be measured photometrically. Assays using rate-blanking minimize interference by bilirubin. Serum and plasma samples contain proteins which react non specifically in the Jaffe method. Serum and plasma results must be reduced by 0.3mg/dl (26µmol/l) to obtain accurate values. This correction causes a measurement error of ≤ 1% in urine specimens because these do not contain non-specific proteins.

Reagent concentration:

R1:	
NaOH (#H5003)	175 mmol/l
NaOH (#H5001)	250 mmol/l
R2:	
Picric acid	25 mmol/l

Preparation and stability:

R1: Ready for use
R2: Ready for use

Unopened kit components: up to the expiration date at +15° to +25°C.

Onboard stability: R1: 28 days opened and refrigerated on the analyser
R2: 56 days opened and refrigerated on the analyser

Specimen:

Serum, heparinized or EDTA- plasma.

Stability: 7 days at +2° to +8°C
Freeze for long-term storage.

Urine
Stability: 4 days at +2° to +8°C
Freeze for long-term storage.

Collect urine without additives.

Roche/Hitachi 911:

Urine samples are determined with a reduced volume of the sample. This reduction in sample volume is considered in the calculation of the results.

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Analytical sensitivity (lower detection limit)

Detection limit: 0.1mg/dl (8.8 µmol/l)

The lower detection limit represents the lowest measurable creatinine concentration that can be distinguished from zero.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 10 (approximate bilirubin: 10 mg/dl)

Hemolysis: No significant interference up to an index H of 1100 (approximate haemoglobin concentration: 1100 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 1074 (approximate triglycerides concentration: 2148 mg/dl)

Negatively biased results have been reported due to a temporary production of turbidity in the early stages of the reaction. This effect correlates with increased triglycerides in the serum sample, but does not correlate as well with sample Lipemic Index values. The effect disappears after overnight storage.

Do not use creatinine Jaffe for the assay of creatinine in hemolyzed samples from neonates, because fetal hemoglobin (HbF) can cause negative values.

High homogentisic acid concentrations in urine samples lead to false results.

In very rare cases, samples of patients with gammopathy might give falsified results.

The effect of drugs or their metabolic products on the test is not known in all cases.

In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl
- Roche/Chimneys

Measuring /reportable range:

Serum/plasma: 0.1-25 mg/dl (8.8 - 2210 µmol/l)

Roche/Hitachi 911:

Urine: 0.1-500 mg/dl (8.8-44200 µmol/l)

Determine samples with higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl or distilled/deionized water (e.g. 1+1). Multiply the result by the appropriate dilution factor (e.g. 2)

Expected values:

Expected values according to Mazzachi, Peake and Ehrhardt

Serum/plasma (adults)

Men: 0.70 - 1.20 mg/dl (62 - 106 µmol/l)

Women: 0.50 - 0.90 mg/dl (44 - 80 µmol/l)

1st morning urine (adults)

Men: 39 - 259 mg/dl (3450 - 22900 µmol/l)

Women: 28 - 217 mg/dl (2470 - 19200 µmol/l)

24-hour urine

600-2000 mg/24 h (5-18 mmol/24 h)

corresponding to 40-133 mg/d** (3300-12 000 µmol/l**)

**based on a urine volume of 1.5 l/24 h

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the creatinine results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Imprecision:

Serum

Reproducibility was determined using controls. The following results were obtained:

Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Control serum 1	1.21	0.025	2.07
Control serum 2	3.86	0.076	1.97
Control serum 3	4.03	0.074	1.84

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Control serum 1	1.34	0.021	1.57
Control serum 2	4.14	0.039	0.94
Control serum 3	4.17	0.028	0.67

Method comparison:

A comparison of the Analyticon Fluitest® CREA (y) with a commercial obtainable assay (x) gave the following result (mg/dl):

$y = 1.011 x - 0.0093$; $r = 0.971$

Quality control:

Human Control Serum:

Contronorm [®] Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath [®] Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Urine Control:

Urine Control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Absorption of atmospheric CO₂ by the opened reagent bottle R1 leads to impaired calibration stability. This kit therefore requires the use of color-coded chimneys, which reduce the uptake of CO₂ by the reagents. The chimneys should be placed directly into the appropriate reagent: black for R2. The chimneys can be reused for reagent bottles within the same kit.

Chimneys are used all systems except the Roche/Hitachi 704 and 911. Calibration stability on these instruments can be achieved by the use of suitable caps for bottle R2.

Standardization: The creatinine Jaffé method was standardized against GC/MS as a reference method.

S1: 0.9% NaCl or distilled/deionized water

S2: Bio Cal [®] E	10 x 3 ml	#1430
Bio Cal [®]	20 x 3 ml	#1420

Calibration frequency:

Two-point calibration is recommended:

- every 14 days if the reagent bottles are onboard the analyzer for more than 14 days.
- after reagent bottle change if the previous reagent bottles were onboard the analyzer for more than 14 days.
- after reagent lot change
- as required following quality control procedures

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