BioLyzer® Order information:

Catalog No.	BioLyzer	Contents	
B5051	200 / 600	R1	6 x 40 ml
		R2	6 x 21 ml
B5053	300 / 600*	R1	6 x 40 ml
		R2	6 x 21 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:


Enzymatic in vitro test for the quantitative determination of creatinine in human serum, plasma and urine.

Summary:

In muscle metabolism, creatinine is synthesized endogeneously from creatine and creatine phosphate. Under conditions of normal renal function, creatinine is excreted by glomerular filtration. Creatinine determinations are performed for the diagnosis and monitoring of acute and chronic renal disease as well as for the monitoring of renal dialysis. Creatinine concentrations in urine can be used as reference values for the excretion of certain analytes (albumin, α -amylase).

Test principle:

Creatinine is determined using creatininase, creatinase and sarcosin oxidase. H_2O_2 formed in final step of the reactions is used in the Trinder for the development of a chinonimin dye.

Reagent concentration:

R1:	
Goods buffer, pH 8.1	25 mmol/l
Creatinase	≥ 30 KU/l
Sarcosine oxidase	≥ 10 KU/l
Ascorbate oxidase	≥ 2.5 KU/l
Catalase	≥ 350 KU/l
HITB (3-Hydroxy 2,4,6-triiodo benzoic acid)	2.3 mmol/l
Detergents, stabilizers and preservatives	

R2:

Goods buffer, pH 8.1	25 mmol/l
Creatininase	≥ 150 KU/l
Peroxidase	≥ 50 KU/l
4-Aminoantipyrine	2 mmol/l
Potassium hexacyanoferrate	0.18 mmol/l
Stabilizers and preservatives	

Preparation and stability:

R1: Ready for use.

R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C. Protect from light. Do not freeze.

Onboard stability:	R1: 28 days
	R2: 28 days

Specimen:

Serum, heparinized or EDTA- plasma.

Stability:	7 days	at +20°C to +25°C
	7 days	at +2°C to +8°C
	3 month	at -20°C

Urine

Stability:	2 days	at +20°C to +25°C
	6 days	at +2°C to +8°C
	6 month	at -20°C

Collect urine without additives.

Pre-dilute urine 1+9 with distilled water.

Centrifuge samples containing precipitate before performing the assay.

Blood samples should only be drawn prior to the administration of Metaminazole.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within $\pm 10\%$ of initial value.

Icterus: No significant interference up to an index I of 20 (approximate conjugated bilirubin: 20 mg/dl)

Hemolysis: No significant interference up to an index H of 400 (approximate haemoglobin concentration: 400 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 750 (approximate triglycerides concentration: 1500 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Ascorbic acid: No significant interference up to 25 mg/dl.

High homogentisic acid concentrations in urine samples lead to false results.

In very rare cases, samples of patients with gammopathy might give falsified results. False low results can occur in patients taking Acetaminophen (Paracetamol), N-Acetyl cysteine (NAC) or Metaminazole (Dipyron). The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

Measuring/reportable range:

Serum/plasma: 0.02 - 30.0 mg/dl (1.77 - 2650 μ mol/l)

Determine samples with higher concentrations via the rerun function.

Expected values:

Expected values according

Serum/plasma

Women: < 0.5 - 0.9 mg/dl (< 44 - 80 μ mol/l)

Men: < 0.6 - 1.2 mg/dl (< 53 - 107 μ mol/l)

24-hour urine

1.0 - 1.5 g/24 h (8.84-13.3 mmol/24 h)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the creatinine results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 0.02 mg/dl (1.77 μ mol/l)

The lower detection limit represents the lowest measurable creatinine concentration that can be distinguished from zero.

Imprecision:

Reproducibility was determined using controls. The following results were obtained:

<u>Within run</u>			
Sample	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	1.25	0.008	0.7
Sample 2	2.98	0.021	0.7
Sample 3	0.69	0.006	0.8
<u>Between day</u>			
Sample	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	1.26	0.033	2.6
Sample 2	3.01	0.096	3.2
Sample 3	3.94	0.102	2.6

Method comparison:

A comparison of the Analyticon Fluitest CREA PAP (y) with a commercial obtainable assay (x) gave the following result (mg/dl):

$$y = 1.066x + 0.0071; \quad r = 0.9976$$

Quality control:

Human Control Serum:

Contronorm® Plus 5 x 5 ml #1205

20 x 5 ml #1220

ControPath® Plus 5 x 5 ml #1305

20 x 5 ml #1320

Human Urine Control:

Urine Control Set 8 x 5 ml #1507

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Calibration Type: Linear

S1: 0.9% NaCl

S2: Bio Cal® E 10 x 3 ml #1430

Calibration frequency:

Two point calibration is recommended

- after reagent lot change
- as required following quality control procedures

Disposal:

Please note the legal regulations.

Text passages with grey background were changed in the latest revision of this package insert.

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