### Order information:

Catalog No.	Contents			
4281	R1	4 x	40 ml	R2 4 x 20 ml
				R4 1 x 5 ml
H5051 Hit I (ILab*)	R1	6 x	40 ml	R2 6 x 21 ml
H5053 Hit 917 (AU*)	R1	6 x	40 ml	R2 6 x 21 ml
AU5053 AU	R1	6 x	40 ml	R2 6 x 21 ml

(\* ) Kit contains only reagent barcode for Hitachi system.

### System information:

Hitachi 911: ACN 063

Hitachi 917: ACN 652

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

### Intended use:


Enzymatic in vitro test for the quantitative determination of creatinine in human serum, plasma and urine.

### Summary:

In muscle metabolism, creatinine is synthesized endogenously from creatine and creatine phosphate. Under conditions of normal renal function, creatinine is excreted by glomerular filtration. Creatinine determinations are performed for the diagnosis and monitoring of acute and chronic renal disease as well as for the monitoring of renal dialysis. Creatinine concentrations in urine can be used as reference values for the excretion of certain analytes (albumin, a-amylase).

### Test principle:

Creatinine is determined using creatininase, creatinase and sarcosine oxidase. H<sub>2</sub>O<sub>2</sub> formed in final step of the reactions is used in the Trinder for the development of a quinoneimine dye.

### Reagent concentration:

<b>R1:</b>	
Goods buffer, pH 8.1	25 mmol/l
Creatinase	≥30 KU/l
Sarcosine oxidase	≥10 KU/l
Ascorbate oxidase	≥2.5 KU/l
Catalase	≥350 KU/l
HTIB (3-Hydroxy 2,4,6-triiodo benzoic acid)	2.3 mmol/l
Detergents, stabilizers and preservatives	
<b>R2:</b>	
Goods buffer, pH 8.1	25 mmol/l
Creatininase	≥150 KU/l
Peroxidase	≥50 KU/l
4-Aminoantipyrene	2 mmol/l
Potassium hexacyanoferrate	0.18 mmol/l
Stabilizers and preservatives	
<b>R4:</b>	
Creatinine	2 mg/dl (176.8 μmol/l)

### Preparation and stability:

R1: Ready for use.

R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to +8°C. Protect from light. Do not freeze.

Onboard stability: R1: 28 days

R2: 28 days

### Specimen:

Serum, heparinized or EDTA- plasma.

Stability:	7 days	at +20°C to +25°C
	7 days	at +2°C to +8°C
	3 month	at - 20°C

Urine

Stability:	2 days	at +20°C to +25°C
	6 days	at +2°C to +8°C
	6 month	at - 20°C

Collect urine without additives.

Pre-dilute urine 1+9 with distilled water.

Centrifuge samples containing precipitate before performing the assay.

Blood samples should only be drawn prior to the administration of Metamizole.

### Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an index I of 20 (approximate conjugated bilirubin: 20 mg/dl)

Hemolysis: No significant interference up to an index H of 400 (approximate haemoglobin concentration: 400 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 750 (approximate triglycerides concentration: 1500 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Ascorbic acid: No significant interference up to 25 mg/dl.

High homogentisic acid concentrations in urine samples lead to false results.

In very rare cases, samples of patients with gammopathy might give falsified results. False low results can occur in patients taking Acetaminophen (Paracetamol), N-Acetyl cysteine (NAC) or Metamizole (Dipyrone). The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

### Testing procedure:

Applications for automated systems are available on request.

#### Materials provided

• Working solutions as described above

#### Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

<b>Manual procedure:</b>		
Wavelength:	Hg 546 nm (500 - 550 nm)	
Temperature:	+37°C	
Cuvette:	1 cm light path	
Zero adjustment:	reagent blank	
	Blank	Sample/ standard
Sample/ standard	---	24 μl
Dist. Water	24 μl	---
Reagent 1	1000 μl	1000 μl
Mix, incubate. 5 min and read absorbance A <sub>1</sub> , then add.		
Reagent 2	500 μl	500 μl
Mix and read absorbance A <sub>2</sub> after exactly after 5 min.		
<b>Calculation:</b>		
$\Delta A = A_2 - A_1$		
$\frac{\Delta A \text{ Sample}}{\Delta A \text{ Calibrator}} \times \text{Calibrator conc.} = \text{Creatinine conc. (mg/dl)}$		

### Measuring/reportable range:

Serum/plasma: 0.03 - 30.0 mg/dl (2.65 - 2650 μmol/l)

Determine samples with higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl or distilled/deionized water (e.g. 1+1). Multiply the result by the appropriate dilution factor (e.g. 2)

### Expected values:

#### Expected values according

##### Serum/plasma

Women: < 0.5 - 0.9 mg/dl (< 44 - 80 μmol/l)

Men: < 0.6 - 1.2 mg/dl (< 53 - 107 μmol/l)

##### 24-hour urine

1.0 - 1.5 g/24 h (8.84-13,3 mmol/24 h)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the creatinine results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

### Analytical sensitivity (lower detection limit)

Detection limit: 0.03 mg/dl (2.65 μmol/l)

The lower detection limit represents the lowest measurable creatinine concentration that can be distinguished from zero.

### Method comparison:

A comparison of the Analyticon Fluitest CREA PAP (y) with a commercial obtainable assay (x) gave the following result (mg/dl):

$$y = 1.000 x + 0.03; \quad r = 0.998$$

### **Imprecision:**

Reproducibility was determined using controls. The following results were obtained:

<b>Within run</b>			
Sample	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	1.29	0.01	0.91
Sample 2	2.81	0.02	0.63
Sample 3	4.28	0.03	0.69
<b>Between day</b>			
Sample	Mean (mg/dl)	SD (mg/dl)	CV %
Sample 1	1.27	0.03	2.35
Sample 2	2.73	0.02	0.80
Sample 3	4.24	0.03	0.72

### **Quality control:**

Human Control Serum:

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Human Urine Control:

Urine control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### **Calibration:**

S1: 0.9% NaCl

S2: Bio Cal® E	10 x 3 ml	#1430
Bio Cal®	20 x 3 ml	#1420

R4: Creatinine calibrator provided in kit

### **Calibration frequency:**

Two point calibration is recommended

- after reagent lot change
- as required following quality control procedures

### **Disposal:**

Please note the legal regulations.

Text passages with grey background were changed in the latest revision of this package insert.

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