Order information:

Catalog No.	Contents		
448	R1 2 x	100 ml	R2 2 x 100 ml
	R4 1 x	5 ml	

Intended use:


Kinetic in vitro assay using rate-blanking and compensation for the quantitative determination of creatinine in human serum, plasma and urine.

Summary:

In muscle metabolism, creatinine is synthesized endogeneously from creatine and creatine phosphate. Under conditions of normal renal function, creatinine is excreted by glomerular filtration. Creatinine determinations are performed for the diagnosis and monitoring of acute and chronic renal disease as well as for the monitoring of renal dialysis. Creatinine concentrations in urine can be used as reference values for the excretion of certain analytes (albumin, α -amylase).

Test principle:

Kinetic colorimetric assay

In alkaline solution, creatinine forms a yellow-orange complex with picrate. The color intensity is directly proportional to the creatinine concentration and can be measured photometrically.

Reagent Concentration:

R1:	
Picric acid	70 mmol/l
R2:	
Sodium hydroxide	0.32 mol/l
R4:	
Creatinine	2 mg/dl (176.8 μ mol/l)

Preparation and stability:

Dilute 1 volume of picric acid solution/R1 with 1 volume of sodium hydroxide solution /R2.

The solution mixture is stable for: 10 days at + 4° to + 8°C or
1 day at +20° to +25°C

Specimen:

Serum, heparinized or EDTA- plasma.

Stability: 7 days at +2° to +8°C
Freeze for long-term storage.

Urine
Stability: 4 days at +2° to +8°C
Freeze for long-term storage.

Collect urine without additives.

Urine samples are diluted manually with 0.9% NaCl or distilled/deionised water (e.g. 1+10). Multiply the result by the appropriate dilution factor (e.g. 11).

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations interference:

Criterion: Recovery within $\pm 10\%$ of initial value.

Icterus: No significant interference up to an index I of 10 (approximate bilirubin: 10 mg/dl)

Hemolysis: No significant interference up to an index H of 1100 (approximate haemoglobin concentration: 1100 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 1074 (approximate triglycerides concentration: 2148 mg/dl)

Negatively biased results have been reported due to a temporary production of turbidity in the early stages of the reaction. This effect correlates with increased triglycerides in the serum sample, but does not correlate as well with sample Lipemic Index values. The effect disappears after overnight storage.

Do not use creatinine Jaffe for the assay of creatinine in hemolyzed samples from neonates, because fetal hemoglobin (HbF) can cause negative values.

High homogenetic acid concentrations in urine samples lead to false results.

In very rare cases, samples of patients with gammopathy might give falsified results. The effect of drugs or their metabolic products on the test is not known in all cases.

In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

Manual Procedure:

Wavelength:	492 nm (480 - 520 nm)
Temperature:	+25° / +30° / +37°C
Cuvette:	1 cm light path
Zero adjustment:	air or distilled water

Since the reaction is highly temperature sensitive, care must be taken to ensure that the solutions are preheated to an exact temperature and the reaction must proceed at a constant temperature.

The reaction temperature of the standard and of the sample must be identical! Provided these conditions are fulfilled, the determination can be carried out at any chosen temperature between +20 to +37°C.

	Semi micro	Micro
Working solution	1000 μ l	500 μ l
Sample/Calibrator	100 μ l	50 μ l

Mix. After exactly 20 seconds read A_1 of sample and calibrator. Exactly 80 seconds after the first reading, read A_2 of sample and calibrator. Calculate ΔA .

Calculation:

$$\Delta A = A_2 - A_1$$

$$\frac{\Delta A_{\text{Sample}}}{\Delta A_{\text{Calibrator}}} \times \text{calibrator conc.} = \text{Creatinine conc.}$$

Measuring /reportable range:

Serum/plasma: 0.1-10 mg/dl (8.8 - 884 μ mol/l)

Urine: 0.1-200 mg/dl (8.8-17680 μ mol/l)

Determine samples having higher concentrations, manually dilute with 0.9% NaCl or distilled/deionized water (e.g. 1+1). Multiply the result by the appropriate dilution factor (e.g. 2)

Expected values:

Expected values according to Mazzachi, Peake and Ehrhardt

Serum/plasma (adults)

Men: 0.70 - 1.20 mg/dl (62 - 106 μ mol/l)

Women: 0.50 - 0.90 mg/dl (44 - 80 μ mol/l)

1st morning urine (adults)

Men: 30 - 259 mg/dl (3450 - 22900 μ mol/l)

Women: 28 - 217 mg/dl (2470 - 19200 μ mol/l)

24-hour urine

600-2000 mg/24 h (5-18 mmol/24 h) corresponding to 40-133 mg/dl** (3300-12 000 μ mol/l**)

**based on a urine volume of 1.5 l/24 h

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the creatinine results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 0.1mg/dl (8.8 μ mol/l)

The lower detection limit represents the lowest measurable creatinine concentration that can be distinguished from zero.

Imprecision:

Serum

Reproducibility was determined using controls. The following results were obtained:

Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Control serum 1	1.21	0.025	2.07
Control serum 2	3.86	0.076	1.97
Control serum 3	4.03	0.074	1.84

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Control serum 1	1.08	0.015	1.39
Control serum 2	2.45	0.039	1.59
Control serum 3	8.01	0.051	0.64

Method comparison:

A comparison of the Analyticon Fluitest CREA Kinetic (y) with a commercial obtainable assay (x) gave the following result (mg/dl):

$$y = 1.075x - 0.124;$$

$$r = 0.997$$

Quality control:

Human Control Serum:

Contronorm [®] Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath [®] Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Human Urine Control:

Urine control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The creatinine Jaffé method was standardized against GC/MS as a reference method.

S1: 0.9% NaCl or distilled/deionized water

S2: Bio Cal [®] E	10 x 3 ml	#1430
Bio Cal [®]	20 x 3 ml	#1420

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