Turbitex® CRP Latex

C-REACTIVE PROTEIN

Order information:

Catalog No.		Contents		
H5401	Hit I (ILab*)	R1	2 x 15 ml	R2 2 x 15 ml
H5403	Hit 917 (AU*)	R1	2 x 15 ml	R2 2 x 15 ml
AU5404	AU	R1	3 x 60 ml	R2 3 x 60 ml

(*) Kit contains only reagent barcodes for Hitachi systems.

System information:

Hitachi 911: ACN 361-400 (user defined method)
Hitachi 917: ACN 901-905 (user defined method)
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of CRP in human serum and plasma.

Summary:

C-reactive protein is the classic acute phase protein to inflammatory reactions. It is synthesized by the liver and consists of five identical polypeptide chains forming a five-membered ring of molar mass 120 000 daltons. CRP is the most sensitive of the acute phase reactants, and its concentration increases rapidly during inflammatory processes. Complexed CRP activates the complement system beginning with C1q. CRP then initiates opsonization and phagocytosis of invading cells, but its main function is to bind and detoxify endogenous toxic substances produced as a result of tissue damage.

CRP assays are used to detect systematic inflammatory processes (apart from certain types of inflammation such as SLE and colitis ulcerosa); to assess treatment of bacterial infections with antibiotics; to detect intrauterine infections with concomitant premature amniorrhexis; to differentiate between active and inactive forms of disease with concurrent infection, e.g. in patients suffering from SLE or Colitis ulcerosa; to therapeutically monitor rheumatic diseases and assess anti-inflammatory therapy; to determine the presence of postoperative complications at an early stage, such as infected wounds, thrombosis and pneumonia; and to distinguish between infection and bonemarrow transplant rejection. A variety of methods, such as nephelometry and turbidimetry are available for the determination of CRP. The Analyticon Turbitex CRP Latex assay is based on the principle of particle-enhanced immunological agglutination.

Test principle:

Particle-enhanced Immunoturbidimetric assay

- Sample and addition of R1 (buffer)
- Addition of R2 (Anti-CRP antibody-latex) and start of reaction: Anti-CRP antibodies coupled to latex microparticles react with antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically.

Reagent concentration:

R1:
TRIS* buffer pH 7,4 16 mmol/l
Preservative

R2:
Latex particles coated with anti-CRP (rabbit): 0.17 %
Glycine buffer: pH 8,0 50 mmol/l
Preservative.

*TRIS= Tris (hydroxymethyl)-aminomethane

Preparation and stability:

R1: Ready for use.
R2: Ready for use.
Unopened kit components: Up to the expiration date at +2°C to +8°C.
Onboard stability at +2°C to +8°C:
R1: 90 days
R2: 90 days

Note: Mix the latex reagent (R2) well before placing on the analyzer and once a week thereafter.

Specimen:

Collect serum using standard sampling tubes
Li-/Na- heparin, EDTA- or citrate plasma.
Stability: 3 days at +20°C to +25°C
8 days at +4°C to +8°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.
Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.
Icterus: No significant interference up to a bilirubin concentration of 60 mg/dl.
Hemolysis: No significant interference up to a haemoglobin concentration of 500 mg/dl.
Lipemia (Intralipid): No significant interference up to a triglycerides concentration of 1500 mg/dl.
There is poor correlation between turbidity and triglycerides concentration.
Rheumatoid factors up to 100 IU/ml do not interfere.
A high-dose hook-effect does not occur at CRP concentrations below 1000 mg/l.
The effect of drugs or their metabolic products on the test is not known in all cases.
In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above
- Additional materials required
- Calibrators and controls as indicated below
- 0.9% NaCl-solution

For performing the assay on Hitachi-analyzers refer to the appropriate operator's manual for user defined methods.

Hitachi 917: Prior to testing, put the barcodes, according to your chosen application code, on the bottles.

e.g. ACN 901 → put R1 barcode 901 on bottle R1 and R2 barcode 901 on bottle R2.

Place R1 (buffer) on reagent plate 1 and R2 (CRP antibody latex) on reagent plate 2.

Measuring /reportable range:

0.1-20 mg/L
When the CRP concentration of the sample is above the measuring range, manually dilute with 0.9% NaCl or CRP-free serum (e.g. 1+2). Multiply the result by the appropriate dilution factor (e.g. 3) or determine with rerun-function.

Expected values:

IFCC/CRM 470**	
mg/dl	mg/l
< 0.5	< 5

**Reference values according to CRM 470 Protein-Standardized

Elevated CRP range:

It is important to monitor the CRP concentration during the acute phase of the illness.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the CRP results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings

Analytical sensitivity (lower detection limit)

Detection limit: 0.01mg/dl (0.1 mg/l)
The lower detection limit represents the lowest measurable CRP- concentration that can be distinguished from zero.

Imprecision:

Reproducibility within run was determined using human samples and controls (n = 21). The following results were obtained:

Sample	Within run		
	Mean mg/l	SD mg/l	CV %
Sample 1	2.74	0.034	1.24
Sample 2	8.64	0.068	0.78
Sample 3	17.43	0.135	0.77

Reproducibility between day was determined using human samples and controls (n = 21). The following results were obtained:

Sample	Between day		
	Mean mg/l	SD mg/l	CV %
Sample 1	2.65	0.111	4.20
Sample 2	8.55	0.092	1.08
Sample 3	17.13	0.204	1.19

Method comparison:

A comparison of the Analyticon Turbitex CRP Latex (y) with a commercial obtainable assay (x) with 70 samples gave the following result:

$$y = 1.014 x - 0.094; \quad r = 0.999$$

Turbitex® CRP Latex

C-REACTIVE PROTEIN

Quality control:

Contronorm ARC	2 x 1 ml	#7562
Protein Control Level 1	3 x 1 ml	#7661

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The CRP method was standardized against the reference material CRM 470 (RPPHS - Reference Preparation for Proteins in Human Serum).

S1:	0,9% NaCl		
S2-S6:	Bio Cal® CRP ultra Calibration Set	5 x 1 ml	#14540

Calibration frequency

Full calibration is recommended:

- after reagent lot changes
- as required following quality control procedures

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