Order information:

Catalog No.	Contents						
7481	R1	3 x	20 ml	R2	1 x	10 ml	
H5301	Hit I (ILab*)	R1	6 x	19 ml	R2	1 x	16 ml
H5303	Hit 917 (AU*)	R1	6 x	60 ml	R2	6 x	9 ml
AU5303	AU	R1	6 x	60 ml	R2	6 x	9 ml

(*) Kit contains only reagent barcodes for Hitachi system.

System information:

Hitachi 911: ACN 066

Hitachi 917: ACN 109

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Immunoturbidimetric assay for the in vitro quantitative determination of CRP in human serum and plasma.

Summary:

C-reactive protein is the classic acute phase protein to inflammatory reactions. It is synthesized by the liver and consists of five identical polypeptide chains forming a five-membered ring of molar mass 120 000 daltons. CRP is the most sensitive of the acute phase reactants, and its concentration increases rapidly during inflammatory processes. Complexed CRP activates the complement system beginning with C1q. CRP then initiates opsonization and phagocytosis of invading cells, but its main function is to bind and detoxify endogenous toxic substances produced as a result of tissue damage.

CRP assays are used to detect systematic inflammatory processes (apart from certain types of inflammation such as SLE and colitis ulcerosa); to assess treatment of bacterial infections with antibiotics; to detect intrauterine infections with concomitant premature amniorrhexis; to differentiate between active and inactive forms of disease with concurrent infection, e.g. in patients suffering from SLE or Colitis ulcerosa; to therapeutically monitor rheumatic diseases and assess anti-inflammatory therapy; to determine the presence of postoperative complications at an early stage, such as infected wounds, thrombosis and pneumonia; and to distinguish between infection and bone marrow transplant rejection.

A variety of methods, such as nephelometry and turbidimetry are available for the determination of CRP. This CRP assay is based on the principle of immunological agglutination.

Test principle:

Immunoturbidimetric assay

Anti-CRP antibodies react with antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically. Addition of PEG allows the reaction to progress rapidly to the end point, increases sensitivity, and reduces the risk of samples containing excess antigen producing false negative results.

Reagent concentration:

R1:

TRIS/HCl buffer pH 7.5 100 mmol/l

Polyethylene glycol ether 1-5%

Preservative

R2:

Anti-human anti-CRP antibody (goat) dependent on titer

TRIS/HCl buffer* pH 7.6 50 mmol/l

NaCl 300 mmol/l

Preservative

*TRIS= Tris (hydroxymethyl)-aminomethane

Preparation and stability:

R1: Ready for use.

R2: Ready for use.

Unopened kit components: Up to the expiration date at +2°C to + 8°C

Onboard stability: R1: 90 days

R2: 90 days

Specimen:

Collect serum using standard sampling tubes. Li-heparin or EDTA-Plasma

Stability: 3 days at +20°C to +25°C

8 days at + 4°C to + 8°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to a bilirubin concentration of 50 mg/dl.

Hemolysis: No significant interference up to a haemoglobin concentration of 400 mg/dl.

Lipemia (Intralipid): No significant interference up to a triglyceride concentration of 1500 mg/dl.

There is poor correlation between turbidity and triglycerides concentration.

Rheumatoid factors up to 1200 IU/ml do not interfere.

A high-dose hook-effect may occur at CRP concentrations >50 mg/dl.

The effect of drugs or their metabolic products on the test is not known in all cases.

In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below

• 0.9% NaCl

Manual procedure:

Wavelength: 340 nm / Hg 334 nm
 Temperature: +37°C
 Cuvette: 1 cm
 Zero adjustment: against reagent blank

	Blank	Sample/Calibrator
Sample/Calibrator	---	40 µl
0.9% NaCl	40 µl	---
R1:	700 µl	700 µl

Mix, incubate 5 min. at +37°C, read absorbance A₁. Then add:

R2: 100 µl 100 µl

Mix, incubate 5 min. at +37°C. Then read absorbance A₂

Calculation:

$$\Delta A = [(A_2 - A_1) \text{ sample or calibrator}] - [(A_2 - A_1) \text{ blank}]$$

The concentration of CRP in patient sera has to be calculated from ΔA using mathematic function as logit/log or can be read from a graph using values of 6 levels of standards in the concentration range of 0 to 250 mg/dl CRP. For zero value is recommended to use saline solution (0.9%)

Measuring /reportable range:

0.3 - 25 mg/dl (3 - 250 mg/l)

Determine samples having higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl solution or distilled/deionized water (e.g. 1 + 2). Multiply the result by the appropriate dilution factor (e.g. 3).

Expected values:

Reference values according to CRM 470 Protein -Standardized
 < 0.5 mg/dl (<5.0 mg/L)

Elevated CRP range:

Concentrations below 1 mg/dl (10 mg/l) rule out numerous acute inflammatory diseases, but cannot be taken to exclude acute inflammatory processes in general. Elevated concentrations up to 5 mg/dl (50 mg/l) in acute illnesses occur in the presence of slight to moderate inflammatory processes. Values 5 above 5 mg/dl (50 mg/l) indicate high and extensive inflammatory activity.

It is important to monitor the CRP concentration during the acute phase of the illness.

The ranges indicated are intended merely for orientation purposes. They can only be interpreted in conjunction with the clinical findings.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the CRP results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 0.3mg/dl (3 mg/l)

The lower detection limit represents the lowest measurable CRP- concentration that can be distinguished from zero. It is calculated as three standard deviations of 21 replicates of the lowest standard.

Imprecision:

Reproducibility within run was determined using controls in an internal protocol (n = 20). The following results were obtained:

Within run			
Sample	Mean mg/l	SD mg/l	CV %
sample 1	18.65	0.47	2.55
sample 2	30.31	0.49	1.61
sample 3	102.07	0.96	0.94

Between day			
Sample	Mean mg/l	SD mg/l	CV %
sample 1	18.74	0.31	1.63
sample 2	33.02	0.65	1.98
sample 3	107.92	1.70	1.57

Method comparison:

A comparison of the Analyticon Turbitex CRP (y) with a commercial obtainable assay (x) with 102 samples gave the following result:

$$y = 1.09x - 1.09; \quad r = 0.997$$

Quality control:

Protein Control Level 1	3 x 1 ml	#7661
Protein Control Level 2	3 x 1 ml	#7662
Contronorm® ARC	2 x 1 ml	#7562
Controptath® ARC	2 x 1 ml	#7561

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The CRP method was standardized against CRM 470.

S1:	0.9% NaCl	
S2-S6:	Bio Cal® CRP Set	5 x 1 ml #14530
S2-S6:	CRP Standard	1 x 1 ml #8120

Factors for calculating the concentrations using the lot-assigned value of the CRP standard:

$$S2: 0.05; \quad S3: 0.10; \quad S4: 0.20; \quad S5: 1.00; \quad S6: 2.50$$

Calibration frequency:

Full calibration is recommended

- after lot change
- as required following quality control procedures

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