Order information:

Catalog No.	Contents
7440 50 Tests	R1 1 x 2.5 ml R2 1 x 0.5 ml R3 1 x 0.5 ml
7400 100 Tests	R1 1 x 5 ml R2 1 x 1 ml R3 1 x 1 ml
7140 Latex reagent	R1 1 x 2.5 ml
7100 Latex reagent	R1 1 x 5 ml
7111 CRPslide control	R2 1 x 0.5 ml R3 1 x 0.5 ml
7310 Reaction cards	1 x 20 pcs

Intended use:

A rapid slide test for the qualitative and semi-quantitative determination of c-reactive protein in human serum.

Summary:

C-reactive protein is the classic acute phase protein to inflammatory reactions. It is synthesized by the liver and consists of five identical polypeptide chains forming a five-membered ring of molar mass 120.000 daltons. CRP is the most sensitive of the acute phase reactants, and its concentration increases rapidly during inflammatory processes. Complexed CRP activates the complement system beginning with C1q. CRP then initiates opsonization and phagocytosis of invading cells, but its main function is to bind and detoxify endogenous toxic substances produced as a result of tissue damage. CRP assays are used to detect systematic inflammatory processes (apart from certain types of inflammation such as SLE and colitis ulcerosa); to assess treatment of bacterial infections with antibiotics; to detect intrauterine infections with concomitant premature amniorrhexis; to differentiate between active and inactive forms of disease with concurrent infection, e.g. in patients suffering from SLE or Colitis ulcerosa; to therapeutically monitor rheumatic diseases and assess anti-inflammatory therapy; to determine the presence of post-operative complications at an early stage, such as infected wounds, thrombosis and pneumonia; and to distinguish between infection and bone marrow transplant rejection. A variety of methods, such as nephelometry and turbidimetry are available for the determination of CRP. This CRP assay is based on the principle of immunological agglutination.

Principle:

The CRP reagent is a suspension of polystyrene latex particles coated with the gamma globulin fraction of antihuman CRP specific serum. When CRP is present in the sample, presence of agglutination indicates a content of CRP equal or greater than 6 mg/l, without previous sample dilution.

Reagent Concentration:

R1:

Latex particles coated with anti-CRP preservative

R2:

Protein containing solution with C-reactive protein preservative > 6 mg/l

R3:

Protein containing solution with C-reactive protein preservative < 6 mg/l

Preparation and stability

The reagents and control sera are stable: Up to the expiration date at +2°C to +8°C Do not freeze!

Specimen:

Collect serum using standard sampling tubes.

Use fresh clear serum specimens which are received from centrifuged coagulated blood. Plasma, lipemic serum or microbial contamination may cause erroneous results. If the test cannot be performed immediately, store the specimen at +2 to +8°C for up to 24 h.

For longer storage, freeze the serum.

Centrifuge samples containing precipitate before performing the assay.

The test does not work with whole blood.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

The positive and negative controls were prepared from human sera which have been tested using FDA-approved methods and found to be non-reactive for HbsAg and HIV antibodies. However as no test method can rule out the potential risk of infection with absolute certainty, the material should be treated just as carefully as a patient sample. In the event of exposure the directives of the responsible health authorities should be followed.

Precautions:

Reagents containing sodium acid may combine with copper and lead plumbing to form highly explosive metal azides. Dispose of reagent by flushing with large amounts of water to prevent azide build up.

Contaminated sera and a longer reaction time than 3 min. may cause false positive agglutination.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Procedure Qualitative Test:

1. Bring reagents and specimens to room temperature before use.
2. Resuspend the Latex Reagent gently to homogeneity.
3. Place one drop (40 µl) of each test specimen and if reference is necessary, one drop (40 µl) of the pos./neg. control on **separate fields** of the reaction slide.
4. Add one drop Latex Reagent (40 µl) to each test field. Use a fresh spatula dispenser for each field to spread reaction mixture over entire test field.
5. Rotate the slide for 2 minutes and read immediately under optimal light conditions.

Procedure Semi-Quantitative Test:

1. Bring reagents and specimens to room temperature before use.
2. Resuspend the Latex Reagent gently to homogeneity.
3. Using saline, dilute the specimens 1:2, 1:4, 1:8, 1:16 or as needed.
4. Place one drop (40 µl) of the test specimen followed by the corresponding prepared dilutions and as reference one drop (40 µl) of the pos./neg. control on **separate fields** of the reaction slide.
5. Add one drop of the Latex Reagent (40 µl) to each sample. Use a fresh spatula dispenser for each field to spread reaction mixture over entire test field.
6. Rotate the slide for 2 minutes and read immediately under optimal light conditions.

Results Qualitative Test:

A negative reaction is indicated by a uniform milky suspension with no agglutination as observed with the negative control (2). A positive reaction is indicated by any observable agglutination in the reaction mixture. The specimen reaction should be compared with the positive control (1).

Results Semi-Quantitative Test:

The CRP concentration of the serum lies between the highest dilution which just shows a positive reaction (e.g. 1:4) and the following dilution which shows a negative reaction (e.g. 1:8). For calculating the serum concentration, multiply the corresponding dilution factor by 6mg/l (cut off of the test).

e.g. Last positive signal 1:4 $4 \times 6 \text{ mg/l} = 24 \text{ mg/l}$

First negative signal 1:8 $8 \times 6 \text{ mg/l} = 48 \text{ mg/l}$

The CRP conc. of the sample is between 24 and 48mg/l.

Interpretation:

positive (1)

agglutination
within 2 minutes

negative (2)

no agglutination
within 2 minutes

Reference values:

Adults: < 6 mg/l

CRP is a protein which serum levels increase significantly in an acute phase status as it happens in the most forms of different tissue injuries. It was described the presence of such protein in normal persons at a titer not exceeding 6 mg/l. However the increase occurs in a non specific way in different aggressions like malignant tumor, rheumatic fever, myocardial infarct, inflammations, etc.

Each laboratory should investigate the transferability of the expected values to its own patients' population and if necessary determine its own reference range. For diagnostic purposes the CRP results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

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