BioLyzer® Order information:

Catalog No.	BioLyzer	Contents		
B3551	200 / 600	R1	3 x	100 ml
		R2	3 x	100 ml
B3553	300 / 600*	R1	6 x	20 ml
		R2	6 x	20 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

In vitro test for the quantitative determination of calcium in human serum, plasma and urine.

Summary:

The calcium content of an adult is somewhat over 1kg (25000mmol), i.e. about 2% of the body weight. Of this, 99% is present as calcium hydroxyl-apatite in bones and less than 1% is present in the extra-osseous ICS (intra-cellular space) or ECS (extracellular space). The calcium level in the ECS (approx. 100mmol) is in dynamic equilibrium with the rapidly exchangeable fraction of bone calcium. Calcium ions affect the contractility of the heart and the skeletal musculature and are essential for the function of the nervous system. In addition, calcium ions play an important role in blood clotting and bone mineralization. In plasma, calcium is bound to a considerable extent to proteins (approx. 40%), 10% is in the form of inorganic complexes and 50% is present as free (ionized) calcium. The body's calcium balance is regulated by the parathyroid hormone (PTH), calcitriol (CT) and calcitonin.


The test is used for the diagnosis and monitoring of hypocalcemia (calcium deficiency) and hypercalcemia (excess calcium) in serum. The characteristic symptom of hypocalcemia is latent or manifest tetany and osteomalacia. Hypocalcemia is due to the absence or impaired function of the parathyroid or impaired vitamin D-synthesis. Hypercalcemia is brought about by increased mobilization of calcium from the skeletal system (osteoporosis) or increased intestinal absorption. The majority of cases are due to primary hyperparathyroidism (pHPT) or bone metastasis of carcinoma of the breast, prostate or thyroid and bronchial carcinoma.

The main significance of determining urinary calcium lies in the differentiation between hypercalciuria and hypocalciuria and the differential diagnosis of nephrolithiasis.

Complexometric methods are used in addition to atomic absorption spectrometry (AAS) for determining calcium. The following calcium determination is based on the reaction of calcium with o-cresolphthalein complexone in alkaline solution. Magnesium is masked with 8-hydroxy-quinoline.

Test principle:

- Sample and addition of R1 (buffer)
- Addition of R2 (chromogen) and start of reaction:

The colour intensity of the purple complex formed is directly proportional to the calcium concentration and is measured photometrically.

Reagent Concentration:

R1: AMP buffer, pH 10.7	350 mmol/l
R2: o-Cresolphthalein complexone	0.16 mmol/l
8-hydroxyquinoline	6.9 mmol/l
Hydrochloric acid detergent	0.06 mmol/l

Preparation and stability:

R1: ready for use
R2: ready for use

Unopened kit components up to the expiration at + 2 to + 25°C.

Onboard stability:	R1	35 days
	R2	35 days

Specimen:

Collect serum using standard sampling tubes.
Heparin plasma. Plasma must be assayed fresh.
Stability: 7 days at +20°C to + 25°C
3 weeks at + 4°C to + 8°C
8 months at - 20°C

Do not use oxalate-, EDTA-, or citrate-plasma.

Urine

24-hour urine: Place 10 ml hydrochloric acid (6mol/l) in a collection bottle, or acidify (pH<2.0) after urine collection, in order to dissolve calcium salts. Before determination dilute urine 1:1 with distilled water, multiply the result by 2.

Stability :	2 days	at +20°C to + 25°C
	4 days	at + 4°C to + 8°C
	3 weeks	at -20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitation interference:

Criterion: Recovery within ± 10% of initial value in serum/plasma.

Icterus: No significant interference up to an index I of 100 (approximate bilirubin concentration: 100 mg/dl).

Hemolysis: No significant interference up to an index H of 1100 (approximate hemoglobin concentration: 1100 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 175 (approximate triglycerides concentration: 350 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

The effect of drugs or their metabolic products on the test is not known in all cases.

In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl
- Chimneys

Measuring range:

Serum/plasma

1 – 5 mmol/l (4 - 20mg/dl)

Determine samples having higher concentrations via the rerun function using 0.9 % NaCl.

Urine

1 – 5 mmol/l (4 - 20mg/dl)

Determine samples having higher concentrations via the rerun function using 0.9 % NaCl.

Reference value:

Serum/plasma

2.15-2.55 mmol/l (8.60-10.2 mg/dl)

24-hour urine

2.5-8.0 mmol/24 h (100-320 mg/24 h), corresponding to 1.7-5.3 mmol/l* (6.7-21 .3 mg/dl*)

*Assuming a urine volume of 1.5 l/24 h

Reference range acc. to Tietz

Serum/plasma

Children (0-10 days):	7.6-10.4 mg/dl (1.90-2.60 mmol/l)
Children (10 days-2 years):	9.0-11.0 mg/dl (2.25-2.75 mmol/l)
Children (2-12 years): ^	8.8-10.8 mg/dl (2.20-2.70 mmol/l)
Adults (12-60 years):	8.4-10.2 mg/dl (2.10-2.55 mmol/l)
Adults (60-90 years):	8.8-10.2 mg/dl (2.20-2.55 mmol/l)
Adults (> 90 years):	8.2-9.6 mg/dl (2.05-2.40 mmol/l)

Urine: 100-300 mg/24 h (with normal food intake)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the calcium results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 0.01 mmol/l (0.04 mg/dl)

The detection limit represents the lowest measurable calcium concentration that can be distinguished from zero.

Imprecision:

Serum

Reproducibility was determined using controls (n= 20). The following results were obtained.

Sample	Run to run		
	MW (mmol/l)	SD (mmol/l)	VK %
Control serum 1	1.69	0.165	9.7
Control serum 2	2.95	0.27	9.1
Control serum 3	2.94	0.225	7.7

Sample	Within run		
	MW (mmol/l)	SD (mmol/l)	VK %
Control serum 1	2.09	0.034	1.6
Control serum 2	3.50	0.073	2.1
Control serum 3	2.59	0.088	3.4

Method comparison:

A comparison of the Analyticon Fluitest[®] Ca-CPC (y) with a commercial obtainable assay (x) gave the following result:

$$y = 1.0223 x + 0.0697; \quad r = 0.9988$$

Quality control:

Human Control Serum:

Contronorm [®] Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Contropath [®] Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Human Urine Control:

Urine Control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: The calcium method was standardized by atomic absorption spectrometry.

Calibration Type: Linear

S1: 0.9% NaCl

S2: Bio-Cal[®] E 10 x 3 ml #1430

Calibration frequency:

Two-point calibration is recommended:

- After reagent lot change
- As required following quality control procedures

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