REF: CG0305 10 x 2 ml
CG0306 10 x 4 ml

Coagnos® APTT

GB

INTENDED PURPOSE:

The Coagnos® APTT kit is intended for carrying out clot based haemostasis assays in human plasma. The Coagnos® APTT kit is indicated for professional use only.

For use in the determination of activated partial thromboplastin times (aPTT), and related coagulation procedures using phospholipid extract and a near-colloidal particle activator. The test system can be used on manual, semi-automated and automated methods. From its origins through the work of Langdell and coworkers¹, and later modified by Proctor and Rapaport², the aPTT is used to detect disorders in the intrinsic coagulation system, which involves coagulation factors VIII, IX, XI, XII, prekallikrein, and high molecular weight kininogen. The aPTT is also used in assays which quantitate these factors and is routinely used for presurgical screening and monitoring of heparin therapy³. Commercially available reagents typically use one of three activators: kaolin, silica, or ellagic acid. In the basic screening test, the aPTT indirectly measures the formation of thrombin by its action on fibrinogen forming the fibrin clot. In the test, citrated test plasma is mixed with aPTT reagent for a specified period of time (typically 5 minutes) at +37°C followed by the addition of pre-warmed (+37°C) calcium chloride (0.025 M). Timing is begun from the time of addition of calcium chloride. The time required for clot formation is the aPTT. Clot detection can be by mechanical, manual (tilt tube), or photo-optical measurement.

WARNINGS AND PRECAUTIONS:

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. The material safety data sheet is available for download from our homepage <http://www.analyticon-diagnostics.com>. Dispose of components in accordance with local regulations. In case any serious incident has occurred in relation to the device, please report to the manufacturer and, if applicable, to the competent authority of the country in which the users and/or the patients established themselves.

COMPOSITION:

Component	Content	Description	Preparation
APTT	10 x 2 mL (REF CG0305) 10 x 4 mL (REF CG0306)	Reagent contains a near colloidal particle activator (magnesium-aluminium-silicate) for optimum sensitivity to factor deficiencies and to heparin. The reagent also contains phospholipids with buffer and stabilisers.	Bring to room temperature prior to use. Mix well by swirling or inversion prior to use.
Calcium Chloride	10 x 2 mL (REF CG0305) 10 x 4 mL (REF CG0306)	The reagent is a 0.025 M solution of calcium chloride.	The reagent is ready for use as packaged.

STORAGE, SHELF-LIFE AND STABILITY:

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Store at +2° – +8°C. DO NOT FREEZE. Stable for 30 days after opening. Avoid prolonged heating.

SAMPLE COLLECTION AND PREPARATION:

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at +2° – +8°C or +18 – +24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at +37°C prior to testing. Do not keep at +37°C for more than 5 minutes. This will minimise the neutralisation of the lupus inhibitor⁴. Erroneous results may be caused by contamination with tissue fluids or stasis. Avoid agitation, air bubbles or foaming. For the effects of commonly administered drugs, refer to Young, *et al.*⁵.

PROCEDURE:**Manual Method**

1. Prewarm well mixed APTT and 0.025 M Calcium Chloride to +37°C.
2. Prewarm 0.1 mL test plasma in duplicate to +37°C for 2 minutes.
3. Forcibly add 0.1 mL prewarmed APTT to plasma and start timer. Incubate for exactly 5 minutes at +37°C.
4. Add 0.1 mL prewarmed 0.025 M Calcium Chloride.
5. Note time for clot formation. Report result as aPTT time (seconds).

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Analyticon Biotechnologies AG for instrument specific application guides.

INTERPRETATION OF RESULTS:

The results of the aPTT test should be reported to the nearest 1/10 of a second. The normal range (usually $X \pm 2$ standard deviations) for each individual laboratory should be established. Results greater than the upper limits of the normal range should be considered abnormal and follow-up testing should be performed. Any aPTT values less than the lower limits of the normal range should be repeated on a new blood sample. Short aPTT values may be seen in association with *in vivo* thrombosis (e.g. deep vein thrombosis and disseminated intravascular coagulation).

Heparin monitoring

When monitoring heparin therapy, it is important to construct an *in vitro* reference curve which reflects the average heparin response, since individual patients respond differently to heparin. In general, one can consider the therapeutic range for heparin to be 0.2 to 0.5 units/mL^{3,6}. The following precautions should be considered when monitoring heparin therapy:

1. Time of collection is important, since heparin has an *in vivo* half-life of only 1.5 hours⁷
2. Release of platelet factor 4 (heparin neutralising factor) caused by platelet aggregation or damage during collection, should be avoided. Careful blood collection, proper centrifugation and prompt removal of the platelet poor plasma from the cells will help minimise the release of platelet factor
3. Baseline data on each patient's aPTT should be established before therapy, to determine the respective patient aPTT as it relates to the normal range established by the laboratory.
4. Different clot detection systems (mechanical, photo-optical, etc.) show variable sensitivities to heparin. The same test system should be used when monitoring heparinised patients.
5. Heparin response curves should be reestablished when lot numbers of reagent change and at periodic intervals with the same lot number.
6. The curve should also be constructed using the same heparin employed in therapy, to eliminate variables connected with heparins from different sources (e.g. porcine mucosa or bovine lung).

LIMITATIONS:

Expected values for the aPTT test will vary from one laboratory to another, depending on the technique used. The method of clot detection, temperature, pH, collection technique, type of anticoagulant and time and method of specimen storage are all very important. Plasma sample collection and storage conditions should be standardised and carefully controlled. Unexpected results should be confirmed by additional tests. Platelet fragments present in a specimen may cause the release of phospholipids, and thus the neutralisation of any lupus inhibitor present in the specimen. The use of specimens with small plasma volumes should be avoided due to possible physiological pH changes. Testing could be affected by several drugs⁵. An increase in the aPTT results may be caused by the administration of diphenylhydantoin, heparin, warfarin and radiographic agents^{5,8}. Decreased aPTT values may be seen during the use of oral contraceptives, or male estrogen therapy^{9,10}. Thus, laboratories should establish their own expected values for patients and well defined performance standards for the control.

QUALITY CONTROL:

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Analyticon Biotechnologies AG supplies the following controls available for use with this product:

REF CG0361 Coagnos® Special Control N
REF CG0362 Coagnos® Special Control P
REF CG0351 Coagnos® Routine Control N
REF CG0352 Coagnos® Routine Control P

REFERENCE VALUES:

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges. Using the Coagulyzer® Auto Pro, normal values ranging from 25.4 – 38.9 seconds; 0.78 – 1.20 (Ratio) are typical.

PERFORMANCE CHARACTERISTICS:

The following performance characteristics have been determined by Analyticon Biotechnologies AG or their representatives using a Coagulyzer® Auto Pro instrument. Each laboratory should establish its own performance data.

Precision

Precision was performed according to CLSI guideline EP05-A3: Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition, 2008.

Sample	Mean	Repeatability		Within-device	
		SD	CV [%]	SD	CV [%]
Clots [s]	33.1	0.66	1.98	1.22	3.69
	54.4	1.53	2.81	2.57	4.73
	69.8	1.64	2.35	3.28	4.70
Ratio	1.16	0.02	2.03	0.07	6.03
	1.91	0.05	2.86	0.15	7.69
	2.45	0.06	2.43	0.19	7.76

Interferences

Coagnos® APTT is insensitive to Heparin levels of up to 0.03 U/ml. Using a 5% interference threshold, there is no significant interference from Haemoglobin at concentrations up to 500 mg/dL. Using a 5% interference threshold, there is no significant interference from Bilirubin at concentrations up to 60 mg/dL for Coagnos® APTT. Lipid interference testing demonstrates that lipid levels do not directly affect the clot time of the reagent up to 250 mg/dL. Lipid concentrations in excess of this prevent clot detection.


REF: CG0305 10 x 2 ml
CG0306 10 x 4 ml

Method comparison

Comparison of clot time in seconds were determined using Coagnos® APTT on Coagulyzer® Auto Pro and Siemens® APTT with Dade® Actin® FS on a Sysmex® CS-5100 system on 104 samples. The following correlations were obtained:

Coagnos® APTT (Seconds) = 0.842x + 9.059 $r^2 = 0.74$ $n = 104$

SYMBOLS

	In vitro diagnostics product
	The product complies with European legislation
	Follow the instructions for use !
	Use by
	Permitted storage temperature range
	Read warnings and precautions in instructions for use !
	Batch identification number
REF	Item number
	Manufacturer

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BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/ LITERATURVERZEICHNIS:

1. Langdell R, Wagner R, Brinkhous K (1953) Effect of Antihemophilic Factor on One Stage Clotting Tests, J. Lab. Clin. Med, 41: 637
2. Proctor R, Rapaport S (1961) The Partial Thromboplastin Time with Kaolin, Am. J. Clin. Path, 36: 212
3. Brandt JT and Triplett DA (1981) Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the Activated Partial Thromboplastin Time, AJCP, 76: 530-537
4. Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
5. Young DS et al. Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACC Press, Washington, D.C., 1990
6. Kelton JG and Hirsh J (1980) Bleeding Associated with Antithrombotic Therapy, Semin Hematol, 17: 259-91
7. Estes J and Poulin P (1974) Pharmacokinetics of Heparin, Thromb Diath Haemorr, 33: 26-35
8. Solomon G, Hilgartner M, Kutt H (1972) Coagulation defects caused by diphenylhydantoin, Neurology, 22: 1165-71
9. Ambrus JL, Schimert G, Lajos TZ, Ambrus CM, Mink IB, Lassman HB, Moore RH, Melzer J (1971) Effect of antifibrinolytic agents and estrogens on blood loss and coagulation factors during open heart surgery, J Med, 2: 65-81
10. Crowell EB Jr, Clatanoff DV, Kiekhofers W (1971) The effect of oral contraceptives on factor VIII levels, Lab. Clin. Med, 77: 551-7