Coagnos® Antithrombin Xa

GB

INTENDED PURPOSE:

The Coagnos® Antithrombin Xa kit is intended for carrying out chromogenic based haemostasis assays in human plasma. The Coagnos® Antithrombin Xa kit is indicated for professional use only.

AT-III is a major inhibitor of blood coagulation, acting to inhibit plasma serine proteases including factors IXa, Xa, XIa, and thrombin. The rate of the inhibition is significantly increased in the presence of heparin. AT-III deficiency may be congenital or acquired and is associated with an increased risk for thrombosis¹⁻³. Deficiencies may occur in liver disease, disseminated intravascular coagulation (DIC), septicemia, pulmonary embolism, nephrotic syndrome, stroke, and thrombophlebitis^{2,3}. Most functional assays for AT-III rely on the ability of plasma to inactivate thrombin in the presence of heparin, such as the assay described by Odegard et al⁴. More recently, it has been shown that a method based on the ability of plasma to inhibit factor Xa in the presence of heparin may be more accurate, since it eliminates interference due to heparin cofactor II⁵ which does not inhibit factor Xa⁶. In the present two-stage method, factor Xa is added to a plasma dilution containing AT-III in the presence of excess heparin and calcium. After an initial incubation period (stage 1), residual factor Xa is determined with a factor Xa-specific chromogenic substrate (stage 2). The residual factor Xa activity is inversely proportional to the AT-III concentration. The Analyticon Biotechnologies AG Coagnos® Antithrombin Xa kit is intended for the quantitative determination of antithrombin III (AT-III) activity in human citrated plasma by chromogenic assay.

WARNINGS AND PRECAUTIONS:

The reagents contained in this kit are for in vitro diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. The material safety data sheet is available for download from our homepage <http://www.analyticon-diagnostics.com>. Dispose of components in accordance with local regulations. In case any serious incident has occurred in relation to the device, please report to the manufacturer and, if applicable, to the competent authority of the country in which the users and/or the patients established themselves.

COMPOSITION:

Component	Content	Description	Preparation
Factor Xa	5 x 2 mL	Contains freeze-dried bovine factor Xa.	Reconstitute each vial with 2 mL of distilled / deionised water. Allow to stand for 10 minutes and mix well before use.
Factor Xa Substrate	5 x 2 mL	Contains freeze-dried CH3OCO-D-CHAGly-Arg-pNA•AcOH.	Reconstitute each vial with 2 mL of distilled / deionised water. Allow to stand for 10 minutes and mix well before use.
Sample Diluent	5 x 3 mL	Contains 5x buffer concentrate with sodium azide at 0.1% as preservative. When fully diluted, buffer contains 0.05 M Tris-HCl, 0.175 M NaCl, 7.5 mM Na ₂ EDTA and sodium heparin at pH 8.4.	Dilute 1 + 4 with distilled / deionised water and mix well.

ITEMS REQUIRED BUT NOT PROVIDED:

Coagulation instrument or spectrophotometer operable at 405 nm.

Depending on Assay Method Used:

- REF CG0370 Coagnos® Calibration Plasma
- Glacial acetic acid
- +37°C water bath or dry bath
- Laboratory timer

STORAGE, SHELF-LIFE AND STABILITY:

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Factor Xa	Reconstituted reagent is stable for 2 months at +2 →+8°C, 1 month at +17°C or 2 days at +37°C.
Factor Xa Substrate	Reconstituted reagent is stable for 2 months at +2 →+8°C, 1 month at +17°C or 7 days at +37°C.
Sample Diluent	Store diluted buffer in a tightly stoppered bottle at +2 →+8°C and use within one month.

SAMPLE COLLECTION AND PREPARATION:

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at +2 →+8 °C or +18 →+24 °C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20 °C for 2 weeks or -70 °C for 6 months. Thaw quickly at +37 °C prior to testing. Do not keep at +37 °C for more than 5 minutes⁷.

PROCEDURE:

Preparation of Standard (REF CG0370), Control and Patient Dilutions:

%AT-III	Plasma	Dilution Buffer
100%	10 µL Standard	990 µL
50%	500 µL of 100% Std	500 µL
25%	500 µL of 50% Std	500 µL
12.5%	500 µL of 25% Std	500 µL
Patient or Control	10 µL Plasma	990 µL

A. Semi-Micro Endpoint Method

- Add 200 µL standard / control or patient plasma dilution to a test tube.
- Incubate at +37°C for 2-4 minutes.
- Add 200 µL Factor Xa Reagent and mix.
- Incubate at +37°C for EXACTLY 1 minute.
- Add 200 µL Factor Xa Substrate and mix.
- Incubate at +37°C for EXACTLY 3 minutes.
- Add 200 µL acetic acid (50%) and mix.
- * Add 200 µL water (optional).

The yellow colour of the final reaction product is stable for at least 4 hours. Read absorbance at 405 nm in a 1 cm semi-micro cuvette against a blank prepared in the following order:

- 200 µL acetic acid.
- 200 µL standard dilution.
- 200 µL Factor Xa reagent.
- 200 µL Factor Xa Substrate.
- * 200 µL water (optional).

* Some spectrophotometers require a minimum of 1 mL volume in the cuvette. If the patient plasma is very icteric, a second blank containing patient plasma dilution instead of standard dilution should be prepared and the absorbance subtracted from the absorbance obtained for the patient's AT-III determination. As many as ten determinations can be performed simultaneously with the same stop watch by staggering pipetting steps at five second intervals.

B. Kinetic Method

A kinetic analyser may be used to measure the initial rate of hydrolysis of the chromogenic substrate. The procedure to be used is as follows:

To a reaction cuvette:

- Add 200 µL diluted standard or patient plasma.
- Incubate at +37°C for 1 minutes.
- Add 200 µL Factor Xa reagent and mix.
- Incubate at +37°C for 1 minute.
- Add 200 µL Factor Xa Substrate.
- Measure rate of change of absorbance at 405 nm for 1 minute.

INTERPRETATION OF RESULTS:

a) Assay Calibration

A commercially prepared plasma standard in which AT-III has been determined eg. Analyticon Biotechnologies AG Coagnos® Calibration Plasma (REF CG0370) should be used as a reference.

b) Calibration Curve

Plot the absorbance obtained with each of the AT-III calibration standards on the y-axis against % AT-III on the x-axis using linear graph paper. The line of best fit should be determined by linear regression analysis. The AT-III in plasma samples can be determined by interpolation from the calibration curve. The AT-III concentration in the patient specimen should be adjusted for the AT-III concentration in the standard:

$$\% \text{ AT-III (adjusted) } = \% \text{ AT-III (patient) } \times \% \text{ AT-III (standard) } / 100$$

LIMITATIONS:

Avoid icteric, lipaemic, and haemolysed samples. To ensure reproducible results, use accurate pipetting devices and observe recommended procedures with emphasis on incubation times and temperature.

QUALITY CONTROL:

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Analyticon Biotechnologies AG supply the following controls available for use with this product:

REF CG0351	Coagnos® Routine Control N
REF CG0352	Coagnos® Routine Control P
REF CG0361	Coagnos® Special Control N
REF CG0362	Coagnos® Special Control P

REFERENCE VALUES:

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range. The normal range of AT-III is 75-125% in plasma. Plasma activity levels of 30-60% can be observed in patients with hereditary AT-III deficiency^{1,2}. Several clinical conditions associated with acquired AT-III deficiency include: liver disease, DIC, nephrotic syndrome, pulmonary embolism, stroke, and thrombophlebitis. In addition, oral contraceptive use may lead to reduced AT-III levels².

REF: CG0381

PERFORMANCE CHARACTERISTICS:

The following performance characteristics have been determined by Analyticon Biotechnologies AG or their representatives using a Coagulyzer® Auto Pro instrument. Each laboratory should establish its own performance data.

Precision

Precision was performed according to CLSI guideline EP05-A3: Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition, 2008.

Sample	Mean	Repeatability		Within-device	
		SD	CV [%]	SD	CV [%]
Antithrombin Xa Concentration [%]	98.6	6.19	6.28	7.06	7.16
	41.5	4.22	10.18	6.97	16.81
	34.4	4.03	11.73	6.26	18.21

Linearity

The method is designed to give a linear standard curve from 2.9 – 134.8 % Antithrombin Xa.

Interferences


Coagnos® Antithrombin Xa is insensitive to Heparin levels of up to 3 U/ml. Using a 5% interference threshold, there is no significant interference from Haemoglobin at concentrations up to 1000 mg/dL. Using a 5% interference threshold, there is no significant interference from Bilirubin at concentrations up to 29 mg/dL for Coagnos® Antithrombin Xa. Lipid interference testing demonstrates that lipid levels do not directly affect the clot time of the reagent up to 672 mg/dL.

Method comparison

Comparison of AT-III concentration [%] were determined using Coagnos® Antithrombin Xa on Coagulyzer® Auto Pro and Siemens® Antithrombin with INNOVANCE® Antithrombin on a Sysmex® CS-5100 system on 104 samples. The following correlations were obtained:

$$\text{Coagnos® Antithrombin Xa [\%]} = 0.885x - 5.241 \quad r^2 = 0.98 \quad n = 104$$

SYMBOLS

	In vitro diagnostics product
	The product complies with European legislation
	Follow the instructions for use !
	Use by
	Permitted storage temperature range
	Read warnings and precautions in instructions for use !
	Batch identification number
	Item number
	Manufacturer

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Hersteller

BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/ LITERATURVERZEICHNIS:

1. Hirsh J et al (1989) Congenital antithrombin III deficiency. Incidence and clinical features, Am J. Med, 87(suppl. 3B): 34S-38S
2. Panicucci F et al (1981) Antithrombin III, heparin cofactor and antifactor Xa in relation to age, sex and pathological condition, Haemostasis, 9: 297-302
3. Conlan MG et al (1994) Antithrombin: Antithrombin III: Associations with age, race, sex and cardiovascular disease risk factors, Thromb. Haemost, 72(4): 551-556
4. Odegard OR et al (1975) Heparin cofactor activity measured with an amidolytic method, Thromb. Res, 6: 287-294
5. Demers C et al (1993) An antithrombin III assay based on factor Xa inhibition provides a more reliable test to identify congenital antithrombin III deficiency than an assay based on thrombin inhibition, Thromb. Haemost, 69: 231-235
6. Tolefsen DM, Blank MK (1981) Detection of a new heparin-dependent inhibitor of thrombin in human plasma, J. Clin. Invest, 68: 589-596
7. Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5