Coagnos® Calibration Plasma

REF: CG0370 10 x 1 ml

Coagnos® Calibration Plasma

GB

INTENDED PURPOSE:

The Coagnos® Calibration Plasma kit is intended for use as a calibration material. The Coagnos® Calibration Plasma kit is indicated for professional use only. Calibration Plasma may be used as a reference plasma when assaying for Factors II, V, VII, VIII, IX, X¹, XI, XII², fibrinogen³, von Willebrand Factor, antigenic and functional Protein C and Protein S (total and free), as well as the chromogenic assays including Antithrombin (AT) III, Protein C, Factor VIII and Plasminogen. Calibration Plasma is used in factor assays and other testing in the same manner as a fresh plasma pool. Factor II, VII, VIII, IX and X values and the chromogenic Factor VIII, AT III and Protein C values are traceable to World Health Organisation standards to ensure the utmost credibility in values stated⁴. The reference plasma should be used to gauge internal factors in each laboratory's system. The reference plasma should not be used to determine normal ranges since normals vary from population to population.

WARNINGS AND PRECAUTIONS:

The reagents contained in this kit are for in vitro diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. The material safety data sheet is available for download from our homepage <http://www.analyticon-diagnostics.com>. Dispose of components in accordance with local regulations. Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of: Hepatitis B Antigen (HbsAg), HIV 1 antibody, HIV 2 antibody, HCV antibody. However, they should be handled with the same precautions as a human patient sample. In case any serious incident has occurred in relation to the device, please report to the manufacturer and, if applicable, to the competent authority of the country in which the users and/or the patients established themselves.

COMPOSITION:

Component	Content	Description	Preparation
Calibration Plasma	10 x 1 mL	Calibration Plasma is prepared from a frozen pool of citrated plasma from healthy donors and is buffered and lyophilised to ensure stability of all plasma constituents ⁵ .	Reconstitute with 1 mL of distilled or deionised water. Swirl gently. Allow product 20 minutes for complete dissolution.

STORAGE, SHELF-LIFE AND STABILITY:

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Values for Factor VIII, von Willebrand factor and ristocetin co-factor are stable for 2 hours at +2 –+8°C. All other factors are stable for 4 hours at +2 –+8°C. The unreconstituted Calibration Plasma must appear as a light yellow, dry plug. If any unusual conditions are noted, the customer should notify Analyticon Biotechnologies AG, before using.

SAMPLE COLLECTION AND PREPARATION:

Not applicable.

PROCEDURE:

Refer to the appropriate instrument Operator Manual for detailed instructions or contact Analyticon Biotechnologies AG for instrument-specific application guides.

INTERPRETATION OF RESULTS:

The assayed percent activities of the various coagulation factors should be taken from the "Reference Value" column when Coagnos® Calibration Plasma is used to determine the laboratory's standard curves. Ensure the lot number printed on this assay sheet is the same as that on the vial of Calibration Plasma to be used.

LIMITATIONS:

The results obtained with Coagnos® Calibration Plasma depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations^{6,7,8}. Each laboratory should establish an expected range for the particular instrument-reagent system.

QUALITY CONTROL:

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Analyticon Biotechnologies supplies the following controls available for use with this product:

REF CG0361 Coagnos® Special Control N
REF CG0362 Coagnos® Special Control P

REFERENCE VALUES:

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

SYMBOLS

	In vitro diagnostics product
	The product complies with European legislation
	Follow the instructions for use!
	Use by
	Permitted storage temperature range
	Read warnings and precautions in instructions for use!
	Batch identification number
REF	Item number
	Manufacturer

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BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/ LITERATURVERZEICHNIS:

1. Babson AL and Flanagan ML (1975) Quantitative One Stage Assays for Factors V and X, AJCP, 64: 817-819.
2. Hardisty RM et al. (1962) A One Stage Factor VIII Assay and Its Use on Venous and Capillary Plasma. Thrombosis et Diathesis Haemorrhagica, 7:215-229.
3. Morse EE et al. (1971) Automated Fibrinogen Determination, AJCP, 55:671-676.
4. Elodi S et al (1978) Some Sources of Error in the One-Stage Assay of Factor VIII, Haemostasis, 7:1-9.
5. Thelin M (1968) Preparation and Standardization of a Stable AHF Plasma, Thrombosis et Diathesis Haemorrhagica, 19:423.
6. Kirkwood TBL et al. (1977) Identification of Sources of Variation in Factor VIII Assay, British Journal of Haematology, 37:559-568.
7. Goldenfarb MD (1971) Reproducibility in Coagulation Assays, AJCP, 55:561-564.
8. Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, AJCP, 59:231-235.