Coagnos® D-Dimer Control H/L **GB**

INTENDED PURPOSE:

The Coagnos® D-Dimer Control H/L kit is intended for use as a quality control material. The Coagnos® D-Dimer Control H/L kit is indicated for professional use only. D-dimer containing moieties are formed by plasmin degradation of Factor XIIIa cross-linked fibrin. Elevated levels of D-dimer are found in clinical conditions such as deep vein thrombosis (DVT), pulmonary embolism (PE) and disseminated intravascular coagulation (DIC)^{1,2}. D-dimer levels rise during pregnancy and high levels are associated with complications³. The negative predictive value of D-dimer for thrombosis is high⁴, although a negative D-dimer does not completely rule out thrombosis. The Analyticon Biotechnologies AG Coagnos® D-Dimer Control H/L kit contains plasmas with low and high levels of D-dimer. These plasmas are for use with the Analyticon Biotechnologies AG latex-enhanced immunoassay methods for D-dimer and should be tested before patient samples to ensure the satisfactory performance of the assay system.

WARNINGS AND PRECAUTIONS:

The reagents contained in this kit are for in vitro diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. The material safety data sheet is available for download from our homepage <http://www.analyticon-diagnostics.com>. Dispose of components in accordance with local regulations. Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of: Hepatitis B Antigen (HbsAg), HIV 1 antibody, HIV 2 antibody, HCV antibody. However they should be handled with the same precautions as a human patient sample. In case any serious incident has occurred in relation to the device, please report to the manufacturer and, if applicable, to the competent authority of the country in which the users and/or the patients established themselves.

COMPOSITION:

Component	Content	Description	Preparation
D-Dimer Control - L	5 x 1 mL	Each vial contains lyophilized D-dimer enriched human plasma.	Reconstitute each vial with 1 mL of saline solution. Allow to stand for 5 minutes and mix well before use. Do not shake.
D-Dimer Control - H	5 x 1 mL	Each vial contains lyophilized D-dimer enriched human plasma.	Reconstitute each vial with 1 mL of saline solution. Allow to stand for 5 minutes and mix well before use. Do not shake.

STORAGE, SHELF-LIFE AND STABILITY:

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Reconstituted vials are stable for 5 days at +2 –+8 °C, or 3 months at -20 °C. Do not freeze / thaw more than once.

SAMPLE COLLECTION AND PREPARATION:

Not applicable.

PROCEDURE:

Each control should be treated in the same manner as the unknown specimen in accordance with the instructions outlined in each particular test protocol.

INTERPRETATION OF RESULTS:

Lot specific expected values are provided with each pack of controls. Recalibration may be required if results fall outside quoted ranges.

LIMITATIONS:

The values listed are for use with Analyticon Biotechnologies AG D-dimer procedures only. The use of other D-dimer procedures and products may lead to erroneous results.

QUALITY CONTROL:

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

REFERENCE VALUES:

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.


PERFORMANCE CHARACTERISTICS:

The following performance characteristics have been determined by Analyticon Biotechnologies AG or their representatives using an opto-mechanical coagulation instrument. Each laboratory should establish its own performance data.

Precision

	D-Dimer Control - H			D-Dimer Control - L		
	Mean	SD	CV [%]	Mean	SD	CV [%]
D-Dimer Concentration [µg/L]	1349.8	61.36	4.55	390.1	12.77	3.28

SYMBOLS

-  In vitro diagnostics product
-  The product complies with European legislation
-  Follow the instructions for use !
-  Use by
-  Permitted storage temperature range
-  Read warnings and precautions in instructions for use !
-  Batch identification number
- REF** Item number
-  Manufacturer

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BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/ LITERATURVERZEICHNIS:

1. Elms MJ et al (1986) Rapid Detection of Cross-Linked Fibrin Degradation Products in Plasma Using Monoclonal Antibody-Coated Latex Particles, J. Clin. Pathol, 85: 360-64.
2. Declerck PV et al (1987) Fibrinolytic Response and Fibrin Fragment D-dimer in Patients With Deep Vein Thrombosis, Thromb. Haemost, 58: 1025-9.
3. Ballegeer V et al (1987) Fibrinolytic Response to Venous Occlusion and Fibrin Fragment D-dimer Levels in Normal and Complicated Pregnancy, Thromb. Haemost, 58: 1030-2.
4. Hansson PO et al (1994) Can Laboratory Testing Improve Screening Strategies for Deep Vein Thrombosis at an Emergency Unit?, J. Intern. Med, 235: 143-51