Coagnos® D-Dimer **GB**

INTENDED PURPOSE:

The Coagnos® D-Dimer kit is a polystyrene, microparticle, agglutination assay for the quantitative determination of fibrin degradation products containing D-dimer in citrated human plasma using automated coagulation analyzers. The Coagnos® D-Dimer kit is indicated for professional use only.

SUMMARY:

Activation of the hemostatic system produces stabilized or cross-linked fibrin as its end product. The hydrolysis of fibrinogen molecules to form fibrin opens sites which allow the formation of affinity driven ternary complexes of fibrin, tissue plasminogen activator, and plasminogen, resulting in the formation of plasmin. Plasmin will cleave the cross-linked fibrin to form fibrin degradation products (FDP). Among these degradation products are species which contain factor XIII mediated transglutaminase (lysine-glutamine) linkages between the carboxy terminal regions of adjacent fibrin(ogen) gamma chains. The globular domains containing these gamma chains are referred to as D-domains and fibrin degradation products containing adjacent D-domains crosslinked by factor XIII are called D-dimers. Antibodies have been developed which recognize epitopes for D-dimer.^{1,2} Higher than normal levels of D-dimer in plasma would appear to indicate the activation of coagulation with the formation of crosslinked fibrin and the subsequent proteolysis to release D-dimer containing breakdown products. Elevated levels have been observed in disseminated intravascular coagulation (DIC),¹ deep venous thrombosis (DVT),^{2,3,4} pulmonary embolism (PE)⁵ and in other disease states such as cancer.⁶ D-dimer levels may rise in a normal pregnancy, but higher levels may indicate complications.⁷ Therefore, the presence of elevated D-dimer levels is not sufficient for the diagnosis of a thrombotic disorder. However, the absence of elevated D-dimer levels may be used to help rule out the presence of a thrombotic disorder, such as DIC, DVT or PE. In the Coagnos® D-Dimer assay, the polystyrene microparticles have monoclonal antibody attached covalently. In the presence of multimeric D-dimer antigen a complex is formed. The formation of the complexes causes an increase in the absorbance of the reaction mixture that is correlated to the antigen level. The Analyticon Biotechnologies AG procedure provides a means to quantitatively measure plasma levels of D-dimer. The D-dimer level may be used to determine patient management strategies.

WARNINGS AND PRECAUTIONS:

Coagnos® D-Dimer reagents are for "in vitro diagnostic use", Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing local, state and federal laws.

Coagnos® D-Dimer Reaction Buffer is an IRRITANT and HARMFUL. Irritating to eyes, respiratory system and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing.

Coagnos® D-Dimer Reagent and Reaction Buffer contain sodium azide, which may react with lead and copper plumbing to form highly explosive metal azides. Avoid azide accumulation.

Coagnos® D-Dimer Reaction Buffer and Calibrator are POTENTIALLY BIOHAZARDOUS MATERIALS. Source materials from which these products are derived were found negative for HBsAg, and for antibodies against HCV, HIV-1 and HIV-2, by approved test methods. Since no test method can offer complete assurance that infectious agents are absent, this product should be handled observing the same safety precautions employed when handling any potentially infectious material.

The material safety data sheet is available for download from our homepage <http://www.analyticon-diagnostics.com>. In case any serious incident has occurred in relation to the device, please report to the manufacturer and, if applicable, to the competent authority of the country in which the users and/or the patients established themselves.

COMPOSITION:

Important: The reagents are lot-specific. Lots are not interchangeable. Immediately before each run, gently agitate the Coagnos® D-Dimer Reagent by repeatedly inverting the vial to disperse the micro particles evenly. DO NOT SHAKE.

Component	Content	Description	Preparation
D-Dimer Latex	2 x 3 mL	Buffered solution of polystyrene microparticle beads coated with mouse monoclonal antibody MA-8D3 directed against D-dimer. Contains stabilizers and 0.095% sodium azide as a preservative.	Ready to use. The latex may sediment during storage. Mix thoroughly before use.
D-Dimer Buffer	2 x 7 mL	Buffer with stabilizers, pH 7.0. Sodium azide, 0.095%, added as preservative.	Ready to use.
D-Dimer Diluent	1 x 7 mL	0.85% saline solution.	Ready to use.
D-Dimer Calibrator	1 x 1 mL	Lyophilized human plasma enriched with fibrin D-dimer. 4624 ng/ml	Reconstitute Coagnos® D-Dimer Calibrator with 1.0 ml of deionized water. Restopper vial, agitate gently for about five (5) minutes, until a homogenous

			solution has been obtained. Allow to stand for 30 minutes at room temperature (18–26 °C). Mix gently, prior to use. DO NOT SHAKE. DO NOT USE STIR BAR.
--	--	--	--

ITEMS REQUIRED BUT NOT PROVIDED:

Coagulation analyzer capable of measuring absorbance at 405nm. Analyzer specific consumables. Analyticon Biotechnologies AG supplies the following controls available for use with this product:
REF CG0363 Coagnos® D-Dimer Control H/L

STORAGE, SHELF-LIFE AND STABILITY:

Coagnos® D-Dimer kit should be stored in the refrigerator at 2–8 °C. Unopened reagents and calibrators are stable until the expiration date. After opening, Coagnos® D-Dimer Reagent, Reaction Buffer and Diluent are stable for 14 days when stored tightly capped and refrigerated at 2–8 °C in the original vial. Do not freeze.

Onboard stability (Coagulyzer® Auto Pro)
Reaction buffer (POS 37) at 37 °C ≥6 hrs
D-Dimer Latex (POS 38) at 37 °C ≥ 6 hrs
Reconstituted Calibrator is stable for 8 hours at 18–26 °C and 3 days in the refrigerator at 2–8 °C in the original vial.
For optimal stability, remove reagents from the analyzer system and store refrigerated at 2–8 °C.

SAMPLE COLLECTION AND PREPARATION:

It is recommended that specimen collection and storage be carried out in accordance with NCCLS guideline H21-A3⁹ or other appropriate guideline (e.g. DIN58 905 Part 1, "Haemostaseology: blood collection; preparation of plasma from citrated blood for coagulation testing.") No known test method can offer complete assurance that human blood samples will not transmit infection, therefore, all blood derivatives should be considered potentially infectious. Venous blood is collected in 3.2% sodium citrate at a ratio of nine (9) parts blood to one (1) part anticoagulant (1:10 ratio). The ratio is critical. If using commercial vacuum tubes, a full draw must be assured. Trauma or stasis during drawing should be avoided. Blood should not be collected through a heparin lock or other heparinized line. The presence of clots in a specimen is cause for rejection. Centrifuge at 1500 x g for 15 minutes at room temperature. Unless samples are to be processed immediately, transfer plasma to a plastic test tube as soon as centrifugation is complete. Store plasma samples refrigerated 2–8 °C and test within 4 hours. Alternatively, plasma may be stored at -20 °C for up to one month. Frozen plasma samples should be rapidly thawed at 37 °C while gently mixing and tested immediately.

PROCEDURE:

Preparation of Standard Curve
Use the lot-specific value for D-Dimer Calibrator to determine the exact D-dimer concentration in each standard dilution. Users must construct a standard curve each time a new kit lot is used, after any instrument change or service, or every 6 months, and if Coagnos® D-Dimer Control H/L (REF CG0363) is assayed out of range.

Automated Method
Refer to the appropriate instrument operator manual for detailed instructions or contact Analyticon Biotechnologies AG for instrument specific application guides.

INTERPRETATION OF RESULTS:

All patient samples should be initially tested undiluted. D-dimer results are reported in ng/ml. The analytical sensitivity is directly linked to the instrument used to perform the test. For more information, refer to the instrument's application protocol, available from Analyticon Biotechnologies AG. The upper end of the reportable range for undiluted samples is a multiple of the assigned value of the calibrator.

The specific multiplier used is dependent upon the instrument application. If the result exceeds the reportable range, retest the sample using an appropriate dilution. No antigen excess hook effect is observed at D-dimer levels below 60.000 ng/ml.

Elevated D-dimer levels can persist for some time after the active process has ceased. Although the absence of elevated D-dimer levels may be used to help rule out the presence of a thrombotic disorder, a negative D-dimer test does not completely rule out the presence of thrombosis. The D-dimer results should be used with other clinical and diagnostic information for forming a diagnosis and for patient management.

LIMITATIONS:

When tested using samples prepared according to NCCLS publication EP7-P,¹⁰ levels of the following do not interfere with the Coagnos® D-Dimer assay: hemoglobin less than 1200 mg/l, bilirubin less than 19 mg/dl, heparin less than 0.6 U/ml, and triglycerides less than 2000 mg/dl. Interferences are instrument dependent, please refer to specific application protocol. The interference criterion is the concentration of the interferent at which the change, if any, in the measured value of D-dimer in ng/ml is less than or equal to 10%. Turbid or opalescent plasma may cause erratic results and should be interpreted with

REF: CG0380

caution. The presence of rheumatoid factor may produce false positive test results. Plasma samples positive for rheumatoid factor should be interpreted with caution. Human anti-mouse antibody testing was not performed with this product. Results from patients with heterophilic antibody should be interpreted with caution since this test kit contains mouse antibodies and interference may occur resulting in falsely elevated or decreased values.

QUALITY CONTROL:

Each laboratory should establish a control range to determine the allowable variation in the day to day performance of the Coagnos® D-Dimer. Failure to meet Quality Control specifications should be investigated and resolved. Recalibration is suggested whenever control materials are not within the established acceptable range. The assay should be repeated. If the problems cannot be resolved, contact Analyticon Biotechnologies AG Technical Services or your local country Analyticon Biotechnologies AG representative.

REFERENCE VALUES:

For Coagnos D-Dimer a cut-off of 115 µg/L DDU, 200 µg/L FEU is recommended. Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS:

The following performance characteristics have been determined by Analyticon Biotechnologies AG or their representatives using a Coagulyzer® Auto Pro instrument. Each laboratory should establish its own performance data.

Specificity

The monoclonal antibody used in this Kit reacts with D-dimer and is reported to react with D-monomer from fibrin.⁸ In studies performed, the antibody did not react with fibrinogen (purified or in plasma), nor did it react with plasmin or streptokinase treated plasma. Some cross reactivity was found with noncrosslinked fibrin and plasmin degradation products from purified fibrinogen. The in vivo relevance of the cross reactivity is minimized by 2 antiplasmin, as demonstrated by the lack of cross reactivity with plasmin treated plasma.

Correlation Studies

The Coagnos® D-Dimer kit (Coagulyzer® Auto Pro) and another commercially available quantitative D-dimer assay (Siemens® D-Dimer / INNOVANCE® D-Dimer / Sysmex® CS-5100) were compared using normal and abnormal samples. Results of these correlation studies are presented in the table below.

Regression equation	Y = 1.061x - 935
Correlation coefficient	0.89
Number of samples	128

Precision


Precision was performed according to CLSI guideline EP05-A3: Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition, 2008.

Sample	Mean	Repeatability		Within-device	
		SD	CV [%]	SD	CV [%]
D-Dimer Concentration (µg/L DDU)	1961	126	6.42	266	13.55
	264	23	8.80	32	12.17

Interferences

Coagnos® D-Dimer is insensitive to Haemoglobin at concentrations up to 500 mg/dL. Using a 5% interference threshold, there is no significant interference from Bilirubin at concentrations up to 15 mg/dL for Coagnos® D-Dimer. Lipid interference testing demonstrates that lipid levels do not directly affect the test up to 37 mg/dL.

SYMBOLS

	In vitro diagnostics product
	The product complies with European legislation
	Follow the instructions for use !
	Use by
	Permitted storage temperature range
	Read warnings and precautions in instructions for use !
	Batch identification number
	Item number
	Manufacturer

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BIBLIOGRAPHY/ BIBLIOGRAFÍA/ BIBLIOGRAPHIE/ LITERATURVERZEICHNIS:

1. Elms, MJ et al: Rapid detection of cross-linked fibrin degradation products in plasma using monoclonal antibody-coated latex particles. J Clin Path 85:360, 1986
2. Declerck, PV et al: Fibrinolytic response and fibrin fragment D-Dimer in patients with deep vein thrombosis, Thromb Haemost 58:1024, 1987
3. Bounameaux H, et al: Measurement of plasma D-Dimer for diagnosis of deep venous thrombosis. Am J Clin Path 91:82, 1989
4. Hansson PO, Eriksson H, Eriksson E, Jegenburg R, Risberg B: Can laboratory testing improve screening strategies for deep vein thrombosis at an emergency unit? J Intern Med 235:143, 1994
5. Demers C, Ginsberg JS, Johnston M, Brill-Edwards P, Panju A: D-Dimer and antithrombin III complexes in patients with clinically suspected pulmonary embolism, Thromb Haemost 67: 408, 1992
6. Nakashima, J et al: Tumor necrosis factor and coagulopathy in patients with prostate cancer. Cancer Res 55: 4881, 1995
7. Bailegeer, V et al: Fibrinolytic response to venous occlusion and fibrin fragment D-Dimer levels in normal and complicated pregnancy. Thromb Haemost 85:1030, 1987

8. Holvoet, P et al: Binding properties of monoclonal antibodies against human fragment D-Dimer of cross-linked fibrin to human plasma clots in an In Vivo model in rabbits. Thromb Haemost 61:307, 1989

9. National Committee for the National Laboratory (NCCLS) Standards Collection, Transport and preparation of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays Approved Guideline, 3rd Edition. NCCLS Publication H21-A3, December 1998

10. National Committee for the National Laboratory (NCCLS) Interference Testing In Clinical Chemistry Proposed Guideline, 3rd Edition; NCCLS Publication EP7-P, Vol. 6, No. 13, 1986

11. National Committee for the National Laboratory (NCCLS) Evaluation of Precision Performance of a Clinical Chemistry Devices Approved Guideline; NCCLS Publication, EP5-A, Vol. 19, No., February 1999