Coagnos® PT Calibration Plasma

REF: CG0371 4 x 1 ml

Coagnos® PT Calibration Plasma

GB

INTENDED PURPOSE:

The prothrombin time (PT) can be used to monitor oral anti-coagulation therapy. The World Health Organisation (WHO) recommends to report PT values as 'International Normalised Ratio' (INR)^{1,2}. INR can be calculated using the international sensitivity index (ISI) and mean normal prothrombin time (MNPT) of the chosen thromboplastin reagent using the following equation:

$$INR = (PT / MNPT)^{ISI}$$

Each lot of Coagnos® PT reagent is supplied with instrument specific ISI values. However, the ISI value and MNPT may be influenced by local, laboratory conditions and the instrument used. Local determination of the ISI and MNPT is therefore recommended¹⁻⁴. The Coagnos® PT kit is indicated for professional use only.

The Coagnos® PT Calibration Plasma kit contains four PT Calibrant plasmas for standardising the PT test. The plasmas are designed to simulate normal human plasma and a range of plasma pathologies. The Calibration Plasma kit is intended for use with Coagnos® PT, and each PT Calibrant plasma is assigned a %PT and INR for the following purposes:

- To establish a %PT Calibration Curve for reporting PT values as a % of normal.
- To establish an INR Reference Curve for the direct INR determination of a patient sample.
- To establish specific ISI and MNPT values for the system, reagent and instrument, used by the laboratory.

WARNINGS AND PRECAUTIONS:

The reagents contained in this kit are for in vitro diagnostic use only - DO NOT INGEST. Wear gloves when handling all kit components. The material safety data sheet is available for download from our homepage <http://www.analyticon-diagnostics.com>. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample. In case any serious incident has occurred in relation to the device, please report to the manufacturer and, if applicable, to the competent authority of the country in which the users and/or the patients established themselves.

COMPOSITION:

Component	Content	Description	Preparation
PT Calibrant 1	1 x 1 mL	The Calibration Plasma kit contains 4 PT Calibrant plasmas for standardising the PT test. The plasmas are prepared from pooled normal human plasmas.	Reconstitute with 1.0 mL of distilled or deionised water, allow the vial to stand for 15 minutes and mix gently before use to allow complete dissolution. DO NOT SHAKE.
PT Calibrant 2	1 x 1 mL	PT Calibrant 1 simulates normal human plasma, PT Calibrant 2, PT Calibrant 3 and PT Calibrant 4 simulate a range of plasma pathologies and are prepared through the absorption of clotting factors.	
PT Calibrant 3	1 x 1 mL	The assigned values for each of the plasma are indicated on the insert sheet expressed as %PT and INR, the values are intended for use with Coagnos® PT and are instrument / instrument series specific.	
PT Calibrant 4	1 x 1 mL		

STORAGE, SHELF-LIFE AND STABILITY:

Unopened vials should be stored at 2 – 8 °C and are stable until the expiration date stated on the labels. Reconstituted plasma must be used within 1 hour.

SAMPLE COLLECTION AND PREPARATION:

Not applicable.

PROCEDURE:

The reconstituted PT Calibrant plasma material should be treated the same as patient samples, following instrument specific protocols.

The PT Calibrant plasmas can be used on specific automated instrumentation as identified in the lot specific reference values insert. It may be possible to input the assigned values of the PT Calibrant plasmas allowing the instrument to interpret %PT and/or INR results automatically. Carefully read and follow the operating procedures for the specific instrument.

INTERPRETATION OF RESULTS:

Establishing a %PT Calibration Curve:

Determine the coagulation times of each of the 4 PT Calibrant plasmas in duplicate or triplicate.

Plot the mean value obtained for each PT Calibrant plasma on double log paper (x-axis = %PT value; y-axis = Coagulation time in seconds).

Join the points (point – point) and determine the PT coagulation time of the patient plasma and directly read from this reference line the corresponding %PT value of patient plasma.

Establishing an INR Reference Curve for the direct INR determination of a patient sample:

Determine the coagulation times of each of the 4 PT Calibrant plasmas in duplicate or triplicate.

Plot the mean value obtained for each PT Calibrant plasma on double log paper (x-axis=INR value; y-axis=Coagulation time in seconds).

Join the points by fitting the best possible straight line through these points.

Determine the PT coagulation time of the patient plasma and directly read from this reference line the corresponding INR value of patient plasma.

Establishing Laboratory Specific ISI and MNPT:

Determine the coagulation times of each of the 4 PT Calibrant plasmas in duplicate or triplicate and calculate the mean value for each of the plasma.

A linear relationship exists between the Log INR value (x-axis) and Log PT (sec; y-axis), expressed by the equation:

$$\text{Log PT (sec)} = [\text{slope} \times \text{Log (INR)}] + \text{intercept}$$

From this equation the laboratory specific ISI and MNPT can be calculated in the following way:

$$ISI = 1/\text{slope}$$

$$MNPT = 10\text{intercept}$$

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Analyticon Biotechnologies for instrument specific application guides.

LIMITATIONS:

A new calibration is required for each batch of Coagnos® PT and for each instrument used. Recalibration is also recommended if software changes are introduced or following a major service or repair of the instrument. The INR and %PT values of the PT Calibrant plasma samples supplied with this kit are lot specific. This kit is for use solely in conjunction with Coagnos® PT (REF: CG0301, CG0302) on the instruments listed on the accompanying assay sheet.

QUALITY CONTROL:

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Analyticon Biotechnologies supply the following controls available for use with this product:

REF CG0351 Coagnos® Routine Control N

REF CG0352 Coagnos® Routine Control P

REF CG0361 Coagnos® Special Control N

REF CG0362 Coagnos® Special Control P

REFERENCE VALUES:

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

SYMBOLS

In vitro diagnostics product

The product complies with European legislation

Follow the instructions for use !

Use by

Permitted storage temperature range

Read warnings and precautions in instructions for use !

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Hersteller

BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/ LITERATURVERZEICHNIS:

1. Babson AL and Flanagan ML (1975) Quantitative One Stage Assays for Factors V and X, AJCP, 64: 817-819.
2. Hardisty RM et al. (1962) A One Stage Factor VIII Assay and Its Use on Venous and Capillary Plasma. Thrombosis et Diathesis Haemorrhagica, 7:215-229.
3. Morse EE et al. (1971) Automated Fibrinogen Determination, AJCP, 55:671-676.
4. Elodi S et al (1978) Some Sources of Error in the One-Stage Assay of Factor VIII, Haemostasis, 7:1-9.

Coagnos[®] PT Calibration Plasma

REF: CG0371 4 x 1 ml

5. Thelin M (1968) Preparation and Standardization of a Stable AHF Plasma, Thrombosis et Diathesis Haemorrhagica, 19:423.
6. Kirkwood TBL et al. (1977) Identification of Sources of Variation in Factor VIII Assay, British Journal of Haematology, 37:559-568.
7. Goldenfarb MD (1971) Reproducibility in Coagulation Assays, AJCP, 55:561-564.
8. Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, AJCP, 59:231-235.

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In Verbindung mit diesem Produkt bietet Analyticon Biotechnologies AG die folgenden Kontrollen an: