REF: CG0301 10 x 2 ml
CG0302 10 x 4 ml

Coagnos® PT

GB

INTENDED PURPOSE:

The Coagnos® PT kit is intended for carrying out clot based haemostasis assays in human plasma. The Coagnos® PT kit is indicated for professional use only. The first standardised one-stage prothrombin time test was developed by Dr. Armand Quick in 1935. It has now become the basic coagulation screening test for the diagnosis of congenital and acquired deficiencies of clotting factors from the extrinsic pathway (factors II, V, VII and X)^{1,2}. It is also used for the induction and monitoring of oral anticoagulant therapy^{3,4} and can be used to assess the protein synthesis capability of the liver in chronic or acute hepatic disorders. Thromboplastin is of rabbit brain origin but resembles human preparations in its low International Sensitivity Index (ISI). The ISI of Thromboplastin is approximately 1.1 and is calibrated against the WHO international reference preparation⁵. Thromboplastin is particularly suited to the monitoring of oral anticoagulant therapy and, in conjunction with the appropriate factor deficient plasma, the measurement of factor activity in the extrinsic pathway. Tissue thromboplastin, in the presence of calcium ions, is an activator which initiates the extrinsic pathway of coagulation. When a mixture of tissue thromboplastin and calcium ions is added to normal citrated plasma, the clotting mechanism is activated, leading to a fibrin clot. If a deficiency exists within the extrinsic pathway, the time required for clot formation will be prolonged depending on the severity of the deficiency.

WARNINGS AND PRECAUTIONS:

The reagents contained in this kit are for in vitro diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. The material safety data sheet is available for download from our homepage <http://www.analyticon-diagnostics.com>. Dispose of components in accordance with local regulations. In case any serious incident has occurred in relation to the device, please report to the manufacturer and, if applicable, to the competent authority of the country in which the users and/or the patients established themselves.

COMPOSITION:

Component	Content	Description	Preparation
Thromboplastin	10 x 2 mL (REF CG0301) 10 x 4 mL (REF CG0302)	Liquid Rabbit Brain Thromboplastin containing Calcium Chloride, stabilisers and preservatives.	The liquid, calcified thromboplastin is ready-for-use. No further calcium is required to carry out standard PT Assays. The contents of the vial should be mixed well before use. (5 minutes on roller).

STORAGE, SHELF-LIFE AND STABILITY:

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Opened vials are stable for 2 months at +2 – +8°C, and 6 hours at +37°C (on-board Coagulyzer® Auto Pro). DO NOT FREEZE. Large clumps of particles or changes in expected values may indicate product deterioration.

SAMPLE COLLECTION AND PREPARATION:

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at +18 – +24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at +37°C prior to testing. Do not keep at +37°C for more than 5 minutes⁶.

PROCEDURE:

For accurate INR reporting, it is recommended to determine the laboratory specific ISI of the reagent with the testing system in use. The Coagnos® PT Calibration Plasma (REF CG0371) is recommended for this purpose^{7,8}. This should be performed for each new reagent batch. The Coagnos® Routine Control N (REF CG0351) and Coagnos® Routine Control P (REF CG0352) should be used to check for shifts in the local system ISI which have been noted with changes in laboratory temperature and post instrument servicing, amongst other local variances.

Manual Method

- Mix sufficient Coagnos® PT to complete the anticipated testing for the day and incubate at +37°C for no more than 4 hours.
- Pre-warm 0.1 mL of the test plasma at +37°C for 2 minutes.
- Add 0.2 mL of freshly mixed thromboplastin reagent to the plasma while simultaneously starting a stopwatch.
- Note the time for clot formation to the nearest 0.1 seconds.

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Analyticon Biotechnologies AG for instrument specific application guides.

INTERPRETATION OF RESULTS:

Results should be reported to the nearest 0.1 seconds and duplicates should agree within 5% of each other. %PT values can be interpolated from the calibration graph (%PT of PT Calibration plasmas versus measured clot time), which should be a straight line when plotted on log-log graph paper.

INR values can be calculated using the following formula:

$$\text{INR} = (\text{PT Time Patient} / \text{Mean Normal PT Time})^{\text{ISI}}$$

Analyticon® Biotechnologies AG
Am Mühlenberg 10, 35104 Lichtenfels / Germany

For clear guidance on the indications for and management of patients on warfarin, please refer to The British Society for Haematology, for their most current edition of 'Guidelines on oral anticoagulation with warfarin'. At time of printing this is the 2011 fourth edition⁹.

LIMITATIONS:

The use of serial dilutions of a reference plasma for the %PT curve is not recommended as this can lead to discrepancies caused by the low fibrinogen in the reference plasma dilutions which are not reflected in patient samples having predominantly normal fibrinogen levels. Analyticon Biotechnologies AG advise use of the Coagnos® PT Calibration Plasma kit (REF CG0371) for this purpose.

QUALITY CONTROL:

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Analyticon Biotechnologies AG supplies the following controls available for use with this product:

REF CG0361 Coagnos® Special Control N
REF CG0362 Coagnos® Special Control P
REF CG0351 Coagnos® Routine Control N
REF CG0352 Coagnos® Routine Control P

REFERENCE VALUES:

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges. This is particularly important for local ISI calibration. Using the Coagulyzer® Auto Pro, normal values ranging from 12.10 - 15.80 seconds; 0.870 - 1.140 INR; 85.2 - 128.50 %PT are typical.

PERFORMANCE CHARACTERISTICS:

The following performance characteristics have been determined by Analyticon Biotechnologies AG or their representatives using a Coagulyzer® Auto Pro instrument. Each laboratory should establish its own performance data.

Precision

Precision was performed according to CLSI guideline EP05-A3: Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition, 2008.

Sample	Mean	Repeatability		Within-device	
		SD	CV [%]	SD	CV [%]
Clots [s]	13.2	0.30	2.24	0.49	3.74
	27.1	0.39	1.45	1.06	3.93
	44.8	1.23	2.74	2.17	4.83
%PT	108.9	3.60	3.31	5.26	4.83
	40.8	0.80	1.97	2.83	6.93
	22.9	0.67	2.90	1.53	6.69
INR	0.99	0.02	2.50	0.04	3.73
	2.10	0.03	1.49	0.08	3.98
	3.55	0.10	2.78	0.20	5.61

Interferences

Coagnos® PT is insensitive to Heparin levels of up to 2 U/ml. Using a 5% interference threshold, there is no significant interference from Haemoglobin at concentrations up to 250 mg/dL. Using a 5% interference threshold, there is no significant interference from Bilirubin at concentrations up to 30 mg/dL for Coagnos® PT. Lipid interference testing demonstrates that lipid levels do not directly affect the clot time of the reagent up to 500 mg/dL. Lipid concentrations in excess of this prevent clot detection.

Method comparison

Comparison of clot time in seconds and INR values were determined using Coagnos® PT on Coagulyzer® Auto Pro and Siemens® PT with Thromborel® on a Sysmex® CS-5100 system on 106 samples. The following correlations were obtained:

$$\begin{aligned} \text{Coagnos® PT (Seconds)} &= 1.180x + 0.902 & r^2 &= 0.96 & n &= 106 \\ \text{Coagnos® PT (INR)} &= 1.026x + 0.165 & r^2 &= 0.96 & n &= 107 \end{aligned}$$

SYMBOLS

	In vitro diagnostics product
	The product complies with European legislation
	Follow the instructions for use !
	Use by
	Permitted storage temperature range
	Read warnings and precautions in instructions for use !

REF: CG0301 10 x 2 ml
CG0302 10 x 4 ml

LOT Batch identification number
REF Item number
 Manufacturer

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Hersteller

BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/ LITERATURVERZEICHNIS:

- Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, Am. J. Med. Sci, 190: 501.
- Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.
- Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, Chest, 95: 5S-11S.
- Poller L (1986) Laboratory Control of Anticoagulant Therapy, Sem. Thromb. Haemostasis, 12: 13-19.
- World Health Organisation (1984) Expert Committee on Biological Standards, Technical Series, 700: 19.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Poller L., Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study, Amer. J. Clin. Pathol, 103: 358-365.
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios, Amer. J. Clin. Pathol, 103: 366-371.
- Keeling D (2011) Guidelines on Oral Anticoagulation with warfarin: Forth Edition, British Journal of Haematology, 154(3): 311-324.